

ROHO - Mosaic



- Designed to provide a comfortable, functional and stable sitting environment for patients who are at risk of skin breakdown
- The Mosaic is also a Dry Floatation cushion but is made from softer plastic than the classic ROHO DF cushions. It is therefore lighter (0.5kg) but less rugged
- Interconnected air cells disperse pressure and promote stability as the pelvis sinks into them
- The Mosaic is also excellent as a lounge-chair cushion and can also be used as a shower chair cushion
- Maximum User Weight: 113 kg
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions

Warranty Information

- 1 Year

✓ Contracted DVA productZ

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Code	Display Name
PTC586100	ROHO Mosaic - 16 x 16in

Code	Display Name
PTC586200	ROHO Mosaic - 18 x 16in