

Aspire Premium Community Care Bed



- Combines a sleek stylish design with superior comfort and functionality
- Designed with personal design attention allowing it to fit subtly into any home environment
- Head and foot boards with adjustable side rails in a beech timber finish form a traditional bed surround when lowered to compliment a homely interior and care facilities alike
- This robust bed can also be easily assembled and dismantled for ease of transport and storage when not in use
- Designed for the home environment with 5 profiling functions to suit all comfort and positioning needs
 - High-low
 - Backrest (60 °)
 - Knee bend (30 °)
 - Trendelenburg (10.5 °)
 - Combined backrest and knee bend
- Beech full length side rails
- Beech head and foot boards
- Integrated bed extension
- Four total-lock castors
- Self-help pole
- Mattress retainers
- TiMotion electronics
- Large lifting range
- Mattress dimensions: 2000 x 900mm
- Mattress platform length: 2020-2120mm



Contract approved

NDIS

Item Code:

05_181210197_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

style='font-size: 13pt;'

style='color: #48a841;'**Accessories**

Click on the cart icon to display the range of optional accessories

Warranty Information

- 3 Years Frame
- 2 Years Electronics
- 1 Year Handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB046050K	Premium Community Care Bed - 5 Function	2140 - 2240 mm	1025 mm	220 - 710 mm	180 kg	5