


## **Aspire Premium Community Care Bed**





- Combines a sleek stylish design with superior comfort and functionality
- Designed with personal design attention allowing it to fit subtly into any home environment
- Head and foot boards with adjustable side rails in a beech timber finish form a traditional bed surround when lowered to compliment a homely interior and care facilities alike
- This robust bed can also be easily assembled and dismantled for ease of transport and storage when not in use
- Designed for the home environment with 5 profiling functions to suit all comfort and positioning needs
  - High-low
  - Backrest (60 º)
  - Knee bend (30 º)
  - Trendelenburg (10.5 º)
  - · Combined backrest and knee bend
- Beech full length side rails
- · Beech head and foot boards
- Integrated bed extension
- Four total-lock castors
- Self-help pole
- · Mattress retainers
- TiMotion electronics
- Large lifting range
- Mattress dimensions: 2000 x 900mm
- Mattress platform length: 2020-2120mm

## Warranty Information



## ndis

**Item Code:** 

05\_181210197\_0103\_1\_2

Power Bed – Composite:

05\_423\_0103\_1\_2

Click **here** for more information

## **Accessories**

Click on the cart icon to display the range of optional accessories



- 3 Years Frame
- 2 Years Electronics
- 1 Year Handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB046050K	Premium Community Care Bed - 5 Function	2140 - 2240 mm	1025 mm	220 - 710 mm	180 kg	5