

ROHO - Mosaic





Code Display Name

PTC586100 ROHO Mosaic - 16 x 16in

PTC586200 ROHO Mosaic - 18 x 16in

- Designed to provide a comfortable, functional and stable sitting environment for patients who are at risk of skin breakdown
- The Mosaic is also a Dry Floatation cushion but is made from softer plastic than the classic ROHO DF cushions. It is therefore lighter (0.5kg) but less rugged
- Interconnected air cells disperse pressure and promote stability as the pelvis sinks into them
- The Mosaic is also excellent as a lounge-chair cushion and can also be used as a shower chair cushion
- Maximum User Weight: 113 kg
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions

Warranty Information

• 1 Year



ndis

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion -

Composite: 05_429_0103_1_2 Click <u>here</u> for more information