

Liko Viking XL Lifter





- Ideal for bariatric lifts and transfers
- Universal TwinBar 670 Slingbar offers additional sling space for patients of size
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 540 1910 mm

Warranty Information

Frames: 10 Years

Electrical Parts: 2 Years

\bigcirc	Cambuant	approved
	Contract	approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist – Composite:

05_420_0105_1_2

Click **here** for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393170	Liko Viking Lifter - XLARGE	1420 mm	815 - 1130 mm	1500 - 2110 mm	300 kg	Aluminium



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393220 Li	ko Viking XI - Leg Protector Set		_	_	_	_