

Aspire ComfiMotion Luxe Bed



- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
 - **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
 - **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
 - **Box-on-box design** facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
 - **USB Charge:** Connect and charge your devices conveniently.
 - German engineered motors for smooth, quiet adjustment.
 - ***Total Height with:**
 - Low Profile Fixed Legs 280 – 580 mm (available upon request)
 - 75mm Standard Castors 350 - 650 mm
 - 125mm Accessibility Castors 400 - 700 mm (optional)
- Mattress choice may increase or decrease overall bed height.



Compliance

IEC 60601-2-52:2009

Warranty Information

- 10 Years Limited Warranty

Code	Display Name
BEB046710	Aspire ComfiMotion Luxe Bed - Portobello Grey - King Single (Headboard Sep)