

Aspire ComfiMotion Luxe Bed





- Massage: Relaxing vibration modes for personalised neck, shoulder and back relief.
- Personalised Positioning: Hand control allows users to find their perfect position to maximise comfort and relaxation.
- Night Assist: Luminous underbed halo for safe entry and exit throughout the night.
- Box-on-box design facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- USB Charge: Connect and charge your devices conveniently.
- German engineered motors for smooth, quiet adjustment.
- *Total Height with:
 - Low Profile Fixed Legs 280 580 mm (available upon request)
 - 75mm Standard Castors 350 650 mm
 - 125mm Accessibility Castors 400 700 mm (optional)

Mattress choice may increase or decrease overall bed height.

Warranty Information

• 10 Years Limited Warranty



Compliance

IEC 60601-2-52:2009

Code

Display Name

BEB046710 Aspire ComfiMotion Luxe Bed - Portobello Grey - King Single (Headboard Sep)