

## Greiner Dialysis & Oncology Chair Multiline - Next DC



- An infinitely adjustable 5 motor treatment chair (seat height, back recline, seat tilt, leg raise, footrest)
- Tailored to meet the specific needs of oncology departments and the demands of infusion and transfusion medicine
- Innovative armrest technology with individual positioning for height, inclination and direction. This makes the Multiline Next DC ideal for procedures such as chemotherapy; dialysis; blood management; plasmapheresis
- Footrest with automatic, patient-friendly length adjuster (removable)
- Soft upholstery provides comfortable seating for patients of all sizes (60cm seat width)
- Easy access height of 54cm
- Its design is compliant with stringent infection control regulations - clean surface finish and minimum seams
- 10 cm twin-wheel castors with central locking system if required
- Backrest Angle: - 5° to 78°
- Seat Angle: 0° to 22°
- Legrest Angle: 0° to 25°



### Compliance

IEC 60601-2-52:2009

## Warranty Information

- 2 Years

Code	Display Name	Seat Width	Safe Working Load
CHT100200	Multiline Next DC (Light Grey)	600 mm	220 kg

Code	Display Name	Seat Width	Safe Working Load
CHT100210	Multiline Next DC (Navy Blue)	600 mm	220 kg