

## **Greiner Dialysis & Oncology Chair Multiline - Next DC**





- An infinitely adjustable 5 motor treatment chair (seat height, back recline, seat tilt, leg raise, footrest)
- Tailored to meet the specific needs of oncology departments and the demands of infusion and transfusion medicine
- Innovative armrest technology with individual positioning for height, inclination and direction. This makes the Multiline Next DC ideal for procedures such as chemotherapy; dialysis; blood management; plasmapheresis
- Footrest with automatic, patient-friendly length adjuster (removable)
- Soft upholstery provides comfortable seating for patients of all sizes (60cm seat width)
- Easy access height of 54cm
- Its design is compliant with stringent infection control regulations clean surface finish and minimum seams
- 10 cm twin-wheel castors with central locking system if required
- Backrest Angle: 5° to 78°
- Seat Angle: 0° to 22°
- Legrest Angle: 0° to 25°

## Warranty Information

220 kg

2 Years



Code	Display Name	Seat Width	Safe Working Load
------	--------------	------------	-------------------

CHT100210 Multiline Next DC (Navy Blue) 600 mm 220 kg