

Aspire ComfiMotion Care Bed






- Night Assist: Luminous underbed halo for safe entry and exit throughout the night.
- Press and hold safety handset: Intuitive, large button interface enhances usability.
- Enclosed frame design: Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo functions only lower under the force of gravity reducing potential for accidents or injuries.
- Box-on-box design facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- German engineered motors for smooth, quiet adjustment.
- *Total Height with:
 - Low Profile Fixed Legs 280 580 mm (available upon request)
 - 75mm Standard Castors 350 650 mm
 - 125mm Accessibility Castors 400 700 mm (optional)

Mattress choice may increase or decrease overall bed height.

Warranty Information

• 10 Years Limited Warranty

Ω		
Ocontracted	DVA	productZ



Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click **here** for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Height	Shoe Size
BEB046720	Aspire ComfiMotion Care Bed - Portobello Grey - Long Single (Headboard Sep)	280 - 580 mm	10

BEB046730 Aspire ComfiMotion Care Bed - Portobello Grey - King Single (Headboard Sep) 280 - 580 mm 11

