

Aspire ComfiMotion Fixed Companion Bed



- Matching single bed without any adjustability designed to sit flush when positioned together with the ComfiMotion CARE Bed and ComfiMotion LUXE Bed.
- It facilitates a side-by-side sleeping arrangement at a fixed height of 370mm, without compromising on the users care and assistance requirements.

Warranty Information

- 1 Year



✓ Contracted DVA productZ

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046800	Aspire ComfiMotion Fixed Companion Bed - Portobello Grey - Long Single (Headboard Sep)	2030 mm	920 mm	370 mm	180 kg