

Aspire ComfiMotion Memory Foam Mattress



- The ComfiMotion Memory Foam mattress contours to the users body, providing personalised comfort and transitional pressure relief.
- Integrated side-walls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit.
- Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.
- Superior breathability: Up to 8x more breathable than conventional memory foams.
- Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.
- Eco-friendly GECA Foam is produced using cleaner chemicals and production technology.

✓ Contracted DVA productZ

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Compliance

BS7177:2008



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM469910	Aspire ComfiMotion Memory Foam Mattress - Portobello Grey (sides) - Long Single	2030 mm	920 mm	200 mm	180 kg
BEM469920	Aspire ComfiMotion Memory Foam Mattress- Portobello Grey (sides) - King Single	2030 mm	1070 mm	200 mm	180 kg
BEM469930	Aspire ComfiMotion Memory Foam Mattress w/Incontinence Cover - Long Single	2030 mm	920 mm	200 mm	180 kg
BEM469940	Aspire ComfiMotion Memory Foam Mattress w/Incontinence Cover - King Single	2030 mm	1070 mm	200 mm	180 kg