

## **Liko Viking L Lifter**





- An everyday flexible lifting solution
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 460 1790 mm

## Warranty Information

• Frames: 10 Years

Electrical Parts: 2 Years

$\bigcirc$	Contract	approved
_	Contract	approved

## ndis

**Item Code:** 

05\_123603161\_0105\_1\_2

**Hoist – Composite:** 

05\_420\_0105\_1\_2

Click **here** for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393180	Liko Viking Lifter - LARGE	1360 mm	695 - 1070 mm	1390 - 2030 mm	250 kg	Aluminium
LSS393210	Liko Viking L - Leg Protector Set	-	-	-	-	-