

Ceiling Hoist - Liko - Multirall 200





- The Multirall 200 hoist is a portable ceiling hoist that can also be used for room-to-room patient transfers
- Can be used to transfer a patient between rooms in a sitting position, easily moving from one rail system to another with the use of a transfer strap
- It's portability makes it an excellent and economical solution for nursing home and hospital care situations where lifting needs are temporary
- Can be easily attached or removed with no tools and minimum effort from ceiling rails to facilitate convenient transport between rooms
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Can be mounted two ways for ease of use

Warranty Information

•	France. 10 fears
•	Electrical Parts: 2 Years

\bigcirc	
⊘ Contract	approved

Code	Display Name	Total Height	Safe Working Load

LIC359100 Multirall 200 - Portable Ceiling Hoist Up to 1600 mm 200 kg