

**OCT-NOV
2020**



\$ALE!



Aspire Assist 2 Wheelchair

Aspire Assist 2 Ward Vinyl Wheelchair
(Folding Wheelchair - Silver - 140kg - 350mm - MWS449340)
(Folding Wheelchair - Silver - 140kg - 400mm - MWS449350)
(Folding Wheelchair - Silver - 140kg - 450mm - MWS449360)
(Folding Wheelchair - Silver - 140kg - 500mm - MWS449370)

Aspire Assist 2 Ward Vinyl HD Wheelchair
(Folding Wheelchair - Silver - 180kg - 570mm - MWS449380)
(Folding Wheelchair - Silver - 180kg - 610mm - MWS449390)



WARD VINYL
\$399
SAVE: \$126

WARD VINYL HD
\$499
SAVE: \$136



SINGLE
\$899
SAVE: \$151



Aspire Active Air 4 Mattress Overlay
(Single - 19 Cells - 2000x850x100mm - 135kg - PTM639100)



SINGLE
\$1,250
SAVE: \$645

KING SINGLE
\$1,699
SAVE: \$551



SUPER SAVER

Aspire Active Air 5 Mattress Overlay
(Single - 16 Cells - 2000x850x130mm - 175kg - PTM639200)
(King Single - 16 Cells - 2000x1050x130mm - 175kg - PTM639210)

**OCT-NOV
2020**



\$ALE!



**Aspire Tilt & Recline
Shower Commode**

**SUPER
SAVER**



**Pictured with optional open commode seat. Sold Separately (See Below)*

460MM
\$2,250
SAVE: \$900

530MM
\$2,450
SAVE: \$900

Aspire Tilt & Recline Shower Commode

Attendant Propelled - Treated Steel - 250kg - 460mm - BTC066030
Attendant Propelled - Treated Steel - 250kg - 530mm - BTC066040



460MM
\$98
SAVE: \$32

530MM
\$125
SAVE: \$41



Aspire Commode Seat
Open - 460mm - BTC066600
Open - 530mm - BTC066750



470MM
\$72
SAVE: \$37



Aspire Shower Chair Regular w/Arms
470mm - Treated Steel - 200kg - BTS118000



525MM
\$92
SAVE: \$57



**WIDER
FIT!**

Aspire Shower Chair Wide w/Arms
525mm - Treated Steel - 200kg - BTS118010

1300 133 120

ENTERPRISE PARTNER PROMOTIONS

www.aidacare.com.au

* All prices displayed exclude GST