

ComfiMotion


# Quick Selection Guide



## Luxe Adjustable Bed




- Hi-Lo
- Back Raise
- Leg Raise
- Zero Gravity
- Underbed lighting
- Head Tilt
- Memory functionality
- Massage functionality
- Twin USB ports (both sides)
- Foot and side mattress retainers (incl.)



## Care Adjustable Bed




- Hi-Lo
- Back Raise
- Leg Raise
- Zero Gravity
- Underbed lighting
- Clinical Positioning
- Auto Seating Position
- Trendelenburg (and Reverse)
- Foot and side mattress retainers (incl.)



# Quick Selection Guide


To build your ComfiMotion configuration, simply select your bed base, mattress, headboard and accessories.



Create your ComfiMotion ...




Bed Function



Headboard Style



Mattress Support



Accessories



Portobello Grey Upholstery



Allure Headboard



Chesterfield Headboard

	Code	Description	Select
BED	BEB046710	LUXE Adjustable Bed - (KS) - Portobello Grey (Headboard Sep.)	<input type="radio"/>
	BEB046720	CARE Adjustable Bed w/ Tilt - (LS) - Portobello Grey (Headboard Sep.)	<input type="radio"/>
	BEB046730	CARE Adjustable Bed w/ Tilt - (KS) - Portobello Grey (Headboard Sep.)	<input type="radio"/>
	BEB046800	FIXED COMPANION Bed - (LS) - Portobello Grey (Headboard Sep.)	<input type="radio"/>
MATTRESS	BEM469910	Memory Foam Mattress - <b>Classic (LS)</b> - Portobello Grey (sides)	<input type="radio"/>
	BEM469920	Memory Foam Mattress - <b>Classic (KS)</b> - Portobello Grey (sides)	<input type="radio"/>
	BEM469930	Memory Foam Mattress - <b>Incontinence (LS)</b> - Portobello Grey (sides)	<input type="radio"/>
	BEM469940	Memory Foam Mattress - <b>Incontinence (KS)</b> - Portobello Grey (sides)	<input type="radio"/>
	BEM469800	Pocket Spring Mattress - (LS) - Portobello Grey (sides)	<input type="radio"/>
	BEM469810	Pocket Spring Mattress - (KS) - Portobello Grey (sides)	<input type="radio"/>
HEADBOARD	BEA830120	Allure Headboard - (LS) - Portobello Grey	<input type="radio"/>
	BEA830130	Allure Headboard - (KS) - Portobello Grey	<input type="radio"/>
	BEA830100	Chesterfield Headboard - (LS) - Portobello Grey	<input type="radio"/>
	BEA830110	Chesterfield Headboard - (KS) - Portobello Grey	<input type="radio"/>
ACCESSORIES	BEA830200	Crook Handle Bed Stick	<input type="radio"/>
	BEA830210	Bed Rail	<input type="radio"/>
	BEA830220	Battery Backup	<input type="radio"/>
	MWT451080	Large Accessibility Castor (x4 required)	<input type="radio"/>