




ComfiMotion Collection

Luxurious adjustable bed collection for an enhanced sleeping experience




aspire

ComfiMotion



▲ Chesterfield Headboard



ComfiMotion collection


We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.


Create your ComfiMotion ...



**Bed
Function**



**Headboard
Style**




**Mattress
Support**




Accessories

Luxe

Adjustable Bed



◀ Chesterfield Headboard



Supporting the modern, independent lifestyle.

The LUXE Adjustable bed embodies style, function and comfort. Multi-positional settings, invigorating massage and USB connectivity allow you to unwind and recharge.

- Life just became a little more LUXE.

Features

- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Box-on-Box Design** facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- **USB Charge:** Connect and charge your devices conveniently.
- **German Engineered Motors** for smooth, quiet adjustment.

Functions

Hi-Lo: Assists with bed entry and exit by raising the mattress surface to the ideal height.

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Leg Raise: Promotes circulation, relieving tired, swollen and/or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.


Memory Function: Allows the user to set and return to their preferred positions at any time with one click of a button.

Upholstery



Portobello Grey

Variable Bed Height Options



LOW HEIGHT
280 - 580mm with Low Profile Bed Legs

STANDARD
350 - 650mm with Standard Castors

ACCESSIBILITY
400 - 700mm with Large Castors (optional)



▲ Companion Sleeping

Sizes

BEBO46710 King Single
2030 x 1070 x 280 - 580mm¹

10 YEARS LIMITED WARRANTY²



SWL

1. With Low Profile Bed Legs. 2. Refer to the User Manual for Warranty T&Cs.

Care

Adjustable Bed



For comfort you deserve.

The CARE Adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.

- Make life easier with a ComfiMotion CARE bed.

Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Press and Hold Safety Handset:** Intuitive, large button interface enhances usability.
- **Enclosed Frame Design:** Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- **Safe-Lowering Technology:** Ensures that back, leg and Hi-Lo functions only lower under the force of gravity reducing potential for accidents or injuries.
- **Box-on-Box Design** facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors** for smooth, quiet adjustment.



Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Clinical Positioning: One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

Leg Raise: Promotes circulation, relieving tired, swollen or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.


*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Upholstery



Portobello Grey

Variable Bed Height Options



▲ Companion Sleeping

Sizes

BEB046720	Long Single 2030 x 920 x 280 - 580mm ¹	BEB046730	King Single 2030 x 1070 x 280 - 580mm ¹
-----------	---	-----------	--

10 YEARS LIMITED WARRANTY²



SWL


1. With Low Profile Bed Legs. 2. Frame only. Refer to the User Manual for Warranty T&Cs.

Mattresses

Engineered for Motion

Pocket Spring

The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality.




Transitional Layered Pocket Springs arranged in weight activated support zones.



Integrated Side-walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



Luxuriously Soft Pillowtop Layer provides a plush immersive feeling complimented by the responsive pocket springs.



▲ Classic Pocket Spring Mattress

Our ComfiMotion Pocket Spring and Memory Foam mattresses are designed with the flexibility to follow the contouring profile of an adjustable bed.

Memory Foam

The ComfiMotion Memory Foam mattress contours to the users body, providing personalised comfort and transitional pressure relief. Integrated side-walls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit.



Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.




Superior Breathability: Up to 8x more breathable than conventional memory foams.



Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.



Eco-friendly GECA Foam is produced using cleaner chemicals and production technology.



▲ Classic Memory Foam Mattress

Also available with Incontinence Cover



Pocket Spring

BEM469800

Long Single - Classic
2030 x 920 x 240mm

BEM469810

King Single - Classic
2030 x 1070 x 240mm

5 YEARS LIMITED WARRANTY



SWL

Memory Foam

BEM469910

Long Single - Classic
2030 x 920 x 200mm

BEM469930

Long Single - Incontinence
2030 x 920 x 200mm

BEM469920

King Single - Classic
2030 x 1070 x 200mm

BEM469940

King Single - Incontinence
2030 x 1070 x 200mm

10 YEARS LIMITED WARRANTY






SWL

Visit comfirmotion.com.au for Warranty T&Cs


Configurations for Companions

Whichever ComfiMotion package you select rest assured there is a companion solution available to suit couples, partners and caregivers wishing to sleep alongside any ComfiMotion Adjustable bed.

What to Order


Configurations	Bed 1		Bed 2
 KING 2030 x 1840	● OR ●	+	● OR ●
 KING PLUS 2030 x 1990	● OR ●	+	● OR ●
 SUPERKING 2030 x 2140	● OR ●	+	● OR ●

● CARE (LS) BEB046720	● CARE (KS) BEB046730
● Fixed Companion (LS) BEB046800	● LUXE (KS) BEB046710



Dual Adjustable Beds ● ● ●

Designed to sit flush when positioned together to enable close monitoring and enhanced care. This configuration offers the flexibility of independent adjustability, without compromising on the sense of inclusivity and intimacy between partners.



Fixed Companion Bed ●


Matching single bed without any adjustability designed to sit flush when positioned together. It facilitates a side-by-side sleeping arrangement at a fixed height of 370mm, without compromising on the users care and assistance requirements.

BEB046800 Long Single - Fixed Companion
2030 x 920 x 370mm

Headboards




Allure
Long Single BEA830120
King Single BEA830130




Chesterfield
Long Single BEA830100
King Single BEA830110


Accessories




Crook Handle Bed Stick
BEA830200
Provides multiple grip points to assist repositioning in bed, and safe entry and exit.



Bed Rail
BEA830210
Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.







Large Accessibility Castors (125mm)
MWT451080 (4x required)
Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access.



Battery Backup
BEA830220
Provides peace of mind in the event of a power outage.

Collection Summary

 Bed	BEB046720	CARE Adjustable Bed w/ Tilt - (LS)	 Headboard	BEA830120	Allure Headboard - (LS)
	BEB046730	CARE Adjustable Bed w/ Tilt - (KS)		BEA830130	Allure Headboard - (KS)
	BEB046710	LUXE Adjustable Bed - (KS)		BEA830100	Chesterfield Headboard - (LS)
	BEB046800	Fixed Companion Bed - (LS)		BEA830110	Chesterfield Headboard - (KS)
 Mattress	BEM469800	Pocket Spring Mattress - (LS)	 Accessories	BEA830200	Crook Handle Bed Stick
	BEM469810	Pocket Spring Mattress - (KS)		BEA830210	Bed Rail
	BEM469910	Memory Foam Mattress - Classic (LS)		MWT451080	Large Accessibility Castors
	BEM469920	Memory Foam Mattress - Classic (KS)		BEA830220	Battery Backup
	BEM469930	Memory Foam Mattress - Incontinence (LS)			
BEM469940	Memory Foam Mattress - Incontinence (KS)				

Key:
LS - Long Single
KS - King Single

Clinical Considerations

Accessibility and Inclusivity: Unexposed bed surround with padded upholstery removes gaps between beds in co-sleeping arrangements as well as promotes safe entry and exit.


User Independence: Height adjustability improves user mobility and assists carers with manual handling.

Flexible Positioning: Caters to various user requirements such as respiratory conditions, snoring, sleep apnoea, pressure redistribution, joint or muscle pain, lymph oedema and swelling.

Vision and Hand Function: Handset with a clear and intuitive button display, assists users with reduced vision or hand function.

Fall Prevention: Integrated lighting and a minimum bed height of 280mm helps to reduce the likelihood and severity of falls in the home.

User Safety: Battery backup (optionally available) provides peace of mind in the event of a power outage.



ASBR28.54 280720

1300 133 120
www.comfimotion.com.au

 **aspire**
ComfiMotion