

# aspire LIFT RECLINE CHAIRS RANGE OVERVIEW

	GOLD				SILVER			BRONZE				
	RAPHAEL QUATTRO LIFT RECLINE CHAIR				DA VINCI QUATTRO LIFT RECLINE CHAIR	SIGNATURE SPACE SAVER LIFT RECLINE CHAIR		SIGNATURE DUAL ACTION LIFT RECLINE CHAIR	PORTO ZERO GRAVITY PETITE LIFT RECLINE CHAIR	VERONA SPACE SAVER LIFT RECLINE CHAIR	MAUI DUAL ACTION LIFT RECLINE CHAIR	MONTANA SINGLE ACTION MAXI LIFT RECLINE CHAIR
												
Description	4 Motor - Quattro				4 Motor - Quattro	1 Motor - Space Saver		2 Motor - Dual Action	Zero Gravity	1 Motor - Space Saver	2 Motor - Dual Action	1 Motor - Heavy Duty
Functions												
Size	Small	Medium	Large	Maxi	Medium	Small	Medium	Large	Petite	Medium	Medium	Maxi
Upholstery	Microsuede Leather-look				Microsuede	Warwick Global Suede & Vinyl			Soft Suede	Quilted Upholstery	Quilted Upholstery	Microsuede Leather-look
Seat Height	470 mm	470 mm	500 mm	510 mm	480 mm	480 mm	480 mm	500 mm	490 mm	500 mm	500 mm	510 mm
Seat Depth	480 mm	500 mm	530 mm	540 mm	500 mm	460 - 485 mm	490 - 515 mm	510 - 535 mm	500 mm	500 mm	500 mm	540 mm
Seat Width	480 mm	520 mm	550 mm	650 mm	530 mm	480 mm	510 mm	540 mm	490 mm	540 mm	540 mm	650 mm
Overall Height	1060 mm	1090 mm	1080 mm	1050 mm	1050 mm	980 mm	1050 mm	1110 mm	1120 mm	1100 mm	1100 mm	1050 mm
Overall Width	740 mm	780 mm	810 mm	930 mm	870 mm	740 mm	780 mm	780 mm	770 mm	780 mm	780 mm	930 mm
Backrest Height	630 mm	680 mm	700 mm	650 mm	630 mm	580 mm	620 mm	650 mm	690 mm	630 mm	630 mm	650 mm
SWL	130 kg			220 kg	159 kg	145 kg	145 kg	155 kg	150 kg	130 kg	130 kg	220 kg
Warranty	2 Years				2-5 Years *	2 Years			2 Years	2 Years	2 Years	2 Years
AUST R ARTG No.	306075				216484	230186			306075	306075	306075	306075

## COLOUR RANGE



**Bronze**  
S - CHP227290  
M - CHP227300  
L - CHP227310  
MAXI - CHP227320




**Latte**  
CHP227500



**Fog**  
CHP227520



**Dove**  
CHP227530




**Gunmetal**  
CHP227510


### SIGNATURE SPACE SAVER



**Mink**  
S - CHP197700  
M - CHP197710  
L - CHP197720



**Lagoon**  
S - CHP197705  
M - CHP197715  
L - CHP197725



**Heritage - Vinyl**  
M - CHP197716  
L - CHP197730

### SIGNATURE DUAL ACTION

**Mink**  
S - CHP197735  
M - CHP197745  
L - CHP197755

**Lagoon**  
S - CHP197740  
M - CHP197750  
L - CHP197760

**Heritage - Vinyl**  
M - CHP197751  
L - CHP197765



**Fawn**  
CHP227200



**Mushroom**  
CHP227100



**Mocha**  
CHP227110



**Bronze**  
CHP227400

### ASPIRE ACCESSORIES

#### Chair Protector

Universal Size  
CHP227600





#### Signature Chairs - Head & Armrest Cover Set

**Mink**  
CHP197770

**Lagoon**  
CHP197771

**Heritage**  
CHP197772



\* 2 Years Handset/Upholstery & Electronics, 5 Years Frame