


MATTRESSES & CUSHIONS



Aspire LIFEComfort BARIATRIC Mattress

- Pressure reducing
- Safe working load of 300kg
- Four layer foam support for maximum support



Aspire Active Air Acute 8 Mattress Replacement

- Breathable multi-stretch cover
- Integrated cable management
- Simple easy-use digital pump
- Patented cable lock protection
- Welded seams & zip protection



ROHO Cushion

- Dry floatation promotes circulation and pressure relief
- Flexible air cells reduce shear friction
- Easily adjustable



ROHO Mattress Section

Interconnected and flexible air cells provide protection and relief for user at risk of skin breakdown or already with early stage ulcers



Gel Pad Cushion

Used to alleviate pressure build up & evenly distribute body weight while seated. Polymer gel does not displace like fluid and is covered in a waterproof film with breathable outer cover



J2 Plus Bariatric Cushion

- Multiple insert selections for customisation
- Anterior tilted base
- Cover choices including waterproof urethane overlay

TRANSFER AIDS



Aspire Slings - Large & Extra Large Sizes

- Versatile sling range offering General Purpose, Hygiene, Hammock & Standing slings
- Shaped for patient comfort & support
- Available in Polyester, Mesh & all day parasilk



AspireHD Lifter A320

- Lifts patients from the floor, or can accommodate transfers to higher surfaces
- 320kg user weight rating
- Electric leg gait and fold down armrests to assist gait training and toileting



Aspire A205 Lifter

- Compact and lightweight mobile patient lifter capable of 205kg patient transfers
- Ideal for transfers from bed to seating and toilet in the home, community or facility environments

AIDACARE BRANCHES

- Brunswick
- Brunswick East
- Cheltenham
- Geelong
- Bendigo
- Mildura
- Warrongong
- Moorebank
- Seven Hills
- Blacktown
- Central Coast
- Albury
- Orange
- Tamworth
- Port Macquarie
- Coffs Harbour
- Ballina
- Acacia Ridge
- Geebung
- Gold Coast
- Toowoomba
- Sunshine Coast
- Hervey Bay
- Hobart
- Launceston
- Regency Park
- Plympton Park



You may also be interested in Aidacare's other Product Application Guides giving product range information for the following areas:



bariatric XXL range

Equipment Selection Made Easy

1300 133 120
www.aidacare.com.au

DVA Aidacare are a DVA contractor. Ask us about DVA funded products.



bariatric XXL range

At Aidacare we've seen an increased need for equipment with a heavier weight-bearing capacity.

We've put together this guide to help our customers and their carers understand what is available from Aidacare's XXL Bariatric range.

A client is deemed to be Bariatric when they weigh over 125kg and/or require a seat width greater than 50cm. This means they need to be careful when selecting medical equipment to ensure that it is appropriate and safe.

Aidacare's range of XXL Bariatric equipment is engineered to be more robust than our standard range. This increased build-quality results in more durable products that help ensure the safety and independence of our customers.

Aidacare's wide range of XXL Bariatric equipment provides our customers with an extensive choice of brands and styles at an affordable price.

Our friendly staff and after-sales service support mean you can rest assured that you'll be taken care of.

Do you need help with product selection and sizes? Please contact Aidacare on 1300 133 120.

For Therapist Solutions please enquire about how else we can help you.

Rental

Did you know that Aidacare offers a comprehensive range of equipment for rent? We've got you covered should you need any products for home care and recovery!



Gutter Frame Walker – HD with brakes

- Encourages an upright gait
- Folding for transportation
- Powder coated mild steel
- Height adjustable

SWL 180KG



Aspire XL HEAVY DUTY Walker

- Available in silver
- Rustproof, robust aluminium frame for durability & stability
- Easy to manoeuvre
- Wide padded seat, 180 kg SWL

SWL 180KG



Aspire Evoke 2 HD Wheelchair

- Lightweight, highly adjustable manual wheelchair with swing away and removable arm rests & foot rests
- Quick release 24" solid rear tyres

SWL 180KG



Super Heavy Duty Wheelchair

- Reinforced wheelchair for bariatric users
- Extra wide seat widths and reinforced construction accommodates larger users

SWL 400KG



Forearm Crutches Custom

- Sturdy reinforced aluminium construction
- Absence of adjustment holes ensures "silent" use
- Forearm cuff for secure fitting

SWL 190KG



Underarm Crutches Bariatric Adult

- Heavy duty crutches to support the bariatric user
- Steel construction
- High density rubber cushions

SWL 249KG



Aspire Bariatric Shower Commode

- Super heavy duty 790mm wide, Anti-bacterial frame
- Height adjustable, swing away footrests
- Padded seat & arms
- 4x lockable castors

SWL 350KG



Aspire MAXI Shower Chair

- Super heavy duty width adjustable frame to 775mm
- Soft padding on seat and arms
- Reinforced aluminium frame is corrosion resistant
- Height adjustable legs

SWL 310KG



Aspire MAXI Over Toilet Aid

- Super heavy duty width adjustable frame to 775mm
- Soft padding on seat and arms
- Reinforced aluminium frame is corrosion resistant
- Height adjustable legs

SWL 310KG



Aspire MAXI Shower Stool

- Super heavy duty width adjustable frame to 775mm
- Soft padding on seat and arms
- Reinforced aluminium frame is corrosion resistant
- Height adjustable legs

SWL 310KG



Bath Transfer Bench Bariatric

- Safe working load of 227kg
- Fully height adjustable with side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs

SWL 227KG



Shower Stool Maxi Compact

- Height adjustable seat
- Reinforced for the larger user
- Lightweight aluminium and compact in design

SWL 227KG



Aspire Montana MAXI Lift Recline Chair

- Single motor lift recline chair
- Head & armrest covers for improved longevity
- Wipe clean upholstery
- Microsuede Leather-look

SWL 220KG



Aspire SIGNATURE Chair

- Recline & tilt functions for comfort and to aid in standing
- Independent leg elevation
- Depth adjustable seat

SWL 155KG



Aspire DELUXE Bedside Commode

- Large bedside chair can be used for toileting
- Safe working load of 160kg
- Fold back padded seat for easy access to commode bowl

SWL 160KG



Aspire AC3 KING SINGLE Bed

- 2 or 3 function electric lifting & reclining bed to increase independence & comfort
- Safe working load of 220kg

SWL 220KG



Aspire MAXI Day Chair

- Height adjustable backrest and legs - easy push button adjustment and non-slip rubber tips for stability
- Seat with anterior tilt

SWL 310KG



Sorrento Chair

- Maximise comfort, postural support and pressure management.
- Highly adjustable, with Tilt in space to redistribute pressure
- Back angle recline

SWL 294KG