

Product Awareness Bulletin – April 2018

PRODUCT UPDATES: Aspire Swift Patient Mover

We are pleased to introduce to you the Aspire Swift Patient Mover which has been in high demand since its market release at the end of 2017.

Key Features

- Soft knee & seat pads – each pad is single moulded PU structure, therefore completely sealed and easily wiped clean.
- Clip sling attachment points which allow for Aspire Standing Slings (with clips) to be attached for extra support
- Width adjustable legs for closer access (to toilets, seats, beds etc) are operated via two foot pedals
- Soft touch grab handle assists the user to move independently from sitting to standing


Additional Resources:

Promo & Instructional Video



Click [here](#) to view the video

Aidacare Website



Click [here](#) to view the product online

Product Brochure



Click [here](#) to view the brochure