

16:20



Hospital Equipment Range



Range Overviews

01	Hospital Beds
02	Bedroom Furniture & Seating
03	Support Surfaces
04	Patient Trolleys
05	Patient Hygiene
06	Patient Transport & Mobility
07	Mobile Pressure Seating
08	CSSD & Pan Room
09	Maternity & Birthing
10	Patient Handling - Mobile
11	Patient Handling - Ceiling Systems
12	Trolleys & Carts
13	Treatment Chairs
14	Examination Couches & Equipment
15	Bariatric

Aidacare Locations

VICTORIA

-  **MELBOURNE NORTH (1)**
35 Tinning St,
Brunswick VIC 3056
Tel: (03) 9981 2100
Fax: (03) 9386 9170
-  **MELBOURNE NORTH (2)**
82 Albion St,
Brunswick East VIC 3057
Tel: (03) 9981 2100
Fax: (03) 9386 9170
-  **MELBOURNE SOUTH**
313 Warrigal Rd,
Cheltenham VIC 3192
Tel: (03) 9584 1644
Fax: (03) 9584 0562
-  **GEELONG**
74-76 Fyans St,
Geelong VIC 3220
Tel: (03) 5223 2263
Fax: (03) 5223 2269
-  **BENDIGO**
51 Bridge St,
Bendigo VIC 3550
Tel: (03) 5441 8213
Fax: (03) 5441 7762
-  **MILDURA**
18 Eighth St,
Mildura VIC 3500
Tel: (03) 5022 0755
Fax: (03) 5021 1845

NEW SOUTH WALES

-  **SYDNEY SOUTH**
4 Noonan Pl,
Ingleburn NSW 2565
Tel: (02) 8706 2300
Fax: (02) 9618 5111
-  **SYDNEY NORTH**
Unit 2/4 Stanton Rd,
Seven Hills NSW 2147
Tel: (02) 8706 2330
Fax: (02) 9624 4256
-  **WOLLONGONG**
140 Church St,
Wollongong NSW 2500
Tel: (02) 8706 2359
Fax: (02) 9618 5111
-  **ALBURY**
347 Urana Rd,
Lavington NSW 2641
Tel: (02) 6040 3444
Fax: (02) 6040 3222
-  **ORANGE**
Unit 1/12-16 Ash St,
Orange NSW 2800
Tel: 0437 134 295
-  **TAMWORTH**
Unit 3/248 Peel St,
Tamworth NSW 2340
Tel: (02) 6761 2018
Fax: (02) 6761 2090

 **PORT MACQUARIE**
31 Jindalee Rd,
Port Macquarie NSW 2444
Tel: (02) 6581 0018
Fax: (02) 6581 0038

 **COFFS HARBOUR**
2/177 Orlando St,
Coffs Harbour NSW 2450
Tel: (02) 6651 5333
Fax: (02) 6652 1272

 **BALLINA**
18 Endeavour Cl,
Ballina NSW 2478
Tel: (02) 6686 6644
Fax: (02) 6686 9383

QUEENSLAND

 **BRISBANE SOUTH**
30 Dulacca St,
Acacia Ridge QLD 4110
Tel: (07) 3086 2900
Fax: (07) 3345 7961

 **BRISBANE NORTH**
354 Bilsen Rd,
Geebung QLD 4034
Tel: (07) 3086 2950
Fax: (07) 3865 7896

 **GOLD COAST**
Unit 1/48 Olympic Circuit,
Southport QLD 4215
Tel: (07) 5532 6070
Fax: (07) 5503 1803

 **TOOWOOMBA**
1/259 James St,
Toowoomba QLD 4350
Tel: (07) 3086 2990
Fax: (07) 4632 6473

 **SUNSHINE COAST**
2/35 Premier Circuit,
Warana QLD 4575
Tel: (07) 3086 2980
Fax: (07) 3345 7961

 **HERVEY BAY**
Shop 3/109 Beach Rd,
Hervey Bay QLD 4655
Tel: (07) 4194 5019
Fax: (07) 3345 7961

SOUTH AUSTRALIA

 **ADELAIDE**
21-23 Tikalara St
Regency Park SA 5010
Tel: (08) 8118 5600
Fax: (08) 8346 3236


TASMANIA

 **HOBART**
18 Mertonvale Circuit,
Kingston TAS 7050
Tel: (03) 6229 4733
Fax: (03) 6229 3676

 **LAUNCESTON**
Unit 8/34 Innocent St,
King Meadows TAS 7250
Tel: (03) 6334 8922
Fax: (03) 6344 7376

Aidacare Branches

-  Brunswick
-  Port Macquarie
-  Brunswick East
-  Coffs Harbour
-  Cheltenham
-  Ballina
-  Geelong
-  Acacia Ridge
-  Bendigo
-  Geebung
-  Mildura
-  Gold Coast
-  Ingleburn
-  Toowoomba
-  Wollongong
-  Sunshine Coast
-  Seven Hills
-  Hervey Bay
-  Albury
-  Adelaide
-  Orange
-  Hobart
-  Tamworth
-  Launceston





VIVO ICU Bed

The VIVO bed was designed with patient and care staff welfare in mind. It offers integrated features to ease care staff burden whilst creating a comfortable and effective patient support surface



Lateral Tilt
Programmable for individual patient needs between +/- 15°

- Lateral tilting patient surface
- Alternating air mattress with full bed integration to maximise therapeutic benefit
- Patient Anti Slipping System (PAS) offsets the patient movement towards the foot end of the bed during backrest adjustment
- Radio transparent patient surface for fluoroscopic examinations
- Touch screen panel in the foot end allows for the bed and integrated air mattress functions to be controlled with ease
- Integrated bed extension to suit a range of patient heights
- Intuitive control panels within the gas assisted folding side rails
- Safety alarm system ensures the base of the bed is clear before any movements, preventing entrapment and damage



BIOCOTE
Antibacterial frame coating for improved infection control

Integrated Nurse Control Panel
Provides fast and effective positioning of the patient support surface



Mattress
Programmable longitudinal air cell in backrest section for improved pulmonary ventilation. Three integrated pressure sensors to ensure maximum accuracy in pressure adjustment for the individual patient



Code	BEB050100
Overall Length	2190 - 2420mm
Overall Width	1080mm
Overall Height	480 - 820mm
Patient Surface Length	2210mm
Patient Surface Width	910mm
Backrest Angle	0° to 68°
Knee Bend Angle	0° to 25°
Trendelenburg	-15° to 15°
Castor Size	150mm
SWL	260kg




Logica Ward Bed

Offering leading functionality, LOGICA Ward Beds are versatile enough for any healthcare environment. The focus on design aesthetics provides every patient room with a comfortable, less clinical atmosphere to promote recovery.



Height Adjustable Side Rails



Low Height for Patient Safety

- Feature rich with low minimum height to meet evolving patient handling procedures
- Modern styling to compliment interior design trends
- Compact footprint and lightweight aluminium frame design
- 150mm Tente castors
- Central locking brake system with integrated directional lock
- X-Ray compatible mattress deck is available in ABS or compact laminate
- Nurse control panel with security drawer for efficient provision of care
- Integrated bed extension
- Full Trendelenburg functionality
- Single touch programming for common functional requirements such as Cardiac Chair and Fowler positioning
- Gas assisted side rails
- Removable Head and Foot Panels

Standard Inclusions

IV Pole, Catheter Hooks and Security Drawer for Nurse Control Panel



Adjustable to Facilitate Safe Patient Repositioning



BACK RAISE



LEG RAISE




BACK/LEG RAISE




KNEE BEND



HEIGHT ADJUSTABLE



TRENDELBURG



REVERSE TRENDELBURG




Code	BEB050110
Overall Length	2175mm
Overall Width	1021mm
Overall Height	375 - 835mm
Patient Surface Length	2000mm
Patient Surface Width	900mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 25°
Trendelenburg	-18° to 18°
Castor Size	150mm
SWL	250kg




Innova Ward Bed


The award winning INNOVA Ward Bed offers a fusion of function and form suitable for the modern healthcare environment. The focus on design aesthetics provides every patient room with a comfortable, less clinical atmosphere to promote recovery.



Adjustable to Facilitate Safe Patient Repositioning



Height Adjustable Drop Down Side Rails



Low Height for Patient Safety

- Feature rich with low minimum height to meet evolving patient handling procedures
- Modern styling to compliment interior design trends
- Compact footprint and lightweight aluminium frame design
- 125mm Tente castors
- Central locking brake system with integrated directional lock
- X-Ray compatible mattress deck is available in ABS or compact laminate
- Nurse control panel with security drawer for efficient provision of care
- Integrated bed extension
- Full Trendelenburg functionality
- Single touch programming for common functional requirements such as Cardiac Chair and Fowler positioning
- Gas assisted side rails
- Removable Head and Foot Panels
- Compatible with many leading bed movers

Standard Inclusions

IV Pole, Catheter Hooks and Security Drawer for Nurse Control Panel



reddot design award winner 2015



BACK RAISE



LEG RAISE



BACK/LEG RAISE



KNEE BEND



HEIGHT ADJUSTABLE



TRENDELENBURG



REVERSE TRENDELENBURG



CARDIAC POSITION



CHAIR POSITION



Code	BEB050120
Overall Length	2230mm
Overall Width	1070mm
Overall Height	365 - 860mm
Patient Surface Length	2000mm
Patient Surface Width	900mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 25°
Trendelenburg	-18° to 18°
Castor Size	125mm
SWL	250kg



Botero Bariatric Ward Bed

The specific bariatric design focus offers width and length adjustable deck to manage changing needs in both ward and respite environments.



Integrated Bed Length Extension
180mm of extension to cater for taller patients



Attendant Control Panel
Single Touch Functions for Speed and Simplicity

- Low height of 380mm
- Antibacterial powder coating for improved infection control management
- Backrest translation (autoregression) reduces compression of the thoracic and abdominal areas as the backrest is raised and alleviates shear and pressure for improved skin integrity
- Compatible with all leading bed movers
- Attendant control panel with intuitive pre-programmed single touch functions for speed and simplicity
- 150mm Tente Castors
- Central locking castors with integrated directional lock
- Full Trendelenburg functionality
- Integrated folding side rails
- Removable width adjustable head and foot panels



S.W.L
400 kg



BIOCOTE
Antibacterial frame coating for improved infection control

Integrated Bed Width Extension
1020mm - 1220mm of tool free extension for changing needs

Optional Accessories



Static Mattress w/ Extension Bolster (BEM050200)



Self Help Pole (BEA825100)



Telescopic IV Pole (BEA825101)



Code	BEB050200
Overall Length	2260 - 2440mm
Overall Width	1020 - 1220mm
Overall Height	380 - 850mm
Patient Surface Length	2000 - 2180mm
Patient Surface Width	1020 - 1220mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 30°
Trendelenburg	-17° to 17°
Castor Size	150mm
SWL	400kg



Aidacare Ward Bed

The robust yet simple design offers a strong value proposition for general ward environments.



Twin 150mm Tente Castors
Central and Directional Locking Functionality



Integrated Bed Length Extension
200mm of extension to cater for taller patients

- Available with ¾ length folding or gas assisted soft drop side rails
- 150mm Tente Castors
- Central Locking with Directional lock functionality
- Attendant Control Panel with function lock-out capability
- Integrated 200mm bed extension
- Easy to clean mattress platform with removable panels
- Removable head and foot board panels to meet patient handling and treatment requirements
- LINAK electronics with battery backup



Soft drop side rails ordered separately (BEA029810)



Attendant Control Panel
Pre Programmed Positioning for Cardiac Chair and CPR Functions

Optional Accessories



Deluxe Soft Drop Side Rails (BEA029810)



Metal Folding Side Rails (BEA029700)



Self Help Pole (BEA029705)



Telescopic IV Pole (BEA029710)

Code	BEB046000
Overall Length	2195 - 2390mm
Overall Width	1020mm
Overall Height	355 - 815mm
Patient Surface Length	2000mm
Patient Surface Width	900mm
Backrest Angle	0° to 77°
Knee Bend Angle	0° to 35°
Trendelenburg	0° to 14°
Castor Size	150mm
SWL	250kg






S.W.L
150kg

Happy Recliner Chair

The Happy Recliner Chair offers a comfortable seating surface in-room or in a mobile configuration

- Intuitive recline mechanism for patient independence
- Three comfort positions
- Adjustable padded neck pillow for postural support
- Mobility castor option for patient transportation
- Comfortable padded arms with tapered front upright for improved table access
- Lightweight and easy to move promoting workplace health and safety



Beige Vinyl

Code	CHP198820
Overall Height	1080mm
Overall Width	690mm
Overall Depth	780mm
Seat Height	420mm
Seat Width	520mm
Seat Depth	580mm
Backrest Height	660mm
S.W.L	150kg




S.W.L
150kg

Happy Sleeper Chair

The Happy Sleeper Chair facilitates the important role that family and friends play in the recovery process

- Durable seat construction preserves shape and comfort in a high use environment
- Intuitive single hand activation allows effortless conversion from chair to a comfortable sleep surface
- Recessed locking castors offer easy positioning and accessibility
- Completely flat sleep surface to maximise rest and relaxation
- Contemporary design compliments a range of care environments



Overall Height	780mm
Overall Width	760mm
Overall Depth	960mm
Seat Height	370mm
Seat Width	600mm
Seat Depth	580mm
Backrest Height	410mm
S.W.L	150kg



Beige Vinyl
(CHP198810)



Teal Vinyl
(CHP198812)



Light Blue Vinyl
(CHP198814)



Navy Blue Vinyl
(CHP198816)



S.W.L
150kg

Happy Dialysis Chair

The Happy Dialysis Chair offers a comfortable seating surface in-room or in a mobile configuration for infusion therapy

- Tapered fold down armrests promote comfort and accessibility during infusion therapy
- Intuitive recline mechanism for patient independence
- Three comfort positions
- IV pole mounted to rear of chair
- Mobility castor option for patient transportation
- Lightweight and easy to move promoting workplace health and safety



Chocolate
Brown Vinyl


Code	CHT100400
Overall Height	1080mm
Overall Width	690mm
Overall Depth	780mm
Seat Height	420mm
Seat Width	520mm
Seat Depth	580mm
Backrest Height	660mm
S.W.L	150kg



S.W.L
160kg

High Back Classic Day Chair

- Contoured backrest for postural support
- Softly padded armrests
- Lightweight yet sturdy design
- Powder coated steel frame
- Easy clean vinyl surfaces
- Height adjustable legs - Easy push button
- Non-slip tips provide stability
- Adjustable chair legrest (optional extra)



Overall Height	1155 - 1315mm
Overall Width	645mm
Overall Depth	600mm
Seat Depth	490mm
Seat Width	490mm
Seat Height	425 - 610mm
Backrest Height	720mm
S.W.L	160kg



S.W.L
200kg

Waterfall Day Chair

- Waterfall backrest provides comfort and support
- Removable cushion sections
- Padded seat and armrests
- Powder coated steel frame
- Non-slip tips provide stability
- Height adjustable
- Sturdy design - high safe working load of 200kg, yet lightweight for easy manoeuvrability
- Available in either vinyl (Ink or Latte) or fabric (Oatmeal or Mocha)




Latte Vinyl
(CHP208071)




Ink Vinyl
(CHP208070)



Overall Height	1045 - 1200mm
Overall Width	720mm
Overall Depth	720mm
Seat Depth	490mm
Seat Width	535mm
Seat Height	420 - 580mm
Backrest Height	685mm
S.W.L	200kg



Mocha Fabric
(CHP208076)



Oatmeal Fabric
(CHP208075)



S.W.L
160kg

Adjustable Day Chair

- Depth adjustable up to 100mm
- Width adjustable up to 100mm
- Height adjustable up to 120mm
- Easy clean vinyl upholstery
- Push handle and transport wheels
- Sturdy design - Generous 160kg weight capacity
- Height adjustable and removable head cushion for extra comfort



Overall Height	1070 - 1220mm
Overall Width	575 - 675mm
Overall Depth	760 - 860mm
Seat Depth	420 - 620mm
Seat Width	510 - 610mm
Seat Height	410 - 530mm
Backrest Height	740 - 790mm
S.W.L	160kg



Latte Vinyl
(CHP208025)



Ink Vinyl
(CHP208030)



S.W.L
160kg

Low Back Classic Day Chair

- Contoured backrest for postural support
- Softly padded armrests
- Lightweight yet sturdy design
- Powder coated steel frame
- Easy clean vinyl surfaces
- Height adjustable legs
- Easy push button
- Non-slip tips provide stability
- Adjustable chair legrest (optional extra)



Overall Height	845 - 1015mm
Overall Width	580mm
Overall Depth	600mm
Seat Depth	440mm
Seat Width	460mm
Seat Height	430 - 605mm
Backrest Height	440mm
S.W.L	160kg



Champagne Vinyl
(CHP208000)




Slate Vinyl
(CHP208005)



S.W.L
125kg

Classic Bedside Commode

- Raised padded armrests
- Removable bowl with lid and carry handle
- Easy clean vinyl surfaces
- Removable, padded seat
- Height adjustable
- Three-sided shroud – discreetly conceals bowl
- Non-slip rubber tips



Overall Height	900 - 1040mm
Overall Width	560mm
Overall Depth	600mm
Seat Depth	440mm
Seat Width	450mm
Seat Height	475 - 615mm
Backrest Height	450mm
S.W.L	125kg



Champagne Vinyl
(BEC045005)



Slate Vinyl
(BEC045006)



S.W.L
160kg

Deluxe Bedside Commode

- Deeply padded comfortable seat, back and armrests
- Removable bowl with lid and carry handle
- Easy clean vinyl surfaces
- Lift-up padded seat
- Height adjustable
- Three-sided shroud – discreetly conceals bowl
- Non-slip rubber tips



Overall Height	865 - 1015mm
Overall Width	600mm
Overall Depth	600mm
Seat Depth	470mm
Seat Width	460mm
Seat Height	410 - 565mm
Backrest Height	440mm
S.W.L	160kg



Champagne Vinyl
(BEC045115)




Slate Vinyl
(BEC045116)



S.W.L
175kg

Perching Stool

- Designed to facilitate regular sitting and standing assistance
- Generously padded seat and back support surfaces
- Raised armrests – Assist the user transferring on and off the stool
- Seat with anterior tilt - Facilitates standing from sitting and provides greater comfort for stiff or painful hips
- Height adjustable legs with easy push button adjustment and non-slip rubber tips for stability



Overall Height	950 - 1100mm
Overall Width	545mm
Overall Depth	400mm
Seat Depth	310mm
Seat Width	375mm
Seat Height	540 - 690mm
Between Arms	495mm
S.W.L	175kg



Champagne Vinyl
(CHP213500)

Transport Legs With Wheels



Code	Description
CHP208090	Transport Legs with Wheels - Aspire High Back, Waterfall & Adjustable Day Chairs

Footrest



Code	Description
CHP208080	Footrest - Champagne Vinyl
CHP208085	Footrest - Slate Vinyl



Standard Option 1
Beige Frame, White Top
(BEA006700)

Laminate & Custom Overbed Tables

- Stable and secure table top to conveniently accommodate meals, drinks, reading materials, tablet devices and other items within easy arms reach from a bed or chair.
- Height adjustment lever is easily accessible to users and carers
- Height fixing knob will lock the height and prevent adjustment by users
- Two locking castors for securing table position by the carer to prevent it rolling away

Standard Option 2
Beige Frame, Beech Top
(BEA006701)



Height Range	730 - 1095mm
Tabletop Width	430mm
Tabletop Length	800mm

Custom Options

Frames	Table Tops	
Gunmetal Grey (BEA006685)	White (BEA006693)	Sublime Teak (BEA006695)
Beige (BEA006690)	Milano Walnut (BEA006696)	Beech (BEA006694)


* We also do bespoke tabletops to order. Please call to discuss your requirements.



Recessed Top Overbed Table

- Recessed top collects and prevents liquids from spilling over the edge
- Thermoformed surface is an impermeable layer that continues over and down the vertical edges so it is easily wiped clean.
- Height adjustment lever is easily accessible to users and carers
- Height fixing knob can be fastened by carer will lock the height and prevent adjustment by users
- Two locking castors for securing table position by the carer to prevent it rolling away in use.

Height Range	730 - 1095mm
Tabletop Width	430mm
Tabletop Length	800mm




Beige Frame, Recessed Cream Top (BEA028000)



Recessed Split Top Overbed Table

- Top can be tilted and fixed in position to suit the user, while the smaller surface remains horizontal
- Recessed top collects and prevents liquids from spilling over the edge
- Thermoformed surface is an impermeable layer that continues over and down the vertical edges so it is easily wiped clean
- Height adjustment lever, height fixing knob and two locking castors for comfort and securing

Height Range	730 - 1095mm
Tabletop Width	430mm
Tabletop Length	800mm
Tilt Section Length	600mm
Fixed Section Length	200mm



Beige Frame, Recessed Split Cream Top (BEA028100)





Premium Overbed Table



- Familiar "C-Shape" base for chair accessibility
- Smooth gas-assisted movement
- Anti-crush mechanism automatically raises table top if bed is raised up underneath
- High density laminate timber top
- Locking castor wheels provide a stable base of support

Code	BEA006710
Height Range	805 - 1050mm
Tabletop Width	400mm
Tabletop Length	800mm



Premium Split Top Overbed Table



- Split top functionality offers a stable drink support surface whilst facilitating reading
- Familiar "C-Shape" base for chair accessibility
- Smooth gas-assisted movement
- Anti-crush mechanism automatically raises table top if bed is raised up underneath
- Edge barriers help prevent accidents
- Locking castor wheels provide a stable base of support
- Document holder accessory

Code	BEA006720
Locker Height	795 - 1040mm
Locker Width	400mm
Locker Length	800mm





Bedside Cabinets

- Available in two different surface finishes, Laminate and Compact
- Laminate finish designed to match modern room decor
- Compact surface is highly durable and robust, catering to acute environments
- Easy to clean for improved infection control
- Four castors for easy positioning and accessibility
- Aluminium corner uprights minimise wear over time (Compact Models)
- Colour contrasting, easy to see handles
- Soft closing drawer runners



1 Drawer + 1 Door
- Laminate



BEA827100

Overall Height	650mm
Overall Width	455mm
Overall Depth	430mm



3 Drawer - Laminate



BEA827150

Overall Height	615mm
Overall Width	600mm
Overall Depth	445mm



2 Drawer - Laminate



BEA827130

Overall Height	450mm
Overall Width	450mm
Overall Depth	450mm



3 Drawer - Laminate



BEA827140

Overall Height	580mm
Overall Width	500mm
Overall Depth	400mm



2 Drawers + 1 Shelf
- Compact



BEA827110

Overall Height	745mm
Overall Width	500mm
Overall Depth	530mm



2 Drawer - Compact



BEA827120

Overall Height	755mm
Overall Width	500mm
Overall Depth	530mm




1 Drawer + 1 Door
- Compact



BEA827090


Overall Height	735mm
Overall Width	530mm
Overall Depth	505mm



Acute Care Static Mattress

- Castellated surface redistributes shear forces whilst increasing airflow for improved microclimate
- Multi-stretch cover with welded seams for stringent infection control
- Side handles allow for easy mattress rotation and movement
- Strong high density base and side walls improve patient transfer efficiency


Code	BEM046615	BEM046616
	(Single)	(King Single)
Overall Length	1980mm	2030mm
Overall Width	880mm	1050mm
Overall Height	150mm	150mm
S.W.L	250kg	250kg
Warranty	5 years foam / 2 years cover	



Bariatric Static Mattress

- Four layer foam design caters for bariatric users
- Multi-stretch cover with welded seams, reinforced waterfall flaps and water resistant zippers for stringent infection control
- Side handles allow for easy mattress movement
- Strong high density base and side walls improve patient transfer efficiency


Code	BEM046700
Overall Length	2030mm
Overall Width	1050mm
Overall Height	200mm
S.W.L	300kg
Warranty	5 years foam / 2 years cover



Botero Bariatric Static Mattress

- Two layer foam design with high density base and castellated surface caters to bariatric patients
- Multi-stretch cover with welded seams for stringent infection control
- Strong high density side walls improve patient transfer efficiency while cut back hinge points allow the mattress to easily move with the mattress platform

Code	BEM050200
Overall Length	2000mm
Overall Width	1220mm
Overall Height	200mm
S.W.L	400kg
Warranty	5 years foam / 2 years cover







S.W.L
230kg

Hybrid Acute Mattress

- Foam-in-cell technology enables flexible Hybrid operation as a static support surface or powered-up alternating pressure support surface.
- The memory foam overlay moves independently of the underlying cells, reducing shear as the bed profiles
- A 12° tapered heel zone redistributes pressure up the lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear
- When not in use the air-hose can be stored internally, reducing infection control and OH&S risks



Code	PTM638020 (Single)	PTM638025 (King Single)
Overall Length	1980mm	1980mm
Overall Width	895mm	1050mm
Overall Height	170mm	170mm
Number of Cells	9 (incl 3x static head cells)	9 (incl 3x static head cells)
S.W.L	230kg	230kg
Digital Pump	PTM638030	PTM638030
Warranty	2 years pump / 1 year mattress	



Additional pump for an active support surface



S.W.L
135kg

S.W.L
175kg

Active Air 4 & 5 Mattress Overlay

- 4 or 5 inch mattress overlay with digital pump suitable for patients at medium to at risk of pressure injury
- Multi-stretch cover with high Moisture Vapour Transmission Rate (MVTR)
- Welded seams, reinforced waterfall flaps and water resistant zippers for stringent infection control
- Rapid deflation CPR valve for use in emergency situations
- Automatic weight detection ensures optimal care and protection



Overlay	Active Air 4	Active Air 5	Active Air 5
Code	PTM639100 (Single)	PTM639200 (Single)	PTM639210 (King Single)
Overall Length	2000mm	2000mm	2000mm
Overall Width	850mm	850mm	1050mm
Overall Height	100mm	130mm	130mm
Number of Cells	19 (incl 3x static head cells)	16 (incl 3x static head cells)	16 (incl 3x static head cells)
S.W.L	135kg	175kg	175kg
Pump	Semi-Digital	Digital	Digital
Warranty	2 years pump / 1 year mattress		

Optional Accessory



Rapid Inflation Turbo Pump (PTM639000)



S.W.L
200kg

Active Air 8 Mattress Replacement

- 8 inch mattress replacement with digital pump suitable for patients at high risk of pressure injury
- Multi-stretch cover with high Moisture Vapour Transmission Rate (MVTR)
- Welded seams, reinforced waterfall flaps and water resistant zippers for stringent infection control
- Rapid deflation CPR valve for use in emergency situations
- Automatic weight detection ensures optimal care and protection
- Automatic pressure increase into backrest when profiling bed to prevent bottoming out



Code	PTM639300 (Single)	PTM639310 (King Single)
Overall Length	2000mm	2000mm
Overall Width	900mm	1050mm
Overall Height	200mm	200mm
Number of Cells	20 (incl 3x static head cells)	20 (incl 3x static head cells)
S.W.L	200kg	200kg
Pump	Digital	Digital
Warranty	2 years pump / 1 year mattress	



Digital Pump

Optional Accessory



Rapid Inflation Turbo Pump (PTM639000)



ONEday Electric Hospital Stretcher

- 4-section fully electric stretcher offering the functions of a bed with the smaller footprint and mobility of a stretcher
- Electric Trendelenburg and reverse-Trendelenburg functions to meet clinical and patient handling needs
- 200mm castors minimise rolling resistance
- Fifth central wheel improves steering and reduces manual handling burden in transportation
- 260kg SWL and inbuilt leg extension facilitate changing patient and clinical needs across multiple ward requirements
- Hand Control flexible arm promotes patient independence and comfort
- Enhanced battery-pack, with day-long operating time, dedicated battery charger for a full and fast recharge
- Audible alarm system sounds when mains power is lost and when brakes are unlocked
- Quick Release CPR for back rest and knee break ensures the stretcher is emergency ready at all times
- Frame is finished with Biocote, a technologically advanced surface coating that promotes infection control via its anti-bacterial properties



BIOCOTE
Antibacterial frame coating for improved infection control

Optional Accessories



Oxygen bottle holder (TSE674110)



Telescopic IV Pole (TSE674130)



Static Mattress w/Extension Bolster (BEM050300)



Tip-over monitor tray/writing desk (TSE674120)



Code	TSE674100
Overall Length	1990 - 2170mm
Overall Width	890mm
Overall Height	400 - 830mm
Patient Surface Length	1880 - 2060mm
Patient Surface Width	710mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 30°
Trendelenburg	-17° to 17°
Castor Size	200mm
SWL	260kg



S.W.L
225kg

Advanced Emergency Stretcher

- Foot operated hydraulic lift, tilt, back raise, Trendelenburg and reverse Trendelenburg functions
- Heavy Duty 225kg Safe Working Load
- Full length sliding X-ray deck ensures complete spinal coverage
- Comfortable 100mm foam mattress with multi-stretch cover
- Safety straps for staff and patient security
- CPR functions integrated on either side for rapid response time
- Collapsible vertical side rails fold away below the patient surface
- Fold-away carer push handles to maximise workspace
- 200mm castors minimise rolling resistance
- Fifth wheel improves steering and reduces manual handling burden in transportation
- Central and Directional Locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Large corner bumpers minimise wear and tear on equipment and facilities
- Intuitive design and function layout from Red Dot design winning team

Standard Inclusions: Mattress, IV Pole, Oxygen Bottle Holder, Utility Basket, Urine Bag Hooks, Accessory Slots



Code	TSE674210
Overall Length	2140mm
Overall Width	840mm
Overall Height	600 - 910mm
Patient Surface Length	1955mm
Patient Surface Width	600mm
Backrest Angle	0° to 90°
Trendelenburg	-12° to 12°
Castor Size	200mm
SWL	225kg



*Pictured without Mattress



S.W.L
175kg


Emergency Stretcher

- Foot operated hydraulic lift, tilt, back raise, Trendelenburg and reverse Trendelenburg functions
- Safe Working Load of 175kg
- Comfortable 100mm foam mattress with multi-stretch cover
- CPR functions integrated on either side for rapid response time
- Collapsible horizontal side rails fold away below the patient surface
- 200mm castors minimise rolling resistance
- Fifth wheel improves steering and reduces manual handling burden in transportation
- Central and Directional Locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Corner bumpers minimise wear and tear on equipment and facilities

Standard Inclusions: Mattress, IV Pole, Oxygen Bottle Holder, Urine Bag Hooks, Accessory Slots



Code	TSE674200
Overall Length	2140mm
Overall Width	840mm
Overall Height	590 - 900mm
Patient Surface Length	1920mm
Patient Surface Width	600mm
Backrest Angle	0° to 70°
Trendelenburg	-12° to 12°
Castor Size	200mm
SWL	175kg



Hydraulic Shower Trolley

- Comfortable, padded and waterproof patient surface with welded seams for infection control
- Internal patient surface 1950mm x 600mm
- Height adjustable from 500mm to 775mm making sideways slide transfers a simple process
- Wall buffers promote durability and limit facility wear
- Durable stainless steel frame
- Central locking mechanism with directional lock functionality
- Hospital grade 125mm castors
- Integrated hydraulic Trendelenburg tilt function aids patient positioning and water drainage

Code	BTS148100
Overall Length	2150mm
Overall Width	840mm
Overall Height	500 - 775mm
S.W.L	150kg

S.W.L
150kg

Electric Shower Trolley

- Epoxy powder coated & ergonomic push handle
- Liner can be secured in the corners to facilitate in bath
- Padded vinyl mattress & corner buffer discs
- Large screw in diameter 60mm flexible waste hose
- Lightweight swing under side rails
- Electric Trendelenburg/tilt 8° degrees
- Electric high-low operation 600-1000mm
- Curved base construction to allow access for carers
- 125mm diameter total frame lock twin wheel castors
- Plug-in rechargeable battery pack
- Stainless steel & P.V.C. Top



Code	BTS147100 (Standard)	BTS147105 (Wide Long)
Overall Length	1890mm	2000mm
Overall Width	730mm	900mm
Overall Height	600mm - 1000mm	600 - 1000mm
S.W.L	150kg	150kg

S.W.L
150kg




S.W.L
200kg

* Attendant Propelled Option

Weight Bearing Platform Commode - 460mm

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Easy-slide weight bearing platform footrest slides completely away for transfers
- *Self propelled option:* Anti-tippers for maximum safety & 2x rear lockable self propelling wheels
- *Attendant propelled option:* 4x total locking stainless steel castors provide brake access from all angles

Propelling Code	Attendant	Self
	BTC066000	BTC066200
Overall Height	995mm	995mm
Overall Width	535mm	665mm
Seat Width	460mm	460mm
Backrest Height	305mm	305mm
Armrest Height	225mm	225mm
Anterior Seat	590mm	590mm
Posterior Seat	560mm	560mm
S.W.L	200kg	200kg



Seats Sold Separately




S.W.L
200kg

* Attendant Propelled Option

Swing Away Legrest Commode - 460mm

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- *Self propelled option:* Anti-tippers for maximum safety & 2x rear lockable self propelling wheels
- *Attendant propelled option:* 4x total locking stainless steel castors provide brake access from all angles

Propelling Code	Attendant	Self
	BTC065900	BTC066100
Overall Height	995mm	995mm
Overall Width	535mm	665mm
Seat Width	460mm	460mm
Backrest Height	305mm	305mm
Armrest Height	225mm	225mm
Anterior Seat	590mm	590mm
Posterior Seat	560mm	560mm
S.W.L	200kg	200kg



Seats Sold Separately



S.W.L
200kg

* Attendant Propelled Option

Weight Bearing Platform Commode - 530mm

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Easy-slide weight bearing platform footrest slides completely away for transfers
- *Self propelled option*
- *Attendant propelled option:* 4x total locking stainless steel castors provide brake access from all angles

Propelling Code	Attendant	Self
	BTC066050	BTC066055
Overall Height	995mm	995mm
Overall Width	605mm	735mm
Seat Width	530mm	530mm
Backrest Height	305mm	305mm
Armrest Height	225mm	225mm
Anterior Seat	590mm	590mm
Posterior Seat	560mm	560mm
S.W.L	200kg	200kg



Seats Sold Separately



S.W.L
200kg

* Attendant Propelled Option

Swing Away Footrest Commode - 530mm

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- *Self propelled option available* -----
- *Attendant propelled option: 4x total locking stainless steel castors provide brake access from all angles*

Propelling Code	Attendant	Self
	BTC066060	BTC066065
Overall Height	995mm	995mm
Overall Width	605mm	735mm
Seat Width	530mm	530mm
Backrest Height	305mm	305mm
Armrest Height	225mm	225mm
Anterior Seat	590mm	590mm
Posterior Seat	560mm	560mm
S.W.L	200kg	200kg



Seats Sold Separately




S.W.L
400kg



Weight Bearing Platform Commode - 600mm

- Reinforced heavy duty platform commode for users of size
- Height adjustable, reinforced and padded armrests with safety ends drop flat and out of the way to assist with positioning or lateral slide transfers -----
- Open back straps provide access for bathing
- Easy-slide weight bearing platform footrest slides completely away for transfers
- 4x total locking stainless steel castors provide brake access from all angles

Propelling Code	Attendant
	BTC061400
Overall Height	995mm
Overall Width	725mm
Seat Width	600mm
Backrest Height	300mm
Armrest Height	290mm
Anterior Seat	590mm
Posterior Seat	560mm
S.W.L	400kg



Seats Sold Separately



S.W.L
200kg

* Attendant Propelled Option

Swing Away Footrest Commode - 600mm

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- *Self propelled option available* -----
- *Attendant propelled option: 4x total locking stainless steel castors provide brake access from all angles*


Propelling Code	Attendant	Self
	BTC061300	BTC061305
Overall Height	995mm	995mm
Overall Width	675mm	805mm
Seat Width	600mm	600mm
Backrest Height	305mm	305mm
Armrest Height	225mm	225mm
Anterior Seat	590mm	590mm
Posterior Seat	560mm	560mm
S.W.L	200kg	200kg



Seats Sold Separately



S.W.L
250kg



Seats Sold Separately


Tilt & Recline Shower Commode

- Backrest may be configured as full length or 3/4 length to accommodate an advanced headrest
- Open back straps provide access for bathing
- 30° backrest angle adjustment in 10° increments for effortless repositioning and maximum user comfort
- Safety arms prevent falls and swing back out of the way for transfers
- Open back straps provide access for bathing
- -5° to 35° gas-assisted tilt-in-space functionality minimises manual handling efforts

Code	BTC066030	BTC066040
Overall Height	1315mm -1485mm	1315mm -1485mm
Overall Width	560mm	630mm
Seat Width	460mm	530mm
Backrest Height	900mm	900mm
Armrest Height	275mm	275mm
Anterior Seat Height	460mm - 630mm	460mm - 630mm
Posterior Seat Height	460mm - 620mm	460mm - 620mm
S.W.L	250kg	250kg



S.W.L
350kg



Bariatric Shower Commode

- Super Heavy Duty Bariatric design with 790mm seat and 350kg SWL
- Padded seat and arms for maximum comfort
- Removable arms to assist with positioning and allow lateral slide transfers
- 4x lockable castors maximise safety by providing carer access to a brake at all angles around the commode
- Height adjustable, swing away footrests for optimal positioning and ease of transfers

Code	BTC061200
Overall Height	990mm
Overall Width	900mm
Seat Width	790mm
Backrest Height	430mm
Armrest Height	170mm
Anterior Seat	570mm
Posterior Seat	570mm
S.W.L	350kg





Standard Seats

Open Seat - 460mm	BTC066600
Closed Seat - 460mm	BTC066700
Open Seat - 530mm	BTC066750
Open Seat - 600mm	BTC066800



Hand Made Pressure Reducing Seats



Open Front - 320mm	BTC089122
Closed Front - 320mm	BTC089124
Open Front - 400mm	BTC089126
Closed Front - 400mm	BTC089128
Open Front - 460mm	BTC089130
Closed Front - 460mm	BTC089140
Open Left - 460mm	BTC089150
Open Right - 460mm	BTC089160
Open Front - 500mm	BTC089170
Closed Front - 500mm	BTC089180
Open Front - 550mm	BTC089190
Closed Front - 550mm	BTC089200
Open Front - 600mm	BTC089210
Closed Front - 600mm	BTC089220
Open Front - 650mm	BTC089230
Closed Front - 650mm	BTC089240
Custom Option	BTC089250



ROHO Padded Seats

380mm x 425mm	BTC067500
425mm x 425mm	BTC067900
425mm x 470mm	BTC067600
470mm x 470mm	BTC067300



Shower Commode Accessories



1. Self Propelling Conversion Kit



2. Elevating Leg Rest



3. Commode Bowl & Pan



4. Pan Carrier



5. Oxygen Bottle Carrier



6. Adjustable IV Pole



7. Adjustable Pelvic Belt



8. Advanced Headrest



9. Stump Support



10. Adjustable Thoracic Supports

1. Self Propelling Conversion Kit	BTC082900
2. Elevating Legrests - Left	BTC067110
2. Elevating Legrests - Right	BTC067115
3. Commode Bowl with Lid & Handle	BTC081610
3. Bed Pan	BTT151800
4. Pan Carrier	BTC066400
5. Oxygen Bottle Carrier	BTC067200
6. Adjustable IV Pole	BTC067103
7. Adjustable Pelvic Belt	BTC067100
8. Advanced Headrest	BTC067050
9. Stump Support	BTC067105
10. Adjustable Thoracic Supports	BTC067101




S.W.L
310kg

MAXI Shower Chair

- Super heavy duty adjustable width shower chair
- Soft padded seat, back and arms provide postural support
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the chair
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability



Code	BTS118020
Overall Height	670 - 1005mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg




S.W.L
310kg

MAXI Shower Stool

- Super heavy duty adjustable width shower stool
- Soft padded seat and arms provide postural support
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability



Code	BTS119030
Overall Height	670 - 790mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg




S.W.L
310kg

MAXI Over Toilet Aid

- Super heavy duty adjustable width over toilet aid
- Soft padded flip up seat for easy cleaning
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
- To create a MAXI bedside commode (BEC045120), simply add a padded backrest (as shown in MAXI Shower Chair above)



Code	BTT145020
Overall Height	670 - 790mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg






S.W.L
200kg

Shower Chair

- Contoured plastic seat with drainage holes
- Height adjustable legs
- Frames individually treated for corrosion resistance
- Raised armrests for assisted user transfer
- Non-slip grip and stability in wet conditions.
- Stainless steel washer within enhances durability
- Heavy duty, reinforced frame accommodates very high SWL




Code	Regular	Wide
Overall Height	735 - 910mm	735 - 910mm
Overall Width	512mm	570mm
Overall Depth	570mm	570mm
Seat Height	375 - 550mm	375 - 550mm
Seat Width	440mm	440mm
Seat Depth	400mm	400mm
Width Between Arms	470mm	525mm
S.W.L	200kg	200kg

Steel
(BTS118000)

Steel - Wide
(BTS118010)

Aluminium
(BTS118005)

Aluminium - Wide
(BTS118015)



S.W.L
175kg

Shower Stool

- Clip on seat or padded back attachment – For easy assembly and cleaning
- Height adjustable legs – Easy push button height adjustment for changing needs
- All frames are individually treated for corrosion resistance



Description	Regular	Wide
Overall Height	670 - 855mm	575 - 765mm
Overall Width	545mm	570mm
Overall Depth	400mm	500mm
Seat Height	470 - 655mm	425 - 615mm
Seat Width	375mm	420mm
Seat Depth	310mm	365mm
Width Between Arms	495mm	525mm
S.W.L	175kg	175kg

Steel - Regular
(BTS119000)

Aluminium - Regular
(BTS119005)

Steel with Padded Seat - Regular
(BTS119020)

Steel with Padded Seat and Back - Regular
(BTS119025)

Steel - Wide
(BTS119010)


Aluminium - Wide
(BTS119015)



S.W.L
175kg

Over Toilet Aid

- Raised armrests - Assist the user transferring on and off the over toilet aid
- Height adjustable legs
- Comfortable, contoured seat with clip on attachment
- Frames individually treated for corrosion resistance
- Removable, moulded splash guard provided
- Durable non-slip feet for grip and stability in wet conditions and a stainless steel washer within enhances durability



Description	Regular	Wide
Overall Height	700 - 880mm	700 - 880mm
Overall Width	512mm	570mm
Overall Depth	625mm	625mm
Seat Height	450 - 630mm	450 - 630mm
Seat Width	370mm	370mm
Seat Depth	435mm	435mm
Width Between Arms	460mm	525mm
S.W.L	175kg	175kg

Steel - Regular
(BTT145000)

Aluminium - Regular
(BTT145005)

Steel - Wide
(BTT145010)

Aluminium - Wide
(BTT145015)



Seats Sold Separately




Toilet Surround Frame



S.W.L
175kg

- Height adjustable frame – Easy push button height adjustment for changing needs
- Raised armrests – Assist the user transferring on and off the toilet
- Durable non-slip feet – Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Additional height increments on one leg – Accommodates uneven surfaces



Code	BTT155100
Overall Height	535 - 660mm
Overall Width	500mm
Overall Depth	435mm
Between Arms	450mm
Unit Weight	3kg
S.W.L	175kg




Raised Toilet Seat



S.W.L
200kg

- Three heights available – helps facilitate sit-to-stand transfers for safe and comfortable toileting
- Contoured plastic seat – Comfortable and durable
- Available with or without lid
- Fully sealed structure – Easy to clean for increased infection control
- Generous cut outs at the front and rear – Facilitates personal care and hygiene
- Two adjustable mounting brackets – Provide an adjustable, secure fit



Overall Height	2" (50 mm) 4" (100 mm) 6" (150 mm)
Unit Weight	1.4kg
S.W.L	200kg

- 50mm - No Lid (BTT149800)
- 100mm - No Lid (BTT150200)
- 150mm - No Lid (BTT150300)
- 50mm - With Lid (BTT149900)
- 100mm - With Lid (BTT150100)
- 150mm - With Lid (BTT150400)



Urinals and Bed Pans



1. Spill Proof Urinal



2. Male Urinal With Lid



3. Male Urinal



4. Male Urinal With Screw Lid



5. Female Urinal Cygnet



6. Female Urinal



7. Bed Pan



8. Bed Pan Slipper



9. Urinal Bottle Holder



10. Commode Bowl

1. Spill Proof Male Urinal	BTT148200
2. Male Urinal with Lid	BTT148450
3. Male Urinal - Green	BTT145700
3. Male Urinal - Clear	BTT160900
4. Male Urinal with Screw Lid	BTT148400
5. Female Urinal Cygnet	BTT148800
6. Female Urinal	BTT148700
7. Bed Pan	BTT151800
8. Bed Pan Slipper - Small	BTT148000
8. Bed Pan Slipper - Large	BTT161500
9. Urinal Bottle Holder	BTT149000
10. Commode Bowl with Lid and Handle	BTC081610



Assist 2 - WARD Vinyl

- Hospital grade fire retardant vinyl is easily sanitised for infection control
- Robust folding steel frame can also be fixed to stop chair folding - optional mod available
- Swing-away removable legrests
- Height adjustable footrests
- Swing-back, removable armrests
- Quick release, solid 24" tyres
- Push-to-lock brakes
- 8" height adjustable castors
- Frame fix mod to stop chair being folded as optional extra (AAL848590)



Seat Width	350 - 610mm
Seat Depth	400mm
Seat Height	460 - 500mm
S.W.L	140kg
HD S.W.L	180kg

- MWS449340 - 350mm (14")
- MWS449350 - 400mm (16")
- MWS449360 - 450mm (18")
- MWS449370 - 500mm (20")
- MWS449380 (HD) - 570mm (22")
- MWS449390 (HD) - 610mm (24")



Assist 2 - Everyday

- Folding steel frame
- Swing-away, removable height adjustable footrests
- Swing-back, removable armrests
- Quick release, solid 24" tyres
- Push-to-lock brakes
- 8" height adjustable castors
- Dacron fire retardant fabric



Seat Width	350 - 500mm
Seat Depth	400mm
Seat Height	460 - 500mm
S.W.L	140kg

- MWS449300 - 350mm (14")
- MWS449310 - 400mm (16")
- MWS449320 - 450mm (18")
- MWS449330 - 500mm (20")



Transit 2 - Transporter

- Folding steel frame
- Swing-away, removable legrests
- Height adjustable footrests
- Swing-away, removable armrests
- Quick release, solid 12" tyres
- Attendant push-to-lock brakes
- 8" height adjustable castors
- Dacron fire retardant fabric



Seat Width	400 - 500mm
Seat Depth	400mm
Seat Height	460 - 500mm
S.W.L	140kg

- MWS449410 - 400mm (16")
- MWS449420 - 450mm (18")
- MWS449430 - 500mm (20")



S.W.L
140kg

HD S.W.L
180kg

EVOKE 2 pictured with optional Contoured Padded Overlay

Evoke 2 - Lightweight & Highly Adjustable

- Highly adjustable lightweight folding wheelchair
- Swing-away & removable legrests
- Height & angle adjustable footrests
- Swing-back & removable armrests
- Quick release, solid 24" tyres, height & depth adjustable
- 8" height & angle adjustable castors
- Easy access push-to-lock brakes
- Depth & height adjustable seat
- Height, tension, angle & depth adjustable backrest
- Dacron fire retardant fabric
- Standard cushion included



Seat Width	350 - 610mm
Seat Depth	400 - 475mm
Seat Height	440 - 480mm
S.W.L	140kg
HD S.W.L	180kg


MWS449500	- 350mm (14")
MWS449510	- 400mm (16")
MWS449520	- 450mm (18")
MWS449530	- 500mm (20")
MWS449540 (HD)	- 570mm (22")
MWS449550 (HD)	- 610mm (24")



S.W.L
140kg

Evoke 2 JNR

- Highly adjustable lightweight paediatric folding wheelchair
- Swing-away & removable legrests
- Height & angle adjustable footrests
- Height adjustable & removable armrests
- Quick release, solid 20" tyres
- 8" height & angle adjustable castors
- Easy access push-to-lock brakes
- Height & depth adjustable seat
- Height, angle & depth adjustable backrest




Code	MWS449560
Seat Width	300mm
Seat Depth	300 - 375mm
Seat Height	360 - 410mm
S.W.L	140kg



S.W.L
180kg

Rehab RX JNR

- Highly configurable tilt-in-space chair for Paediatric users in the rehab environment
- Advanced articulating elevating legrests mirror the anatomical range of movement at the knee
- Quick release, solid 20" tyres, height & depth adjustable
- Attendant operated drum brakes & push-to-lock handbrakes
- Adjustable backrest, armrests, footrests, 6" castors and headrest
- PU breathable multi-stretch fabric covers
- Multilayer castellated back & seat cushions with gentle lumbar build up and lateral supports



Code	MWS449630
Seat Width	320 - 370mm
Seat Depth	440 - 490mm
Seat Height	400 - 475mm
S.W.L	180kg





Rehab RX - Advanced Tilt-In-Space




S.W.L
180kg



- Highly configurable tilt-in-space wheelchair for users with complex needs in the rehab environment
- Advanced articulating elevating legrests mirror the anatomical range of movement at the knee
- Quick release, solid 24" tyres, height & depth adjustable
- Attendant operated drum brakes & push-to-lock handbrakes
- Adjustable backrest, armrests, footrests, 6" castors and headrest
- Breathable PU multi-stretch fabric covers for pressure redistribution
- Multilayer castellated back & seat cushions with gentle lumbar build up and lateral supports

Seat Width	410 - 560mm
Seat Depth	440 - 560mm
Seat Height	400 - 475mm
S.W.L	180kg



MWS449710
410 - 460mm (16 - 18")

MWS449720
460 - 510mm (18 - 20")

MWS449740
510 - 560mm (20 - 22")



Rehab RS - Classic Tilt-In-Space



S.W.L
180kg



- Tilt and recline configurable wheelchair for the rehab environment
- Swing-away, removable, padded & elevating legrests
- Quick release, solid 24" tyres, height & depth adjustable wheels
- Attendant operated drum brakes & push-to-lock handbrakes
- Adjustable backrest, armrests, footrests, 6" castors and headrest
- PU breathable multi-stretch fabric covers
- Standard foam comfort cushions on back & seat

Seat Width	410 - 560mm
Seat Depth	440 - 560mm
Seat Height	400 - 475mm
S.W.L	180kg



MWS449600 - 410 to 460mm (16 - 18")

MWS449610 - 460 to 510mm (18 - 20")

MWS449620 - 510 to 560mm (20 - 22")





Heavy Duty Wheelchair

- Reinforced wheelchair for bariatric users
- Extra wide seat widths and reinforced construction
- Moulded composite mag wheels
- Swing-away detachable footrests
- Removable padded armrests
- Durable chrome-plated steel construction
- Vinyl wipe clean seating surfaces

Description	Heavy Duty	Super Heavy Duty
Seat Width	500 - 600mm	650 - 700mm
Seat Depth	400 - 450mm	500mm
Seat Height	500mm	500mm
S.W.L	136kg	318kg
HD S.W.L	204kg	318kg



Heavy Duty

- MWS442910 - S.W.L 136kg - 500mm (20")
- MWS442911 (HD) - S.W.L 204kg - 550mm (22")
- MWS442912 (HD) - S.W.L 204kg - 600mm (24")

Super Heavy Duty

- MWS442901 - S.W.L 318kg - 650mm (26")
- MWS442902 - S.W.L 318kg - 700mm (28")



Ward Transfer Chair




- Lightweight rigid design indoor ward transfer chair with comfort seating
- Clean and accessible design is easily sanitised for infection control
- Stackable design to maximise utilisation of space and maintain safe workflows - *see Park Station*
- Automatic and attendant braking ensures patient safety
- Frame colour differentiation allows rapid identification of hospital assets
- Swing back arms and generous seat offer comfort and support to a wide range of users
- Contoured back and seat promote correct sitting posture
- Molded PU padded surfaces are durable and easy to clean
- Optional security token access improves asset management
- Sprung flip up footplate for easy entry and egress
- IV Pole with height adjustment














Code	MWS449900
Seat Width	470mm (18")
Seat Depth	430mm
Seat Height	455mm
S.W.L	140kg



Park Station - MWS449910
























Aspire Wheelchairs Accessories Matrix

CODE	DESCRIPTION	SUITS							
		 Evoke 2	 Evoke 2 JNR	 Rehab RX JNR	 Rehab RX	 Rehab RS	 Assist 2	 Transit 2	
ELEVATING LEGRESTS									
WCA848100	Elevating Legrest - Left - 300/320/350/400/410mm		●	●	●	●	○	●	●
WCA848110	Elevating Legrest - Right - 300/320/350/400/410mm		●	●	●	●	○	●	●
WCA848120	Elevating Legrest - Left - 450/460/500/510mm		●			●	○	●	●
WCA848130	Elevating Legrest - Right - 450/460/500/510mm		●			●	○	●	●
WCA848140	Elevating Legrest - Left - 570/610mm		●					●	
WCA848150	Elevating Legrest - Right - 570/610mm		●					●	
ADVANCED ARTICULATING ELEVATING LEGRESTS									
WCA848160	Elevating Legrest - Articulating - Left - 300/320/350/400/410mm		●	○	○	●	●	●	
WCA848170	Elevating Legrest - Articulating - Right - 300/320/350/400/410mm		●	○	○	●	●	●	
WCA848180	Elevating Legrest - Articulating - Left - 450/460/500/510mm				○	●	●	●	
WCA848190	Elevating Legrest - Articulating - Right - 450/460/500/510mm				○	●	●	●	
WCA848200	Elevating Legrest - Articulating - Left - 570/610mm		●				●		
WCA848210	Elevating Legrest - Articulating - Right - 570/610mm		●				●		
STUMP SUPPORTS									
WCA848220	Stump Support - Left - Silver - 300/320/350/400/450mm		●	●	●	●	●	●	
WCA848230	Stump Support - Right - Silver - 300/320/350/400/450mm		●	●	●	●	●	●	
WCA848240	Stump Support - Left - Silver - 500/550/570/610mm		●			●	●	●	
WCA848250	Stump Support - Right - Silver - 500/550/570/610mm		●			●	●	●	
WCA848260	Stump Support - Extension Bracket - Silver		●	●			●	●	
OXYGEN BOTTLE HOLDERS									
WCA848290	Oxygen Bottle Holder		●	●			●	●	
WAA695600	Oxygen Bottle Bag				●	●	●		
BRAKE EXTENSIONS									
WCA848310	Brake Extension Lever - Left		●	●			●	●	
WCA848320	Brake Extension Lever - Right		●	●			●	●	
LATERAL SUPPORTS									
WCA848430	Lateral Chest Support Assembly - Left		●	●					
WCA848440	Lateral Chest Support Assembly - Right		●	●					

○ - Standard Feature Included

Aspire Wheelchairs Accessories Matrix

CODE	DESCRIPTION	SUITS						
		 Evoke 2	 Evoke 2 JNR	 Rehab RX JNR	 Rehab RX	 Rehab RS	 Assist 2	 Transit 2
CONTOURED PADDED OVERLAYS								
WCA848330	Contoured Padded Overlay - 350mm Wide		●					
WCA848340	Contoured Padded Overlay - 400mm Wide		●					
WCA848350	Contoured Padded Overlay - 450mm Wide		●					
WCA848360	Contoured Padded Overlay - 500mm Wide		●					
WCA848370	Contoured Padded Overlay - 570mm Wide		●					
WCA848380	Contoured Padded Overlay - 610mm Wide		●					
DEPTH & HEIGHT ADJUSTABLE ARMRESTS								
WCA848390	Depth & Height Adjustable Armrest - Left		●					
WCA848400	Depth & Height Adjustable Armrest - Right		●					
GUTTER ARMRESTS								
WCA848410	Gutter Armrest Pad - Left		●					
WCA848420	Gutter Armrest Pad - Right		●					
DRUM BRAKES								
WCA848510	Rear Wheels w/Drum Brakes Kit - 16in		●	●	●	●		
WCA848520	Rear Wheels w/Drum Brakes Kit - 20in		●	○	●	●		
WCA848500	Rear Wheels w/Drum Brakes Kit - 24in		●		○	○		
DEPTH EXTENSIONS								
WCA848530	Base Frame Extension Kit - 50mm - Black				●	●		
WCA848540	Base Frame Extension Kit - 100mm - Black				●	●		
OTHER								
WCA848550	Therapist-Friendly Tool-Free Adjustment Kit				●	●	●	
WCA848560	Meal Tray - 320mm to 410mm				●	●	●	
WCA848570	Meal Tray - 460mm to 510mm					●	●	
WCA848270	Anti-Tip Bars		●	●	○	○	○	●
WCA848280	IV Pole		●	●	●	●	●	●
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle		●	●	●	●	●	●
WCA848470	Extended Back Cane - Pair - 19/21in		●	●				
WCA848480	Headrest Assembly - Deluxe Style		●	●				
WCA848580	Calf Strap		●	●	●	●	●	●
WCA848490	Camber Kit		●	●				
AAL848590	Security Mod - ASSIST 2 WARD Only						●	


○ - Standard Feature Included



Mobile Air Chair



S.W.L
180kg



World's leading high stretch, durable and breathable healthcare fabric



Designed specifically for pressure care needs



Removable anti-bacterial/anti-microbial cover for infection control management



Waterproof yet breathable for improved skin microclimate

- Unique design promotes skin immersion and envelopment, redistributing skin pressure, reducing the risk of pressure injury associated with prolonged sitting or lying
- Lie-flat position converts chair into a bed for prolonged periods of use
- Swing away wings and armrests for improved patient transfer access
- Adjustable foam matrix in backrest promotes skin micro-climate and postural support
- Gas assisted tilt, recline and leg raise
- Removable and washable covers for infection control
- Lockable castors with directional control
- Meal tray and padded footrest optional accessories

Code	CHP198850 (Small)	CHP198855 (Large)
Overall Height	1290mm	1290mm
Overall Width	610mm	700mm
Overall Depth	690mm	690mm
Seat Height	650mm	650mm
Seat Width	410mm	500mm
Seat Depth	560mm	560mm
Recline Angle	0° to 50°	0° to 50°
Tilt Angle	0° to 25°	0° to 25°
Leg Raise Angle	0° to 90°	0° to 90°
S.W.L	180kg	180kg



Shell Chair



S.W.L
135kg



World's leading high stretch, durable and breathable healthcare fabric



Designed specifically for pressure care needs



Removable anti-bacterial/anti-microbial cover for infection control management



Waterproof yet breathable for improved skin microclimate

- Comfortable seating solution designed to maximise user access and social interaction within the healthcare environment
- Anterior and posterior tilting
- Ergonomic grip handles
- Adjustable seat height
- Safety lap belt
- Ergonomically reclining legrest
- Removable seat cushion allowing pressure care to be stepped up with after market seat cushions
- Removable washable covers for improved infection control

Code	CHP198910 (Medium)	CHP198915 (Large)
Overall Width	700mm	720mm
Overall Depth	890mm	890mm
Seat Height	450-530mm	450-530mm
Seat Width	450mm	510mm
Seat Depth	440mm	470mm
Forward Tilt Angle	0° to -3°	0° to -3°
Backward Tilt	0° to 27°	0° to 27°
Leg Raise Angle	0° to 90°	0° to 90°
S.W.L	135kg	135kg



Swish Ward Transport Chair

Comfortable multi-position patient transport chair offers both leg elevation and backrest recline

- Padded PU swing back armrests for comfort and accessibility
- Patient accessible recline function for independence
- Adjustable headrest and integrated lateral supports promote midline sitting posture
- Ergonomically designed carer push handle and operation levers promote sound workplace health and safety practices
- Highly manoeuvrable design with 300mm solid front wheels and 125mm rear castors to negotiate the busy ward environment with minimal carer effort
- Rear mounted full width foot brake lever
- Slide out durable compact laminate footplate promotes patient safety and mobility
- Safety spring in footplate allows for standing transfers without tipping the chair
- Colour differentiation allows rapid identification of hospital asset



Code	MWS450100
Overall Height	1130mm
Overall Width	643mm
Overall Depth	1000mm
Seat Height	410mm
Seat Width	500mm
Seat Depth	440mm
Recline Angle	0° to 90°
Leg Raise Angle	0° to 90°
S.W.L	180kg



Swish HD Ward Transport Chair

Comfortable multi-position heavy duty patient transport chair offers both leg elevation and backrest recline

- Generous 600mm wide seat with open arms facilitates mobilization of larger users without moving to a larger wheel base and reducing accessibility
- 300kg Safe Working Load to support bariatric patients
- Padded PU swing back armrests for comfort and accessibility
- Patient accessible recline function for independence
- Slide out durable compact laminate footplate promotes patient safety and mobility
- Safety spring in footplate allows for standing transfers without tipping the chair
- Colour differentiation allows rapid identification of hospital assets
- Adjustable headrest and integrated lateral supports promote midline sitting posture
- Ergonomically designed carer push handle and operation levers promote sound workplace health and safety practices
- Highly manoeuvrable design with 300mm solid front wheels and 125mm rear twin castors to negotiate the busy ward environment with minimal carer effort
- Rear mounted full width foot brake lever



Standard Inclusion - IV pole

Code	MWS450200
Overall Height	1130mm
Overall Width	743mm
Overall Depth	1000mm
Seat Height	410mm
Seat Width	600mm
Seat Depth	440mm
Recline Angle	0° to 90°
Leg Raise Angle	0° to 90°
S.W.L	300kg



S.W.L
160kg

Phoenix Chair

- Effective for seating clients with the following postural conditions: Pressure ulcers, Reduced upper body tone, Kyphosis, Scoliosis, Muscle contractures and Oedema
- Adjustable head support
- Arm height adjustment
- Backrest recline
- Seat depth adjustment
- Removable arms
- Elevating legrest
- Manual & Powered options



Large - Manual
(CHP197805)

Large - Electric Leg & Tilt
(CHP197810)

Large - Electric Leg, Recline & Tilt
(CHP197815)

Seat Depth	390 - 480mm
Seat Width	400 - 500mm
Total Height	1320mm
Backrest Height	850mm
Seat to Floor	600mm
Unit Weight	67kg
S.W.L	160kg



S.W.L
160kg

Phoenix Chair - Paediatric

- Effective for seating clients with the following postural conditions: Risk of pressure ulcers, Reduced upper body control, Kyphosis, Scoliosis, Muscle contractures, Swelling
- Tilt in space
- Back recline and elevating legrest
- Helps manage pressure
- Forward tilt for stand transfers
- Unique adjustable back and head support
- Fully adjustable in under 1 minute, without tools
- Shown with electric adjustment option



Code	CHP197800
Seat Depth	290 - 380mm
Seat Width	250 - 400mm
Total Height	1220mm
Backrest Height	680 - 830mm
Seat to Floor	540mm
Unit Weight	57kg
S.W.L	160kg

*Seating Matters Chairs are available in other colour options. Please call 1300 133 120 to discuss your requirements.




Sorrento Chair

- Effective for seating clients with the following postural conditions: Cerebral palsy, spinal cord injury, head injury, stroke, multiple sclerosis, Parkinson's disease, motor neurone disease, dementia, Alzheimer's disease
- Tilt in space
- Adjustable back angle recline
- Seat depth adjustment
- Elevating legrest
- Optional lateral supports and supportive pillows can help to maintain or improve posture



Code	CHP197822	CHP197823
Description	Standard	Bariatric
Seat Depth	380 - 480mm	410 - 660mm
Seat Width	350 - 600mm	550 - 850mm
Total Height	1320mm	1320mm
Backrest Height	750mm	750mm
Seat to Floor	600mm	600mm
Unit Weight	62kg	120kg
S.W.L.	160kg	294kg



Bariatric Sorrento Chair



Atlanta Chair

- Effective for seating clients with the following postural conditions: Huntington's disease, Alzheimer's, lower limb oedema, cellulitis, advanced dementia, Parkinson's disease, brain injury
- Variable angle legrest
- Arm height adjustment
- Sliding footplate
- Removable seat cushion
- Tilt and recline
- Push handle

Code	CHP197820
Seat Depth	480mm
Seat Width	380 - 500mm
Total Height	1135mm
Backrest Height	750mm
Seat to Floor	540mm
Unit Weight	67kg
S.W.L.	160kg



Monaco Chair

- Effective for seating clients with the following postural conditions: Minor stroke reduced, mobility, poor sitting balance, lower limb swelling and oedema, post operative rehab needs, pressure management needs, postural problems
- Ideal for users with reduced mobility or sitting balance
- Unique angle adjustable seat
- Adjustable back angle recline
- Padded Waterfall back
- Optional lateral supports
- Pressure management materials

Code	CHP197827
Seat Depth	350 - 550mm
Seat Width	380 - 600mm
Total Height	1205mm
Backrest Height	800mm
Seat to Floor	350 - 500mm
Unit Weight	57kg
S.W.L.	160kg



*Seating Matters Chairs are available in other colour options. Please call 1300 133 120 to discuss your requirements.

Aspire Mobile Air Chair Accessories

Meal Tray



Code	Description
CHP198861	Tray - Mobile Air Chair

Footrest



Code	Description
CHP198860	Footrest - Mobile Air Chair

Aspire Shell Chair Accessories

Pillow Headrest



Code	Description
CHP198930	Pillow Headrest - Shell Chair

Footplate



Code	Description
CHP198925	Footplate - Shell Chair

Meal Tray



Code	Description
CHP198920	Tray - Shell Chair

Seating Matters Accessories

Seat Cushion



Code	Description
CHP197910	350mm - Phoenix Paed Chair
CHP197915	450mm - Phoenix Chair
CHP197920	500mm - Atlanta Chair
CHP197925	Custom - Phoenix, Paed Phoenix & Atlanta Chairs

Head Cushion




Code	Description
CHP197930	Concave - Phoenix, Paed Phoenix & Atlanta Chairs
CHP197935	Horse Shoe - Phoenix, Paed Phoenix & Atlanta Chairs
CHP197940	Standard - Phoenix, Paed Phoenix & Atlanta Chairs

Armrest Cover

Code	Description
CHP197830	Cover - Phoenix Chair
CHP197835	Cover w/1in Padding - Phoenix Chair
CHP197840	Cover w/2in Padding - Phoenix Chair
CHP197845	Cover w/1in Padding - Phoenix Paed Chair
CHP197850	Cover w/2in Padding - Phoenix Paed Chair
CHP197855	Cover w/1in Padding - Atlanta Chair
CHP197860	Cover w/2in Padding - Atlanta Chair

Seating Matters Accessories Continued

Twin Castors



Code	Description
CHP197890	Twin Castors - Phoenix, Paed Phoenix & Atlanta Chairs

Variable Leg Rest




Code	Description
CHP197825	Variable Leg Rest - Atlanta Chair

Tray Table



Code	Description
CHP197875	Tray Table - Phoenix, Paed Phoenix & Atlanta Chairs

Removable Seat Cover



Code	Description
CHP197895	Removable Seat Cover - Phoenix, Paed Phoenix & Atlanta Chairs

Padded Footplate



Code	Description
CHP197900	Padded Footplate - Phoenix, Paed Phoenix & Atlanta Chairs

Tilt Angle Display



Code	Description
CHP197905	Tilt Angle Display - Phoenix Chair

Leg Laterals



Code	Description
CHP197864	Leg Laterals - Paed Phoenix Chair
CHP197865	Leg Laterals - Phoenix Chair

Lateral Supports



Code	Description
CHP197869	Lateral Supports - Paed Phoenix Chair
CHP197870	Lateral Supports - Phoenix Chair

Seat Belt



Code	Description
CHP197880	4 Point Seat Belt - Phoenix, Paed Phoenix & Atlanta Chairs
CHP197885	Single Seat Belt - Phoenix, Paed Phoenix & Atlanta Chairs



Compliant with ISO15883, AS4187:2014

Panamatic Optima 3

The Optima 3 is a powerful Washer Sanitiser with both hot and cold water performance. The large wash chamber accommodates a range of bed pans, urine bottles and utensils.

- Hot & cold water wash performance
- Steam disinfection cycle
- Hands-free opening
- Large capacity chamber
- Easy loading



Greater Versatility



Code	UPM100250
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output /input	Motor - 0.75 KW, Heater - 2.5 KW
Dimensions (w x d x h)	600 x 600 x 900mm
Water supply	15 mm Hot & Cold water 1 - 6 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins - 30ml Ecowash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	3 bedpans + 3 urine bottles/cycle
Warranty	1 Year




Compliant with ISO15883, AS4187:2014

Panamatic Optima 2

The Optima 2 is an efficient Washer Sanitiser for both hot and cold water performance. Hands free open, close and start helps you break the chain of infection in the Hospital environment.

- Hot & cold water wash performance
- Steam disinfection cycle
- Hands-free open, close and start
- Efficient water and power use
- MicrobeSafe antimicrobial lid
- Medium capacity chamber



Easier loading



Hands free cycle start

Code	UPM100200
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output /input	Motor - 0.75 KW, Heater - 2.5KW
Dimensions (w x d x h)	500 x 600 x 1042 mm
Water supply	15 mm Hot & Cold 1 - 6 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins / 30ml EcoWash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	2 bedpans + 2 urine bottles/cycle
Warranty	1 Year



Compliant with ISO15883, AS4187:2014

Panamatic Maxi

The Maxi Washer Sanitiser is the cold water only version of the Optima 3. The large wash chamber accommodates a range of utensils and urine bottles per cycle.

- Cold water wash performance
- Steam disinfection cycle
- Hands-free opening
- Large capacity chamber
- Easy loading



Code	UPM100350
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output /input	Motor - 0.75 KW, Heater - 2.5 KW
Dimensions (w x d x h)	600 x 600 x 900 mm
Water supply	15 mm Cold water 0.5 - 10 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins - 30ml Ecowash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	3 bedpans + 3 urine bottles/cycle
Warranty	1 Year



Compliant with ISO15883, AS4187:2014

Panamatic Midi

The Midi Washer Sanitiser is the cold water only version of the Optima 2. Hands free open, close and start helps you break the chain of infection in the Hospital environment.

- Cold water wash performance
- Steam disinfection cycle
- Hands-free open, close and start
- Efficient water and power use
- MicrobeSafe antimicrobial lid
- Easy loading



Hands free cycle start

Code	UPM100300
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output /input	Motor - 0.425KW, Heater - 2.5KW
Dimensions (w x d x h)	500 x 600 x 1042 mm
Water supply	15 mm Cold water 0.5 - 10 bar
Soil connection	110 mm connections supplied
Operation cycle time	Approx 5 mins / 30ml EcoWash
Steam disinfection	Exceeds Ao 700, 85° x 75 secs
Typical capacity	2 bedpans + 2 urine bottles/cycle
Warranty	1 Year



Compliant with ISO15883, AS4187:2014

Panamatic Mini

The Mini Washer Sanitiser is a hands free cold water only system designed for relatively low workload situations.

- Cold water wash performance
- Steam disinfection cycle
- Efficient water and power use
- Hands-free open, close & start
- MicrobeSafe antimicrobial lid
- Easier loading



Hands free cycle start

Code	UPM100290
Electrical supply	Single Phase 230 V, 13 AMP
Electrical output /input	Pump - 0.425kW
Dimensions (w x d x h)	500 x 600 x 1040 mm
Water supply	15mm Cold water 0.5 – 10 bar
Soil connection	75 mm connections supplied
Operation cycle time	Approx 6 minutes / 30ml Ecowash
Steam disinfection	Exceeds Ao 700, 85 degrees x 75 secs
Typical capacity	1 bedpan & 2 urine bottles / cycle
Warranty	1 Year



Pulpmatic Ultima

The Pulpmatic Ultima is a Medical Pulp Macerator designed to optimize cross-infection control in CSSD.

- Blockage free processing
- Hands-free open, close and start
- Large chamber capacity
- Efficient water and power consumption
- Automatic sanitation cycle



Hands free cycle start

Greater Capacity

Code	UPM100120
Electrical supply	Single Phase 230 V, 15 AMP
Electrical output /input	0.75kW, Noise: 59.5 dba
Dimensions (w x d x h)	500 x 600 x 1042 mm
Water supply	15 mm Cold water 0.5 - 10 bar
Soil connection	50 mm PTrap to 110 mm waste
Operation cycle time	115 seconds / 8ml EcoCleanse
Typical capacity	3 large or 6 small pulp products
Warranty	1 Year



Antibacterial Wipes
UPW300210

Wipes, Detergent & Deodoriser

- **Anti bacterial Wipes** - Ecocleanse+ Antibacterial Wipes clean and decontaminate frequently touched items and surfaces. The skin safe wipes are impregnated with a highly effective sanitiser and provide a long lasting impact for improved infection prevention
- **Macerator Deodoriser** - KILLS 99% OF BACTERIA AND VIRUSES. Maximum sanitisation and germ kill. Not classified as hazardous. For use in DDC Pulpmatic Macerators
- **Sanitiser Detergent** - KILLS 99% OF BACTERIA AND VIRUSES. Maximum sanitisation and germ kill. For use in DDC Washer Sanitisers



Sanitiser Detergent
UPD300110

Macerator
Deodoriser
UPD300120



4 Slot Drying Rack
UPA200270


Drying Racks

- **4 Slot Drying Rack** - stainless steel drying rack fits 4 bedpans / 4 urine bottles
- **6 Slot Drying Rack** - stainless steel drying rack fits 6 bedpans / 6 urine bottles
- **Wall Mounted Basin Rack** - stainless steel



6 Slot Drying Rack
UPA200280


Wall Mounted
Basin Rack
UPA200290



Optima 3 / Maxi
UPA200250

Utensil Trays

- Available in two different sizes to accommodate varying machine sizes



Optima 2 / Midi
UPA200240



1. Bed Pan Slipper



2. Bed Pan



3. Wash Bowls



4. Kidney Dishes



5. Urine Bottle - Female



6. Urine Bottles - Male



7. Bowl with Handle & Lid



8. Bed Pan Liner Support



9. Slipper Pan Liner Support

Pans & Bowls

1. Bed Pan Slipper - Small	BTT148000
1. Bed Pan Slipper - Large	BTT161500
2. Bed Pan	BTT151800
3. Wash Bowl - 600ml	WAT769800
3. Wash Bowl - 1500ml	WAT770000
3. Wash Bowl - 3500ml	WAT761200
3. Wash Bowl - 5500ml	WAT769600
4. Kidney Dish - 800ml	WAT770700
4. Kidney Dish - 700ml - clear	WAT770400
4. Kidney Dish - 600ml	WAT770600
4. Kidney Dish - 200ml	WAT770500
5. Urine Bottle - Female	BTT148700
6. Urine Bottle - Male - green	BTT145700
6. Urine Bottle - Male - clear	BTT160900
7. Bowl with Handle & Lid	BTC081610
8. Bed Pan Liner Support	UPP100110
9. Slipper Pan Liner Support - midi	UPP100115
9. Slipper Pan Liner Support - maxi	UPP100120



1. Pulp - Slipper Pan - Midi



2. Pulp - Slipper Pan - Maxi



3. Pulp - Bed pan



4. Pulp - Kidney Dish



5. Pulp - Male Urine Bottle



6. Pulp - GP Bowl - 1000ml



7. Pulp - GP Bowl - 3000ml



8. Pulp - Commode Bowl


Medical Pulp Shapes (Disposable) for Macerators



1. Pulp Slipper Pan - Midi - 1300ml	UPP100215
2. Pulp Slipper Pan - Maxi - 2000ml	UPP100220
3. Pulp Bedpan - 2000ml	UPP100210
4. Pulp Kidney Dish - 700ml	UPP100230
5. Pulp Male Urine Bottle - 900ml - square	UPP100245
5. Pulp Male Urine Bottle - 875ml - bulbous	UPP100250
6. Pulp GP Bowl - 1000ml	UPP100235
7. Pulp GP Bowl - 3000ml	UPP100240
8. Pulp Commode Bowl - 1700ml	UPP100225
9. Pulp Commode Bowl - Lid	UPP100226



Gynaecology Examination Couch




- Transforms from specialised gynaecology, colposcopy, urology position to normal male/female exam couch without modification
- Integrated gynae cut-out promotes accessibility and effective positioning
- Electrically adjustable backrest for patient comfort and clinically correct positioning
- Comfortable padded leg supports are also easily removable for general examination and procedures
- Healthcare grade upholstery and fully covered base frame to promote efficient infection control
- Intuitive Linak hand control and electronics for reliability and simplicity of servicing
- Integrated Trendelenburg capabilities support clinical outcomes and facilitate manual handling
- Stainless steel bowl with easy access positioning
- Central locking castors with directional functionality



BACK RAISE



HEIGHT
ADJUSTABLE



TRENDLENBURG



Code	BEB050300
Overall Length	2035mm
Overall Width	600mm
Overall Height	550 - 850mm
Patient Surface Length	1440mm
Patient Surface Width	600mm
Backrest Angle	0° to 60°
Trendelenburg	-12° to 12°
Castor Size	125mm
S.W.L	175kg



Gynaecology Delivery Bed

Designed to deliver safe, fast, effective care in the busiest of ward environments



- Retractable leg section, easily and quickly slides under the bed and out of the way for delivery
- Single action gas assisted side rails minimise transfer gaps for patient transfer
- Integrated gynae cut-out promotes accessibility and effective positioning
- Ergonomically designed patient handles
- Removable padded leg supports maximise patient postural support and can be safely stored when not in use
- Electrically adjustable for patient comfort and clinically correct positioning
- Comfortable 2-section mattress includes multi-stretch cover for patient comfort and improved infection control
- Intuitive Linak hand control and electronics for reliability and simplicity of servicing
- Integrated Trendelenburg capabilities support clinical outcomes and facilitate manual handling
- Stainless steel bowl with easy access positioning
- Tente 125mm Central locking castors with directional functionality


Standard Inclusions: Mattress, IV Pole, wall buffers, urine bag holder, belt hooks and accessory slots



BACK RAISE



HEIGHT ADJUSTABLE



TRENDLENBURG



Code	BEB050310
Overall Length	2175mm
Overall Width	1060mm
Overall Height	620 - 890mm
Patient Surface Length	2160mm
Patient Surface Width	955mm
Backrest Angle	0° to 63°
Trendelenburg	-11° to 11°
Castor Size	125mm
S.W.L	175kg



BIOCOTE
Antibacterial frame coating for improved infection control



Oxygen Bottle Holder
(PED490120)

S.W.L
150kg



Smile Paediatric Ward Bed

- Three function electrically adjustable hospital bed designed with growing people in mind
- Colourful side rails and decals to meet interior design themes and add to a more cheerful patient environment
- Capacity to lengthen at the press of a button to transition from child to young adult size proportions
- Biocote antibacterial surface treatment for improved infection control
- Durable blow molded side rails with gas spring dampeners for smooth operation
- Onboard battery provides ward mobility and emergency backup
- Integrated control panel in foot end section for carer and parent access

Code	PED490100
Overall Length	1830 - 2180mm
Overall Width	890mm
Overall Height	550 - 950mm
Patient Surface Length	1580 - 1980mm
Patient Surface Width	890mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 36°
Trendelenburg	-15° to 15°
Castor Size	150mm
S.W.L	150kg



BIOCOTE
Antibacterial frame coating for improved infection control



Clothes Basket
(PED490210)



Matris Rooming-In Crib


- Transparent side panels allow full visibility
- Ventilation sections in either end to promote flow of clean air to the baby
- Drop down side sections allow reinforced padded side sections to lay across the gap between crib and bed. This creates a:
 - More natural and intimate maternal bonding experience
 - Safe lateral transfer of the baby
- Automatically locking wheels ensure stability and safety with active handle compression required for brake release
- Integrated Trendelenburg tilt function to assist with clinical requirements such as relieving reflux
- Gas assisted height adjustment for ergonomically sound manual handling for mother and baby
- Biocote antibacterial surface treatment for improved infection control

Code	PED490200
Overall Length	890mm
Overall Width	520mm
Overall Height	490 - 710mm
Patient Surface Length	660mm
Patient Surface Width	340mm
Trendelenburg	-6° to 12°
Castor Size	125mm




Paediatric Cot & Bath - Small

- Compact and lightweight design moves easily and suits even the smallest room layouts
- Gas assisted height adjustment makes patient handling and lifting bath water a simple and ergonomic task
- Lockable mobility castors on base for ease of transportation
- Hygienic removable bath supports infection control management processes
- Convenient storage (accessory) basket on base ensures essential items are always at hand
- Powder coated frame for easy cleaning and infection control management
- Manufactured to AS2130 standard




Clothes Basket (PED490240)

Code	PED490220
Overall Length	700mm
Overall Width	500mm
Overall Height	610 - 860mm
Trendelenburg	-20° to 15°
Castor Size	50mm



Paediatric Cot & Bath - Large

- Compact and lightweight design moves easily and suits even the smallest room layouts
- Gas assisted height adjustment for safe patient handling
- Convenient storage basket (accessory) on base ensures essential items are always at hand
- Lockable mobility castors on base for ease of transportation
- Powder coated frame for easy cleaning and infection control management
- Manufactured to AS2130 standard



Clothes Basket (PED490240)

Code	PED490230
Overall Length	870mm
Overall Width	610mm
Overall Height	610 - 860mm
Trendelenburg	-20° to 15°
Castor Size	50mm





S.W.L
150kg



Beige Vinyl
(CHP198820)

Happy Recliner Chair

The Aspire Reclining Chair offers a comfortable seating surface in-room or in a mobile configuration

- Intuitive recline mechanism for patient independence
- Three comfort positions
- Adjustable padded neck pillow for postural support
- Mobility castor option for patient transportation
- Comfortable padded arms with tapered front upright for improved table access
- Lightweight and easy to move promoting workplace health and safety

Overall Height	1080mm
Overall Width	690mm
Overall Depth	780mm
Seat Height	420mm
Seat Width	520mm
Seat Depth	580mm
Backrest Height	660mm
S.W.L	150kg



S.W.L
150kg



Beige Vinyl
(CHP198810)



Light Blue Vinyl
(CHP198814)



Teal Vinyl
(CHP198812)



Navy Blue Vinyl
(CHP198816)

Happy Sleeper Chair

The Aspire Sleeper Chair facilitates the important role that family and friends play in the recovery process

- Durable seat construction preserves shape and comfort in a high use environment
- Intuitive single hand activation allows effortless conversion from chair to a comfortable sleep surface
- Recessed locking castors offer easy positioning and accessibility
- Completely flat sleep surface to maximise rest and relaxation
- Contemporary design compliments a range of care environments



Overall Height	780mm
Overall Width	760mm
Overall Depth	960mm
Seat Height	370mm
Seat Width	600mm
Seat Depth	580mm
Backrest Height	410mm
S.W.L	150kg



S.W.L
50kg

Change Table

- Compact footprint for efficient space utilisation
- Durable medical grade vinyl upholstery for easy cleaning
- Medical grade foams for durability and patient comfort
- Dual cupboards offer large internal storage area
- Three dump drawers provide additional segmented spaces for consumables management
- Laminate finish is easy to clean and durable
- Anterior barrier rail holds mattress securely for a static treatment surface



Code	WAT822650
Overall Length	1100mm
Overall Width	550mm
Overall Height	970mm
Patient Surface Length	1064mm
Patient Surface Width	514mm
S.W.L	50kg



Cot
S.W.L
75kg

Bed
S.W.L
175kg

Paediatric Cot with Trundle Bed

- Integrated adult bed allows parents to room-in with children to create a safer and more homely care environment
- Comfortable mattresses improve sleep quality for children and carers
- Cot and trundle both offer durable powder coated steel frame with mobility castors
- Integrated drawer and cupboard for essentials
- Quick release vertical side rails slide up and down for simple patient handling
- Manufactured to AS2130 standard

Code	PED490250	
Description	Paediatric Cot	Trundle Bed (included)
Overall Length	2060mm	1960mm
Overall Width	760mm	760mm
Overall Height	1160mm	300mm
Patient Surface Length	1650mm	1960mm
Patient Surface Width	700mm	760mm
Castor Size	125mm	50mm
S.W.L	75kg	175kg



S.W.L
200kg

Golvo 9000 Patient Lifter

- Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor
- Parallel base-width adjustment contributes to easier lifting of patients to/from chairs, beds and the floor
- Armrests help to support patients during ambulation, as well as making it easier for caregivers to manoeuvre the lift
- User display highlights device usage, including lift count, service required indicator and overloads
- Optional quick-release hook system means caregivers can easily and rapidly change a sling bar or install a scale
- Ergonomic hand control for lift operation displays both battery usage and overload indicators
- An extensive assortment of accessories is available including slings, slingbars and scales



Code	LSS393100
Overall Height	2090mm
Overall Length	1185mm
Overall Width	735-1020mm
Lifting Height	571-1816mm
Material	Aluminium
Leg Function	Electric
S.W.L	200kg



S.W.L
200kg

Golvo 9000 LowBase Patient Lifter

- LowBase version allows for even greater clearance under low beds or stretchers versus the standard model
- Armrests help to support patients during ambulation, as well as making it easier for caregivers to manoeuvre the lift
- Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor
- Parallel base-width adjustment contributes to easier lifting of patients or residents to/from chairs, beds and the floor as well as making it easier for caregivers to manoeuvre the lift
- User display highlights device usage, including lift count, service required indicator and overloads
- Optional quick-release hook system means caregivers can easily and rapidly change a sling bar or install a scale
- Ergonomic hand control for lift operation displays both battery usage and overload indicators
- An extensive assortment of accessories is available including slings, slingbars and scales



Code	LSS393105
Overall Height	2090mm
Overall Length	1185mm
Overall Width	745-1028mm
Lifting Height	571-1816mm
Material	Aluminium
Leg Function	Electric
S.W.L	200kg




Viking M Patient Lifter



S.W.L
205kg

- Ideal for meeting the lifting needs in small spaces
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales

Code	LSS393190
Overall Height	1340-2020mm
Overall Length	1230mm
Overall Width	690-1110mm
Lifting Height	420-1790mm
Material	Aluminium
Leg Function	Electric
S.W.L	205kg



Viking L Patient Lifter



S.W.L
250kg

- An everyday flexible lifting solution
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales

Code	LSS393180
Overall Height	2030mm
Overall Length	1360mm
Overall Width	695-1070mm
Lifting Height	460-1790mm
Material	Aluminium
Leg Function	Electric
S.W.L	250kg



Viking XL Patient Lifter



S.W.L
300kg

BARIATRIC

- Ideal for bariatric lifts and transfers
- Universal TwinBar 670 Slingbar offers additional sling space for patients of size
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales

Code	LSS393170
Overall Height	2110mm
Overall Length	1420mm
Overall Width	815-1130mm
Lifting Height	540-1910mm
Material	Aluminium
Leg Function	Electric
S.W.L	300kg





Sabina II EE Standing Lifter



S.W.L
200kg

- Conventional patient lifting to a limited extent, lifting in the sitting position
- Powerful lifting motor ensures the smooth raising and lowering of patients
- Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility
- The lower leg support and the foot rest can easily be removed for cleaning
- Adjustable height settings accommodate a variety of different patient sizes

Code	LSS393250
Overall Height	900-1000mm
Overall Length	1060mm
Overall Width	695-1115mm
Lifting Height	1050-1725mm
Material	Steel
Leg Function	Electric
S.W.L	200kg



RollOn Standing and Positioning Aid



S.W.L
160kg

- Sit to stand mobile aid supporting patients in mobilising themselves
- Padded rear supports provide security and are easily positioned and adjusted
- Safety lock ensures rear supports remain in an ideal position
- Width-adjustable base for optimal adaptability
- Crafted from aluminium for a lightweight design
- Lockable dual castors ensure stable positioning and usage

Code	LSS393360
Overall Height	1410mm
Overall Length	890mm
Overall Width	620-900mm
Material	Aluminium
Leg Function	Manual
S.W.L	160kg





Original Highback Sling with Head Support

- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester or Mesh



Size	XS	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	200kg	300kg



Polyester (Loops)


XS
(LSS394610)

SML
(LSS394620)

MED
(LSS394630)

LRG
(LSS394640)

XL
(LSS394650)



Mesh (Loops)

XS
(LSS394660)

SML
(LSS394670)

MED
(LSS394680)

LRG
(LSS394690)

XL
(LSS394700)



Soft Original Highback Sling with Head Support

- Soft pliable plate provides comfortable support when lifting from the recumbent position
- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester

Size	MED	LRG	XL	XXL
S.W.L	200kg	200kg	500kg	500kg



Polyester (Loops)

MED - Slim
(LSS394560)

MED
(LSS394570)

LRG
(LSS394580)

XL
(LSS394590)

XXL
(LSS394600)




ComfortSling with Head Support

- This sling provides comfortable sitting posture with good back support
- Adapts to the body's contours without the need for individual adjustment
- Allows the patient to remain in the sling after application and transfer
- Enhanced width increases patient comfort
- Available in Mesh



Size	MED	LRG	XL
S.W.L	200kg	200kg	300kg




Mesh (Loops)

MED
(LSS394710)

LRG
(LSS394720)

XL
(LSS394730)




Liko

S.W.L 200kg

SilhouetteSling with Head Support

Size	XS	SML	MED	LRG
S.W.L	200kg	200kg	200kg	200kg

- Pliable sling that adapts to the body and requires minimal wheelchair space
- Recommended for lifting to sitting to/from close-fitting wheelchairs
- Provides a slightly reclined sitting posture suitable for patients with poor trunk stability
- Adjustable head support
- Suitable to be used during bathing
- Available in Mesh



Mesh (Loops)

XS (LSS394450)

SML (LSS394460)

MED (LSS394470)

LRG (LSS394440)



Liko

S.W.L 200kg S.W.L 300kg

HygieneSling

Size	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	300kg

- Easy to apply sling that enables quick, smooth and secure transfers
- Equipped with a large opening that facilitates dressing and undressing
- Works well with patients regardless of what clothing they are wearing
- Recommended for lifting to sitting to/from close-fitting wheelchairs
- Suitable to be used during bathing
- Available in Polyester




Polyester (Loops)

SML (LSS394480)

MED (LSS394490)

LRG (LSS394500)

XL (LSS394510)



Liko

S.W.L 200kg

AmputeeSling

Size	SML	MED	LRG
S.W.L	200kg	200kg	200kg

- Designed for safe and secure lifting of patients with amputated legs
- Supports patient with single and double leg amputations as well as patients with contractures
- Can also be used for non-amputee patients who are slender and require a smaller sling
- Offers a sitting posture that is slightly reclined with leg supports applied parallel to each under both thighs
- Available in Polyester



Polyester (Loops)

SML (LSS394514)

MED (LSS394516)

LRG (LSS394518)



AmputeeSling with Head Support

- Designed for safe and secure lifting of patients with amputated legs
- Supports patient with single and double leg amputations as well as patients with contractures
- Can also be used for non-amputee patients who are slender and require a smaller sling
- Offers a sitting posture that is slightly reclined with leg supports applied parallel to each under both thighs
- Adjustable head support
- Available in Polyester



Size	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	300kg



Polyester (Loops)
 SML (LSS394520)
 MED (LSS394530)
 LRG (LSS394540)
 XL (LSS394550)



UltraSling Bariatric Sling

- Designed for safe and comfortable lifting of bariatric patients
- Supports the entire back up to the shoulders
- Leg supports are installed individually on the sling bar, which gives the stomach space and the legs optimal support
- Available in Short, Regular and Long lengths
- Available in Mesh & Polyester



Size	SHORT	REGULAR	LONG
S.W.L	500kg	500kg	500kg



Polyester (Loops)
 SHORT (LSS394400)
 REGULAR (LSS394410)
 LONG (LSS394420)



Mesh (Loops)
 REGULAR (LSS394430)



Solo Highback Disposable Sling with Head Support

- Designed for single patient use
- Can be used in most lifting situations
- Designed to adapt to the patient without any individual adjustments
- Made of strong non-woven polypropylene fabric that is non-washable and engineered to be disposed of when it has become soiled or when the patient no longer requires it
- Supplied in boxes of 20 (Med-Lrg) & 8 (XL)



Size	MED	LRG	XL
S.W.L	200kg	200kg	500kg



Polyester (Loops)
 MED - Slim (LSS394740)
 MED (LSS394750)
 LRG (LSS394760)
 XL (LSS394770)



SupportVest Standing Sling

Size	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	200kg



- Comfortable sling which lifts behind the back and under the arms
- Easy to use and requires only minimal individual adjustment so suit any patient
- Soft padding redistributes pressure evenly
- Supportive belt holds the vest in place at all times
- Available in Polyester



Polyester
(Loops)
SML
(LSS394120)
MED
(LSS394130)
LRG
(LSS394140)
XL
(LSS394150)

S.W.L
200kg



ComfortVest Standing Sling

Size	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	200kg



- Sit-to-stand sling designed especially for patients who must avoid underarm pressure
- Outer vest supports the patient behind the back and outside the arms, while the inner waist belt prevents gliding and keeps the patient from sliding out of the vest
- Self-regulating tensioning function
- Ideal choice for patients who are paralysed on one side or have other conditions that inhibit underarm pressure
- Available in Polyester



Polyester
(Loops)
SML
(LSS394180)
MED
(LSS394190)
LRG
(LSS394200)
XL
(LSS394210)

S.W.L
200kg




SafetyVest Standing Sling

Size	MED	LRG
S.W.L	200kg	200kg



- Sit-to-stand sling designed to assist patients with poor stability
- Self-regulating tensioning function ensures both caregiver and patient feel secure
- Padded fabric provides a high degree of friction and effectively prevents gliding
- Optional Sheepskin padding is available for even more comfort
- Available in Polyester



Polyester
(Loops)
MED
(LSS394160)
LRG
(LSS394170)

S.W.L
200kg



S.W.L
200kg

MasterVest Gait Training Sling

- Versatile training aid that facilitates toilet visits, standing training, gait training, mobility training and placing patients into standing positions
- Reduces strain on caregivers while helping them deliver physical training
- Transfers to an from the toilet can be accomplished easily and smoothly, reducing strain
- Lifts around the upper body and partly under the arms, and is best suited for patients that can support some of their body weight
- Optional leg harness for improved support
- Available in Polyester



Size	SML	MED	LRG	XL	XXL
S.W.L	200kg	200kg	200kg	200kg	200kg
	Polyester (Loops)			Leg Harness (Optional Accessory)	
	SML (LSS394260)			SML (LSS394310)	
	MED (LSS394270)			MED/LRG (LSS394320)	
	LRG (LSS394280)			XL/XXL (LSS394330)	
	XL (LSS394290)				
	XXL (LSS394300)				



S.W.L
200kg

HygieneVest Toileting Sling

- Designed for safe and comfortable transfers to and from the toilet
- Offers almost fully upright sitting, providing functional and safety advantages
- Works well in complicated situations, such as when lifting amputees or patients with poor muscle tone
- Vest section firmly supports the upper body during transfer, while a generous opening leaves the lower body free for dressing and undressing
- Available in Polyester



Size	SML	MED	LRG	XL
S.W.L	200kg	200kg	200kg	200kg
	Polyester (Loops)			
	SML (LSS394220)			
	MED (LSS394230)			
	LRG (LSS394240)			
	XL (LSS394250)			




S.W.L
200kg

S.W.L
500kg

LiftPants Gait Training Sling

- Enables caregivers the ability to provide gait training to their patients while providing a sense of security for patients
- Offer freedom of movement while relieving the burden of body weight
- Provide necessary support to allow patient to move on his or her own with confidence
- Can be used on overhead or mobile patient lifters
- Available in Mesh



Size	SML	MED	LRG	XL	XXL	XXXL
S.W.L	200kg	200kg	200kg	500kg	500kg	500kg
	Mesh (Loops)					
	SML (LSS394340)			XL (LSS394370)		
	MED (LSS394350)			XXL (LSS394380)		
	LRG (LSS394360)			XXXL (LSS394390)		



S.W.L
50kg

TeddySling Paediatric Sling with Head Support

- Enable a relaxed, reclined seated position
- Can be used for all transfers to and from the horizontal and seated positions
- Perfect for transfers to and from the floor
- Equipped with a reinforced head support allowing a more upright seated position
- Ideal for children with spasticity problems
- Available in Teddy Patterned Polyester

Code	LSS394100	LSS394110
Size	XS	SML
S.W.L	50kg	50kg



Polyester
(Loops)



Liko Accessories



Liko Scale



Liko SlingBar



LikoStretch



Liko Clips

Please contact Aidacare to discuss delivery, installation and accessory options



A205 Aluminium Lifter

- Compact yet versatile mobile patient lifter ideal for high-use facility environments
- Ergonomic handset - Linak technology fits neatly in the hand with an intuitive control panel
- Powered leg function - Opens and closes the legs for lifting stability or accessibility
- Quick Release Yoke with Safety Clips - One touch release and attach mechanism. Padded yoke for added patient safety
- Safety stop button - Emergency stop button, emergency lowering and anti-entrapment technology
- High quality castors - 2x 75mm front castors, 2x 100mm rear locking castors

Code	LSS390610
Overall Height	2045mm
Overall Length	1135mm
Overall Width	620mm - 1095mm
Lifter Weight	38.1kg
S.W.L	205kg



A320 Aluminium Bariatric Lifter

Heavy duty yet manoeuvrable bariatric patient lifter ideal for the most demanding clinical environments.

- Ergonomic handles & support arms - provide multiple grip points. Folding support arms assist with stability or mobility
- Safety stop button - emergency stop button, emergency lowering and anti-entrapment technology
- Ergonomic handset - Linak technology fits neatly in the hand with an intuitive control panel
- Powered leg function - opens and closes the legs for lifting stability or accessibility
- Quick release 4 point yoke - one touch release and attach mechanism
- High quality castors - 2x 100mm front castors, 2x 125mm rear locking castors

Code	LSS390630
Overall Height	2165mm
Overall Length	1505mm
Overall Width	765mm - 1190mm
Lifter Weight	63kg
S.W.L	320kg



Pivot Frame

- Designed to be fitted on Aidacare Aspire Patient Lifters
- Pivot frame allows carers to pivot and manoeuvre a patient during transfers
- Powered pivot function option ensures that carers reduce the amount of manual handling
- Allows a more comfortable user seating position and more control for the carer
- Designed to be used in conjunction with Clip/Pivot Slings



Manual Pivot Frame (LSS390670)



Compatible Lifters	A205 Lifter, A320 Lifter
Compatible Slings	Aspire CLIP Slings
S.W.L	227kg



In-Line Weigh Scale

- Designed to be fitted on Aidacare Aspire Patient Lifters, however adaptable to most Lifter brands
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display



Code	LSS390680
Compatible Lifters	A150F Lifter, A205 Lifter, 200 Lifter, A320 Lifter
S.W.L	320kg



Deluxe General Purpose Sling


Size	S.M.L to M.E.D	LRG to XL
S.W.L	200kg	300kg



- Versatile sling offering body contoured fabric shaped for maximum patient comfort and support
- Available in Polyester or Mesh
- For FIM score 1 or 2
- Sizing colour reference on edging of slings for easy visual identification



Material	Size	S.W.L
Polyester (Loops)	S.M.L	(LSS390700)
	M.E.D	(LSS390705)
	LRG	(LSS390710)
	XL	(LSS390715)



Material	Size	S.W.L
Mesh (Loops)	S.M.L	(LSS390716)
	M.E.D	(LSS390717)
	LRG	(LSS390718)
	XL	(LSS390719)



Deluxe General Purpose Sling with Head Support

Size	S.M.L to M.E.D	LRG to XXXL
S.W.L	200kg	300kg



- This sling has a comfortable padded head support with rigidisers for posture management
- Body contoured fabric
- Available in Polyester or Mesh for longer use periods and Loop or Clip options.
- Sizing colour reference on edging of slings for easy visual identification



Material	Size	S.W.L
Polyester (Clips)	S.M.L	(LSS390900)
	M.E.D	(LSS390905)
	LRG	(LSS390910)
	XL	(LSS390915)



Material	Size	S.W.L
Polyester (Loops)	S.M.L	(LSS390740)
	M.E.D	(LSS390745)
	LRG	(LSS390750)
	XL	(LSS390755)
	XXL	(LSS3907551)
	XXXL	(LSS3907552)



Material	Size	S.W.L
Mesh (Clips)	S.M.L	(LSS390916)
	M.E.D	(LSS390917)
	LRG	(LSS390918)
	XL	(LSS390919)



Material	Size	S.W.L
Mesh (Loops)	S.M.L	(LSS390756)
	M.E.D	(LSS390757)
	LRG	(LSS390758)
	XL	(LSS390759)
	XXL	(LSS3907591)
	XXXL	(LSS3907592)




Deluxe Hygiene Access Sling

Size	S.M.L to M.E.D	LRG to XXXL
S.W.L	200kg	300kg



- This sling offers generous underarm padding, quilted leg straps and a reinforced wide-cut back.
- Padded polyester fabric for increased patient comfort
- Security strap features velcro and buckle fastening for easy positioning and added safety
- Sizing colour reference on edging of slings for easy visual identification



Material	Size	S.W.L
Polyester (Loops)	S.M.L	(LSS390760)
	M.E.D	(LSS390765)
	LRG	(LSS390770)
	XL	(LSS390775)
	XXL	(LSS390776)
	XXXL	(LSS390777)






Deluxe Hygiene Access Sling with Head Support


Size	S.M.L to M.E.D	L.R.G to X.L
S.W.L	200kg	300kg



- Adjustable head support for additional patient safety and comfort.
- Padded poly fabric and generous under arm padding.
- Security belt with velcro and buckle
- For FIM score 2 or 3
- Available in padded polyester or polyester with clips
- Sizing colour reference on edging



Polyester (Loops)
SML (LSS390780)
MED (LSS390785)
LRG (LSS390790)
XL (LSS390795)



Polyester (Clips)
SML (LSS390940)
MED (LSS390945)
LRG (LSS390950)
XL (LSS390955)



Deluxe Hammock Sling with Head Support


Size	S.M.L to M.E.D	L.R.G to X.L
S.W.L	200kg	300kg



- This full body sling provides complete comfort and support for high needs patients.
- Available in Mesh for general use or Supersoft for long periods of use
- Parasilk leg strap lining for easy positioning
- For FIM score 1 or 2
- Sizing colour reference on edging



Mesh (Loops)
SML (LSS390860)
MED (LSS390865)
LRG (LSS390870)
XL (LSS390875)



Supersoft (Loops)
SML (LSS390800)
MED (LSS390805)
LRG (LSS390810)
XL (LSS390815)




Deluxe Low Friction Parasilk Sling


Size	S.M.L to M.E.D	L.R.G to X.L
S.W.L	200kg	300kg



- Parasilk material is low friction for easy positioning
- Supports the whole body and is suitable for general transfers, including: seated and bed to chair transfers
- Clip option available
- For FIM score 1 or 2
- Sizing colour reference on edging



Parasilk (Loops)
SML (LSS390880)
MED (LSS390885)
LRG (LSS390890)
XL (LSS390895)





Parasilk (Clips)
SML (LSS390960)
MED (LSS390965)
LRG (LSS390970)
XL (LSS390975)



Deluxe Standing Sling


- Generous underarm padding and a reinforced wide-cut back for comfort
- Quilted and padded under arm and on back to enhance patient comfort, stability and rigidity
- Security strap with velcro and buckle fastening for easy positioning and added safety.
- Available in padded polyester or polyester with clips
- For FIM score 3 or 4
- Sizing colour reference on edging



Size	SML to MED	LRG to XXXL
S.W.L	200kg	300kg



Polyester (*Loops*)
 SML
 (LSS390720)
 MED
 (LSS390725)
 LRG
 (LSS390730)
 XL
 (LSS390735)
 XXL
 (LSS390736)
 XXXL
 (LSS390737)



Polyester (*Clips*)
 SML
 (LSS390920)
 MED
 (LSS390925)
 LRG
 (LSS390930)
 XL
 (LSS390935)




Evacuation Sheet

Code	TSE672526
Length	1900mm
Width	900mm
S.W.L	300kg



- Fits discreetly underneath the resident's mattress secured by cross and corner straps
- Allows bed making to occur as normal
- Pulling straps located at head and foot ends
- Remains readily accessible in case of emergency evacuation
- Is fire retardant and highly durable
- Low friction surface
- Suitable for use by one or two carers




Slide Sheet

Size	Small	Large
Code	TSE673070	TSE673080
Length	1000mm	2000mm
Width	1450mm	1450mm



- Nylon material designed to assist in sliding and transferring patients
- Low friction fabric slides on itself to facilitate patient bed mobility
- High quality Nylon fabric and stitching for improved durability
- Aspire slide sheets can be laundered for ongoing patient hygiene




Non-Slip One Way Slide

Code	TSE673050
Length	440mm
Width	450mm



- Assists in the correct positioning of a person in the back of a chair
- Helps maintain posture for the user, preventing them from slipping forward
- Reduces the risk of strain or injury from a fall or slip
- Works in conjunction with most chairs whether a gap exists between the back and seat or not
- Machine washable





Foot Operated One Way Slide

Code	TSE673052
Length	1300mm
Width	450mm



- Non-slip mesh grips securely to smooth surfaces such as vinyl upholstery
- Soft seating surface can remain under patient throughout the day and facilitate multiple repositions
- Foot strap threads through back of the chair to assist the carer with effortless positioning into seat
- Machine washable



Positioning Slide - Foot Operated

Code	TSE673054
Length	1375mm
Width	460mm



- Designed with Care staff in mind, low friction surfaces allow for easy patient manoeuvring with minimal lift
- For patients with a lower risk of slipping while sitting in a chair
- Manoeuvre the patient easily to the back of the chair or sideways for transfer assistance
- Foot strap threads through back of chair to assist the carer with positioning
- Machine Washable



Swivel Cushion

Code	TSE673056
Width	440mm



- Low friction pivot allows easy sit and swivel entry and exit from vehicles
- Minimises exertion and the likelihood of injury, and improves the angle of access to or from mobility equipment
- Plush velour fabric and generous padding for maximum comfort on long car trips
- Machine washable



Ceiling Hoist Considerations


Aidacare is an experienced supplier and installer of ceiling hoist systems.

We understand the risks involved in the installation and use of ceiling hoists and can help during early-stage specification and design to reduce risk and improve user-acceptance

Our systems are world-class and are accompanied by accredited installation guidelines which allow us to work with a range of installers. We are an ISO 9001 accredited company and consistently provide high-quality outcomes for Australia's largest Aged Care and Hospital providers

A wide range of motor & track options are available to suit client specifications

Some common examples are provided below




- A. Structural Attachment Point
- B. Threaded Rod
- C. Ceiling Bracket
- D. Ceiling Track
- E. Traverse Boom (Track)
- F. Charging Station
- G. X-Y Gantry
- H. Ceiling Lift Unit
- I. Carry Bar
- J. Lift Strap
- K. Sling
- L. Gait Assembly
- M. Linear Track
- N. Modified Door Header


Consulting Solutions

Aidacare has extensive experience with installation and user-acceptance of ceiling systems. Simple things can be identified during planning to reduce risk in the installation and improve user and patient outcomes. We can assist in developing the best solution. Some examples of installations are provided below.


Linear Track with Turn Table




Straight Monorail



Turntable




J Track



Consulting Solutions

No two installations are alike. Aidacare has extensive experience dealing with many different installation substrates and scenarios. Aidacare utilises manufacturer-certified installation guidelines. Some of the more complex examples are shown below. This leads to a consistent high-quality outcome for our clients.



Track Options



Maintenance Safety

Aidacare fully-supports its installations with ongoing maintenance programs. Aidacare also assists our clients to develop their own internal maintenance regime.

Training & User Acceptance

Aidacare Equipment Specialists offer user training and ongoing support. We also provide WHS signage with all equipment.

Quality & Compliance

All systems are tested and installed in compliance with Australian Standard AS/NZ ISO 10535-2011 Hoists for the transfer of disabled persons.

All systems are maintained in accordance with manufacturer's guidelines and in accordance with AS/NZ 3551-2012 Management Programs for Medical Equipment.

Aidacare is an ISO9001 compliant company.

Aidacare uses an accredited network of trained installers, engineers, and load testers.





P200 Portable Ceiling Hoist



- With its unique innovative design, the P200 is one of the lightest portable hoists on the market, yet it is capable of lifting up to 200kg
- Suitable for use with fixed ceiling tracks or mobile gantry systems
- The design and features make it perfect for use in community or nursing environments
- Quick Release System (QRS) allows for the P200 to be safely and easily attached or removed from the gantry or track system with an easy reacher
- Pneumatic handset is ideal for use in wet environments
- Emergency lowering allows for the safe lowering of the in emergency situations
- Client safety is assured with locking latches to secure slings in place for all lifts
- Room to Room transfer allows patients to be transferred without the need to be disconnected from the hoist

Manual Traverse
(LIC350120)

Manual Traverse
& Quick Release System
(LIC350110)

Quick Release System reacher bar
- sold as optional accessory -
(LIC354170)



Height Range	Up to 2070mm
Weight	4.8kg
S.W.L	200kg



CP136 Portable Ceiling Hoist



- Slimline, low profile, padded carry bar that has six loop attachments as standard (four-point hook optional). Different width carry bars also available
- Low profile hoist unit is unobtrusive and provides greater lifting height
- All lifting and traversing operations have soft start/stop to ensure smooth transfers for the patient
- A host of safety features
- Smooth running tape with large lifting range
- State of the art display with lift counter
- Ergonomically designed handset with hold to run clearly marked buttons and directional arrows requires minimum pressure to operate. Pneumatically controlled, making it safe for use in wet conditions
- Fire retardant cover completely encapsulates hoist
- Neat in track transformer, IP 67 rated
- Coloured direction indicators on hoist unit match hand control to ensure correct operation first time
- Emergency stop and manual emergency window



Code	LIC350100
Height Range	Up to 1960mm
Weight	3.8kg
S.W.L	136kg

PRISM MEDICAL

Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



Transactive Xtra Ceiling Hoist

- The Transactive XTRA fixed ceiling hoist allows for a large array of configurations, catering to changing needs, making it a suitable fixed lifting solution
- Digital display screen provides essential information for the user and caregiver
- Return to Charge (RTC) option automatically returns the hoist to its charger when required
- Manual wind down and raising allows the hoist to continue to be used during a power outage
- Low profile hoist unit is unobtrusive and provides greater lifting height
- Pneumatic handset is ideal for use in wet environments
- Smooth lifting/running tape with large lifting range
- Emergency Lowering allows for the safe lowering of the in emergency situations
- Client safety is assured with locking latches to secure slings in place for all lifts
- Available in manual or powered traverse options

Manual

- Manual 130kg - (LIC350130)
- Manual 160kg - (LIC350150)
- Manual 200kg - (LIC350170)
- Manual 270kg - (LIC350200)
- Manual Curtain Jumper 270kg - (LIC350220)

Powered

- Powered 130kg - (LIC350140)
- Powered 160kg - (LIC350160)
- Powered 200kg - (LIC350180)
- Powered RTC 200kg - (LIC350190)
- Powered 270kg - (LIC350210)
- Powered RTC 270kg - (LIC350240)
- Powered Curtain Jumper 270kg - (LIC350230)

S.W.L
130-270kg
 Dependant on Model

Height Range Up to 2080mm
Weight 9.7kg (manual) / 10.4kg (powered)
S.W.L 130 - 270kg

PRISM MEDICAL

Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



C380 Bariatric Fixed Ceiling Hoist

- The C380 is specifically designed to lift clients that exceed a standard hoists maximum weight capacity
- As the strongest Ceiling Track Hoist in its class, the C380 safely transfers bariatric clients with respect and dignity
- Low profile hoist unit is unobtrusive and provides greater lifting height
- Digital display screen provides essential information for the user and caregiver
- Return to Charge (RTC) option automatically returns the hoist to its charger when required
- Pneumatic handset is ideal for use in wet environments
- Smooth lifting/running tape with large lifting range
- Emergency Lowering allows for the safe lowering of the in emergency situations
- Client safety is assured with locking latches to secure slings in place for all lifts
- Available in manual or powered traverse options



S.W.L
380kg

Manual Traverse (LIC350250)

Powered Traverse (LIC350260)

Powered Traverse & RTC (LIC350270)

Height Range Up to 2080mm
Weight 10.4kg
S.W.L 380kg



Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



S.W.L 272kg S.W.L 363 kg



LikoGuard L/XL Fixed Ceiling Hoist

- LikoGuard L & XL hoists combine powerful performance and enhanced safety design features
- Newly designed 450 SlingGuard slingbar encourages the correct placement of sling loops for greater confidence in safe patient mobilisation
- Innovative emergency stop and manual emergency lowering mechanism that is easy to reach, engage and operate with no tools, ladder or force needed
- Two choices of IPX7* certified hand controls, one LED and one LCD backlit display. Both easy-to-use hand controls features battery, service and overload-indication. Advanced and intuitive LCD hand control offers graphical user interface with tracking of user and service data
- The LikoGuard lift system offers lift intervals of 230cm and supports quick release hooks, fixed sling bars and stretchers
- Choice of on-wall or in-rail charger

Code	LIC359300	LIC359310
Size	Large	X-Large
Height Range	Up to 2300mm	Up to 2300mm
Weight	14kg	14kg
S.W.L	272kg	363kg



Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



S.W.L 200kg

Multirall 200 Portable Ceiling Hoist

- The Multirall 200 hoist is a portable ceiling hoist that can also be used for room-to-room patient transfers
- Can be used to transfer a patient between rooms in a sitting position, easily moving from one rail system to another with the use of a transfer strap
- It's portability makes it an excellent and economical solution for nursing home and hospital care situations where lifting needs are temporary
- Can be easily attached or removed with no tools and minimum effort from ceiling rails to facilitate convenient transport between rooms
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Can be mounted two ways for ease of use

Code	LIC359100
Height Range	Up to 1600mm
Weight	8.7kg
S.W.L	200kg



Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



Likorall 200 Fixed Ceiling Hoist

- The Likorall 200 is an economical fixed ceiling hoist needed to manage everyday lifting situations
- The patented safety drum with speed limiter facilitates electrical emergency lowering
- Ideal for floor lifts in most environments
- Can be equipped with a special carriage for planned transfers between different rooms
- Ergonomically designed hand control
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Simple wall charging via hand control



Code	LIC359490
Height Range	Up to 1800mm
Weight	8.7kg
S.W.L	200kg



Exclusively Distributed by:
 aidacare
 HEALTHCARE EQUIPMENT



Likorall 242 Fixed Ceiling Hoist

- The Likorall 242 ceiling hoist is designed to meet demanding lifting requirements
- The system is an ideal choice for those who want to move the lift motor along a track electronically with a transfer motor, or to control the system with a wireless hand control option
- The patented safety drum with speed limiter facilitates electrical emergency lowering
- Can be equipped with a special carriage for planned transfers between different rooms
- Wireless hand control
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Simple wall charging via hand control or in rail charging



Likorall 242 Manual (LIC359535) Likorall 242 Powered Room to Room (LIC359500) Likorall 242 Manual Room to Room (LIC359510)

Height Range	Up to 2000mm
Weight	8.7kg
S.W.L	200kg



Exclusively Distributed by:
aidacare
 HEALTHCARE EQUIPMENT



Likorall 250 Fixed Ceiling Hoist

- The Likorall 250 ceiling hoist is designed to meet demanding lifting requirements
- The system is an ideal choice for those who want to move the lift motor along a track electronically with a transfer motor, or to control the system with a wireless hand control option
- The patented safety drum with speed limiter facilitates electrical emergency lowering
- Can be equipped with a special carriage for planned transfers between different rooms
- Wireless hand control
- Quick sling bar change thanks to Q-link and Quick-release Hook
- Simple wall charging via hand control or in rail charging



Likorall 250 Powered
(LIC359520)

Likorall 250 Manual
(LIC359530)

Height Range	Up to 1800mm
Weight	8.7kg
S.W.L	250kg


Prism Ceiling Hoist Accessories & Installation



Clax Transport Cart
(WAT758020)



Free Standing Gantry
(LIC352100)



Fitout and Installation Solutions

Please contact Aidacare to discuss delivery, installation and accessory options

Liko Ceiling Hoist Accessories & Installation



Liko SlingBars



Liko Gate Training Sling



LikoStretch

Please contact Aidacare to discuss delivery, installation and accessory options



Single

Linen Skip - No Lid

- Stainless steel frame
- Open top (no lid)
- Fully welded joints assists with infection control
- Heavy duty 100mm swivel castors for safe manual handling
- Fixed uprights for swift bag removal



Double



Triple

Code	LAC500800	LAC500810	LAC500820
Size	Single	Double	Triple
Length	430mm	430mm	430mm
Width	430mm	860mm	1280mm
Height	940mm	940mm	940mm



Soiled Linen Skip Bag Support - Stainless Steel (LAC500920)



Buffer Disc and Reducing Bush to Suit Linen Skips (LAC500930)

Soiled Linen Skip Push Handle - Single - Stainless Steel (LAC500925)



Single

Linen Skip - Foot Operated Lid

- Stainless steel frame
- Foot operated lid operation
- Fully welded joints assists with infection control
- Heavy duty 100mm swivel castors for safe manual handling
- Fixed uprights for swift bag removal
- Coloured lid labels included



Double



Triple

Code	LAC500860	LAC500870	LAC500880
Size	Single	Double	Triple
Length	430mm	430mm	430mm
Width	460mm	920mm	1360mm
Height	945mm	945mm	945mm



Lid for Linen Skip White - With stickers (LAC500636)

Buffer Disc and Reducing Bush to Suit Linen Skips (LAC500930)

Soiled Linen Skip Push handle - Single - Stainless Steel (LAC500925)

Soiled Linen Skip Bag Support - Stainless Steel (LAC500920)



Single

Economy Linen Skip - No Lid

- Powder coated steel frame
- Open top (no lid)
- Fully welded joints assists with infection control
- 75mm swivel castors for safe manual handling
- Fixed uprights for swift bag removal



Double



Triple

Code	LAC500830	LAC500840	LAC500850
Size	Single	Double	Triple
Length	430mm	430mm	430mm
Width	430mm	860mm	1290mm
Height	890mm	890mm	890mm



Buffer Disc and Reducing Bush to Suit Linen Skips (LAC500930)



Single

Economy Linen Skip - Foot Operated Lid

- Powder coated steel frame
- Foot operated lid operation
- Fully welded joints assists with infection control
- 75mm swivel castors for safe manual handling
- Fixed uprights for swift bag removal
- Coloured lid labels included



Double



Triple

Code	LAC500890	LAC500900	LAC500910
Size	Single	Double	Triple
Length	430mm	430mm	430mm
Width	430mm	860mm	1290mm
Height	890mm	890mm	890mm




Lid for Linen Skip
White - With stickers
(LAC500636)

Buffer Disc and Reducing
Bush to Suit Linen Skips
(LAC500930)



Soiled Linen Bags

- Replacement soiled linen bags to suit linen skips
- Available in full size, half size and waterproof variations
- Available in a range of colours



Code	Type	Colour	Capacity
LAC353130	Standard	White	15 kg
LAC353140	Standard	Blue	15 kg
LAC353160	Standard	Green	15 kg
LAC353170	Standard	Red	15 kg
LAC353180	Standard	Yellow	15 kg
LAC353190	Standard	Purple	15 kg
LAC353060	Waterproof	White	15 kg
LAC353020	Waterproof	Blue	15 kg
LAC353040	Waterproof	Green	15 kg
LAC353050	Waterproof	Red	15 kg
LAC353070	Waterproof	Yellow	15 kg
LAC353030	Waterproof	Black	15 kg
LAC353105	Half Size (OH&S)	White	10 kg
LAC353106	Half Size (OH&S)	Blue	10 kg
LAC353107	Half Size (OH&S)	Green	10 kg
LAC353108	Half Size (OH&S)	Red	10 kg
LAC353109	Half Size (OH&S)	Yellow	10 kg
LAC353110	Half Size (OH&S)	Purple	10 kg

Soiled Linen Skip Lids

- Available in a range of colours



Green Plastic
(LAC354110)



Red Plastic
(LAC354130)



White Plastic
(LAC354140)



Yellow Plastic
(LAC354150)




Blue Plastic
(LAC354100)

Moist Linen Trolley



- Mobile trolley for easy transport
- Designed to transport and dispose soiled linen in hospital environments
- Strong galvanised base
- Durable polyethylene linen tub
- 4x heavy duty 125mm swivel castors for safe manual handling



Code	LAC354025	LAC354030
Size	Small	Large
Length	1000mm	1360mm
Width	710mm	750mm
Height	790mm	790mm
Capacity	320L	450L

Push Handle
LAC500935

Moist Linen Trolley Lid
Small - LAC500940
Large - LAC500945



Spring Base Option
Small - LAC354062
Large - LAC354064

Wet & Dry Linen Trolley



Small

- Lightweight yet robust mesh design
- Spring assisted adjustable base
- Spring covers reduce linen entrapment
- Small has 4x HD 125mm rubber swivel castors.
- Medium has 2x HD 150mm rubber fixed castors + 2x HD 150mm rubber swivel castors.
- Large has 4x HD 200mm rubber swivel castors
- Powder coated steel frame
- Self lubricating anti rattle platform guides
- 4x corner rotating wall buffers



Push Handle - Bolt on
Medium - LAC354220
Large - LAC354225

Code	LAC353380	LAC354190	LAC354210
Size	Small	Medium	Large
Length	800mm	1100mm	1400mm
Width	600mm	680mm	800mm
Height	800mm	850mm	850mm



Medium



Large

Economy Clean Linen Trolley




Large - open
(3 Shelf)

- Durable lightweight design
- Raised top tray edging prevents spillages
- Powder coated steel frame
- 4x corner buffers prevents damage to walls
- Multi height push handle (large model only)
- Directional castor (large model only)

Code	LAC353640	LAC353610	LAC353620	LAC353590
Size	Small - Open	Large - Open	Small - Closed	Large - Closed
Length	920mm	1100mm	920mm	1100mm
Width	525mm	550mm	525mm	550mm
Height	1040mm	1400mm	1040mm	1400mm



Small - closed
(2 shelf)




Front Cover
Small - LAC353630
Large - LAC353600

Deluxe Clean Linen Trolley



- Durable lightweight design
- Stainless steel tow hitch connection
- Raised top tray edging prevents spillages
- Two tone powder coated steel frame
- 4x corner buffers prevents damage to walls
- Multi height push handle
- 2x directional castors

Code	LAC353560
Length	1285mm
Width	625mm
Height	1430mm



Front Cover
Large - LAC353570





Clean Linen Trolley
Tow Hitch
Large - LAC353575

Linin Folding Table



- 4x heavy duty 125mm swivel castors
- Raised PVC work surface
- Durable robust construction
- Powder coated steel frame
- Promotes quick and efficient linin storage time

Code	LAC501000	LAC501010
Operation	Fixed Height	Electric
Length	1200mm	1200mm
Width	650mm	650mm
Height	900mm	900mm




Garment Dispensing Trolley



Garment Dispensing
Trolley (15 Trays)

- Durable lightweight design
- Powder coated steel frame
- 4x corner buffers prevents damage to walls
- Push handle on both ends
- Name tags to identify patient belongings

Code	LAC353650	LAC353660	LAC353680	LAC353470
Capacity	15 Trays	18 Trays	30 Trays	Tray
Length	1100mm	1100mm	1450mm	430mm
Width	490mm	490mm	700mm	325mm
Height	1300mm	1500mm	1375mm	130mm




Garment Rack



- Epoxy Powder coated frame
- Flat pack base and upright design
- Nestable 'Z' Base design

Code	LAC353462
Type	Z Base
Length	1425mm
Width	650mm
Height	1750mm




Garment Rack - Z base

Stock Trolley



- Lightweight robust construction
- Fixed perforated basket
- Push Handle on both ends
- Heavy duty 125mm polyurethane swivel castors
- Swivel buffer discs for wall protection
- Epoxy powder coated finish

Code	WAT825510
Length	700mm
Width	450mm
Height	950mm



Food Tray Dispensing Trolley



8 Tray Trolley

- Stainless steel top shelf
- Push handle at both ends
- Stainless steel frame
- 2x directional lock castors
- 2x directional swivel castors
- Corner buffers
- Can be manufactured to suit any tray size or configuration



16 Tray Trolley



30 Tray Trolley

Code	WAT823500	WAT823600	WAT824000
Capacity	8 Tray	16 Tray	24 Tray
Length	460mm	575mm	575mm
Width	615mm	900mm	1330mm
Height	1490mm	1490mm	1490mm

Food Tray - 450 x 350 mm - Plastic Small - WAT824700

Cover for 8 Tray Trolley Pockets & 4 Grips - WAT824716

Cover for 16 Tray Trolley Pockets & 4 Grips - WAT824717

Cover for 24 Tray Trolley Pockets & 4 Grips - WAT824715

Food Tray Clearing Trolley



2 Shelf Trolley




3 Shelf Trolley

- Lightweight stainless steel, fully welded construction
- 40mm safe edge on top of tray to prevent spills
- Push handle at both ends
- Stainless steel frame
- 4x 125mm grey rubber swivel castors

Code	WAT824400	WAT824600
Size	2 Shelf	3 Shelf
Length	800mm	800mm
Width	500mm	500mm
Height	950mm	950mm

Food Tray - 450 x 350mm - Plastic Small - WAT824700




File Trolley



- 25mm diameter push handle at both ends
- Lightweight
- Swivel wall buffers
- Storage Tray
- File trolley lid available

Code	WAT900910
Size	File Trolley
Length	950mm
Width	420mm
Height	980mm

File Trolley Lid WAT900920



Flow Chart Trolley



- Lightweight aluminium work platform 670mm x 800mm
- Spring assisted high / low adjustment
- 3x flow chart clips
- Pen & implement retainer
- Stainless steel uprights
- Epoxy powder coated top
- Stable 4x point pre-zinc treated epoxy powder coated contoured base
- 100mm total lock grey rubber castors

Code	WAT900940
Width	800mm
Height	850-1300mm



Arm Operating Table

Code	WAT820190
Length	850mm
Width	450mm
Height	670-825mm

- Robust stainless steel construction with no cross bars allowing access from both sides
- X-ray translucent top
- Height adjustable 670-825mm
- 100mm total lock grey rubber castors



Instrument / Dressing Trolley (Stainless Steel)



Instrument Trolley - Square Frame
WAT900800



Instrument Trolley - Square Frame
+ 1 Drawer - WAT900810

- Stainless steel construction
- Total lock castors on front wheels
- Guard rail on 3 sides
- Drawers fitted with high quality runners with soft close magnetic latch
- 140mm deep drawers



Instrument Trolley - Square Frame
+ 2 Drawers - WAT900820



Instrument Trolley - Square Frame
+ 4 Drawers - WAT900830



Dressing Trolley - Square Frame
WAT900840



Dressing Trolley - Square Frame
+ 1 Drawer - WAT900850



Dressing Trolley - Square Frame
+ 2 Drawers - WAT900860



Dressing Trolley - Square Frame
+ 4 Drawers - WAT900870




Dressing Trolley - Square
Frame + 2 Drawers side by
side - WAT900861

Trolley	Instrument	Dressing
Length	500mm	800mm
Width	500mm	500mm
Height	1040mm	1040mm

Instrument / Dressing Trolley - Flat Top (Stainless Steel)



- Stainless steel construction
- Flat top fully welded
- 100mm total lock castors



Length	500/800/1000mm
Width	500mm
Height	950mm

Instrument Trolley - Flat Top
500 x 500 x 950mm - WAT900880

Instrument Trolley - Flat Top
500 x 800 x 950mm - WAT900890

Instrument Trolley - Flat Top
500 x 1000 x 950mm - WAT900900

eCart Range

- Lightweight composite aluminium/plastic panels, reduces cart noise while in motion and panels can be replaced economically if damaged
- Drawers fitted with soft close runners and dampeners - prevents drawers from opening while in motion
- Height adjustable push handles
- Raised drawer handles for grip and easy cleaning
- Recessed ABS plastic top and storage compartment
- Slide out shelf fitted with runners for smooth, stable operation
- Keyed mechanical lock or optional digital lock
- 3x 125mm anti-static total lock castors and 1x directional lock castor

Economy Cart - 2 Drawers Fixed Push Handle (Webster)



WAT901100	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1075mm

Economy Cart - 3 Drawers Fixed Push Handle & Dividers (12)



WAT901110	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1075mm

Medication Cart - 3 Drawers (Webster)



WAT901120	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

Medication Cart - 4 Drawers Dividers (12)



WAT901130	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

Medication Cart - 5 Drawers Dividers (12)



WAT901140	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

Med/Theatre Wide Cart 4 Drawers, Dividers (12 & 14)



WAT901150	Overall Length	570mm
	Overall Width	870mm
	Overall Height	1065mm

Med/Theatre Wide Cart 5 Drawers, Dividers (12 & 14)



WAT901160	Overall Length	570mm
	Overall Width	870mm
	Overall Height	1065mm

Anaesthetic Cart - 5 Drawers
Blue, Stainless Steel Top, Guardrail,
Tilt Containers (17)



WAT901170	Overall Length	545mm
	Overall Width	760mm
	Overall Height	1580mm

Emergency Cart - 5 Drawers
Red, Oxy carrier, IV Pole,
Cardio Board



WAT901180	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1075mm

Theatre Equip Cart - 5 Drawers
Stainless Steel Top, Tilt Containers



WAT901190	Overall Length	541mm
	Overall Width	732mm
	Overall Height	1574mm

Theatre Equip Cart - 5 Drawers
Equip Shelves, Accessory Module



WAT901200	Overall Length	537mm
	Overall Width	997mm
	Overall Height	1724mm

Sachet Cart - 6 Drawers
Dividers (10)



WAT901210	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

Procedure Cart - 5 Drawers
Grey, IV Pole, Storage Basket



WAT901220	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

Isolation Cart - 4 Drawers
Yellow, Glove Dispenser Holder,
Dividers (12)



WAT901230	Overall Length	520mm
	Overall Width	760mm
	Overall Height	1045mm

eCart Accessories



Accessory Module



Basket



Bottle Cabinet - Fitted with Key Lock



Cardiac Board






Chart Holder




Dividers (1x set) (4x3=12)



Equipment Shelf - Double



Equipment Shelf - Single



Guard Rail



Intravenous Pole with Corner Basket



Oxygen Bottle Carrier




Sharps Holder



Stainless Steel Top




Tablet Arm




Tilt Container - 4



Tilt Container - 5



Tilt Container - 6



Waste Bin



Webster Clips (Pair) to Suit Wire Webster Rack



Glove Dispenser



Corner Bracket for Universal IV Pole



Anti-Tamper Lock



Tray

Accessory Module	WAT902100
Basket	WAT902110
Bottle Cabinet - Fitted with Key Lock	WAT902120
Cardiac Board	WAT902130
Chart Holder	WAT902140
Dividers (1x Set) (4x3 = 12)	WAT902150
Divider Name Tabs	WAT902160
Defibrillator Arm	WAT902170
Digital Lock for eCart	WAT902180
Digital Lock Key - Spare - Blue	WAT902190
Digital Lock Key - Master - Black	WAT902200
Digital Lock Key - Program - Yellow	WAT902210
Equipment Shelf - Double	WAT902220
Equipment Shelf - Single	WAT902230
Folding Side Shelf	WAT902240
Guard Rail	WAT902250
Intravenous Pole with Corner Bracket	WAT902260
Lock Rail	WAT902270
Oxygen Bottle Carrier	WAT902280
Suture Cabinet	WAT902290
Sharps Holder	WAT902300
Stainless Steel Top	WAT902310
Tablet Arm	WAT902320
Tilt Container - 4	WAT902330
Tilt Container - 5	WAT902340
Tilt Container - 6	WAT902350
Tilt Container - 9	WAT902360
Mount Arms for Tilt Container (Pair)	WAT902370
Waste Bin	WAT902380
Webster Clips (Pair) to Suit Wire Webster Rack	WAT902390
Glove Dispenser	WAT902400
Corner Bracket for Universal IV Pole	WAT902261
Anti-Tamper Lock	WAT902105
Drive Unit	WAT902440
Tray	WAT902315



Multiline NEXT IT Treatment Chair

- 3 motor treatment chair (back recline, seat tilt and leg raise)
- Innovative five position padded armrests
- Ergonomic push handle with paper roll holder
- Infusion pole holder can be pivoted left and right with just one hand
- Vinyl covering with no patient facing seams, promoting stringent infection control
- Soft upholstery promotes comfort whilst a 60cm wide seat caters for a large range of patients
- Market leading side accessibility, promoting safe manual handling
- Battery backup for use in power outages or ward transport
- Full +/- Trendelenburg functionality & foot control for assistance in emergency situations
- Genuine chair position with leg section at 90° angle.
- Included IV Pole Holder can be pivoted left and right with just one hand
(IV Pole is an optional accessory)



Product safety
GREINER treatment
couches are independently
tested and certified



Seated position



Relax position



Bed/flat position



Trendelenburg & shock position

Code	CHT100100
Overall Height	1400mm
Overall Width	800mm
Overall Depth	1100mm
Flat Position Length	1920mm
Seat Height	540mm
Seat Width	600mm
Seat Depth	600mm
Backrest Angle	-5° to 78°
Footrest Angle	0° to 90°
S.W.L	220kg



Multiline NEXT DC Treatment Chair

- 5 motor treatment chair (seat height, back recline, seat tilt, leg raise, footrest)
- Ergonomic push handle with paper roll holder
- Footrest with automatic, patient friendly length adjuster
- Market leading side accessibility, promoting safe manual handling
- Premium vinyl covering with no patient facing seams, promoting stringent infection control
- Innovative armrest technology with individual positioning for height, inclination, direction and volume, with adjustment to width of torso – ideal for infusion therapy
- Soft upholstery promotes comfort whilst a 60cm wide seat caters for a large range of patients
- Battery backup for use in power outages or ward transport
- Full +/- Trendelenburg functionality & foot control for assistance in emergency situations
- Protective footrest cover extends the life of vinyl coverings around the foot area
- Included IV Pole Holder can be pivoted left and right with just one hand
(IV Pole is an optional accessory)



Product safety
GREINER treatment
couches are independently
tested and certified



Seated position



Relax position



Bed/flat position



Trendelenburg & shock position



Code	CHT100210 CHT100200
Overall Height	1400 - 1610mm
Overall Width	800mm
Overall Depth	1480mm
Flat Position Length	2200mm
Seat Height	540 - 750mm
Seat Width	600mm
Seat Depth	600mm
Backrest Angle	-5° to 78°
Footrest Angle	-25° to 0°
S.W.L	220kg





Multiline XXL Bariatric Treatment Chair




- 3 motor bariatric treatment chair (back recline, seat tilt and leg raise)
- True bariatric transport and treatment solution
- Full sit-to-stand functionality
- Easy access with fold down armrests and adjustable access height of 55 to 80cm
- Padded fold down armrests
- Adjustable push handles for ease of movement
- Soft upholstery promotes comfort whilst a large 70cm wide seat caters for bariatric users
- Battery backup for use in power outages or ward transport
- Full +/- Trendelenburg functionality
- Included IV Pole Holder can be pivoted left and right with just one hand (IV Pole is an optional accessory)

Optional Accessories



LED reading light



A range of custom colours are available upon request (Lead times may vary)



IV Pole



Sit-to-stand position



Seated position



Relax position



Trendelenburg & shock position



Bed/flat position



Code	CHT100300
Overall Height	1340 - 1600mm
Overall Width	890mm
Overall Depth	1300mm
Flat Position Length	2200mm
Seat Height	540 - 750mm
Seat Width	700mm
Seat Depth	700mm
Backrest Angle	-5° to 78°
Footrest Angle	0° to 90°
S.W.L	300kg



Tilting Examination Couch

- 2 section height adjustable patient surface also tilts from horizontal to vertical to assist with examinations and promote good posture for the therapist and practitioner
- Secure patient safety straps assist with positioning and support
- Gas assisted independent head section movement
- Durable medical grade vinyl upholstery with robust powder coated steel frame for easy cleaning and sanitising in support of infection control processes
- Medical grade foams for patient comfort and durability
- Linak electronics with intuitive hand control to provide simple and efficient patient handling
- Base frame allows access for lifting devices
- 125mm locking castors


Code	WAT822630
Overall Length	1865mm
Overall Width	600mm
Overall Height	560 - 820mm
Backrest Angle	0° to 82°
Headrest Angle	0° to 30°
Castor Size	125mm
S.W.L	175kg



Bobath Therapy Couch

- Specifically designed for rehabilitation involving patients with neurological conditions
- Generous width therapy surface for mobilisation of patients up to 220kg SWL
- Height adjustable patient surface assists with therapy and examinations
- Gas assisted independent head section movement
- Durable medical grade vinyl upholstery with robust powder coated steel frame for easy cleaning and sanitising in support of infection control processes
- Medical grade foams for patient comfort and durability
- 125mm locking castors
- Linak electronics with intuitive hand control to provide simple and efficient patient handling


Code	WAT822640
Overall Length	2400mm
Overall Width	1200mm
Overall Height	560 - 820mm
Backrest Angle	0° to 69°
Castor Size	125mm
S.W.L	220kg



Examination Couch

- 2 section adjustable height examination couch
- Head section with gas assisted movement for +/- angle adjustment
- Lightweight and compact design with 200mm head end castors for easy manoeuvring
- Paper towel roll holder for convenience and hygiene management to avoid cross contamination
- Linak electronics with intuitive hand control
- Durable medical grade vinyl upholstery for easy cleaning
- Medical grade foams for durability and patient comfort

Code	WAT822620
Overall Length	2055mm
Overall Width	650mm
Overall Height	585 - 935mm
Headrest Angle	-62° to 62°
Castor Size	105mm (foot) / 200mm (head)
S.W.L	175kg





Examination Table

- 2 section fixed height examination table
- Gas assisted head section movement with integrated comfort pillow
- Durable medical grade vinyl upholstery for easy cleaning
- Medical grade foams for durability and patient comfort
- Robust powder coated steel frame can be cleaned easily in support of infection control processes

Code	WAT822600
Overall Length	1860mm
Overall Width	600mm
Overall Height	665mm
Headrest Angle	0° to 30°
S.W.L	150kg



Examination Table with Drawers

- 2 section fixed height examination table with drawers
- Gas assisted head section with integrated comfort pillow
- Integrated drawers for management of consumable items
- Durable medical grade vinyl upholstery for easy cleaning
- Medical grade foams for durability and patient comfort
- Stylish timber laminate surround trim to compliment the interior environment
- Robust powder coated steel frame can be cleaned easily in support of infection control processes

Code	WAT822610
Overall Length	1860mm
Overall Width	600mm
Overall Height	665mm
Headrest Angle	0° to 30°
S.W.L	150kg



Examination Stool

- Practitioners swivelling Examination Stool with backrest
- Gas assisted height adjustment
- 5 legs with rolling castors for stable manoeuvrability
- Comfort padded vinyl seat
- Contoured backrest

Code	WAT803850
Overall Height	560 - 750mm





Botero Bariatric Ward Bed

The specific bariatric design focus offers width and length adjustable deck to manage changing needs in both ward and respite environments.



Integrated Bed Length Extension
180mm of extension to cater for taller patients



Attendant Control Panel
Single Touch Functions for Speed and Simplicity

- Low height of 380mm
- Antibacterial powder coating for improved infection control management
- Backrest translation (autoregression) reduces compression of the thoracic and abdominal areas as the backrest is raised and alleviates shear and pressure for improved skin integrity
- Compatible with all leading bed movers
- Attendant control panel with intuitive pre-programmed single touch functions for speed and simplicity
- 150mm Tente Castors
- Central locking castors with integrated directional lock
- Full Trendelenburg functionality
- Integrated folding side rails
- Removable width adjustable head and foot panels



S.W.L
400 kg



BIOCOTE
Antibacterial frame coating for improved infection control

Integrated Bed Width Extension
1020mm - 1220mm of tool free extension for changing needs

Optional Accessories



Static Mattress w/ Extension Bolster (BEM050200)




Self Help Pole (BEA825100)



Telescopic IV Pole (BEA825101)



Code	BEB050200
Overall Length	2260 - 2440mm
Overall Width	1020 - 1220mm
Overall Height	380 - 850mm
Patient Surface Length	2000 - 2180mm
Patient Surface Width	1020 - 1220mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 30°
Trendelenburg	-17° to 17°
Castor Size	150mm
SWL	400kg



Bariatric Static Mattress

- Four layer foam design caters for bariatric users
- Multi-stretch cover with welded seams, reinforced waterfall flaps and water resistant zippers for stringent infection control
- Side handles allow for easy mattress movement
- Strong high density base and side walls improve patient transfer efficiency



Code	BEM046700
Overall Length	2030mm
Overall Width	1050mm
Overall Height	200mm
S.W.L	300kg
Warranty	5 years foam / 2 years cover




Botero Bariatric Static Mattress

- Two layer foam design with high density base and castellated surface caters to bariatric patients
- Multi-stretch cover with welded seams for stringent infection control
- Strong high density side walls improve patient transfer efficiency while cut back hinge points allow the mattress to easily move with the mattress platform



Code	BEM050200
Overall Length	2000mm
Overall Width	1220mm
Overall Height	200mm
S.W.L	400kg
Warranty	5 years foam / 2 years cover




Weight Bearing Platform Commode - 600mm

Reinforced heavy duty platform commode for users of size

- Height adjustable, reinforced and padded armrests with safety ends drop flat and out of the way to assist with positioning or lateral slide transfers
- Open back straps provide access for bathing
- Easy-slide weight bearing platform footrest slides completely away for transfers
- 4x total locking stainless steel castors provide brake access from all angles




Propelling Code	Attendant BTC061400
Overall Height	995mm
Overall Width	725mm
Seat Width	600mm
Backrest Height	300mm
Armrest Height	290mm
Anterior Seat	590mm
Posterior Seat	560mm
S.W.L	400kg



Seats Sold Separately



Bariatric Shower Commode



- Super Heavy Duty Bariatric design with 790mm seat and 350kg SWL
- Padded seat and arms for maximum comfort. Removable arms to assist with positioning and allow lateral slide transfers
- 4x lockable castors maximise safety by providing carer access to a brake at all angles around the commode
- Height adjustable, swing away footrests for optimal positioning and ease of transfers

Code	BTC061200
Overall Height	990mm
Overall Width	900mm
Seat Width	790mm
Backrest Height	430mm
Armrest Height	170mm
Anterior Seat	570mm
Posterior Seat	570mm
S.W.L	350kg



Swish HD Ward Transport Chair

Comfortable multi-position heavy duty patient transport chair offers both leg elevation and backrest recline



- Generous 600mm wide seat with open arms facilitates mobilization of larger users without moving to a larger wheel base and reducing accessibility
- 300kg Safe Working Load to support bariatric patients
- Padded PU swing back armrests for comfort and accessibility
- Patient accessible recline function for independence
- Slide out durable compact laminate footplate promotes patient safety and mobility
- Safety spring in footplate allows for standing transfers without tipping the chair
- Colour differentiation allows rapid identification of hospital assets
- Adjustable headrest and integrated lateral supports promote midline sitting posture
- Ergonomically designed carer push handle and operation levers promote sound workplace health and safety practices
- Highly manoeuvrable design with 300mm solid front wheels and 125mm rear twin castors to negotiate the busy ward environment with minimal carer effort
- Rear mounted full width foot brake lever

Standard Inclusion - IV pole



Code	MWS450200
Overall Height	1130mm
Overall Width	743mm
Overall Depth	1000mm
Seat Height	410mm
Seat Width	600mm
Seat Depth	440mm
Recline Angle	0° to 90°
Leg Raise Angle	0° to 90°
S.W.L	300kg




MAXI Shower Chair

- Super heavy duty adjustable width shower chair
- Soft padded seat, back and arms provide postural support
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the chair
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability



Code	BTS118020
Overall Height	670 - 790mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg




MAXI Shower Stool

- Super heavy duty adjustable width shower stool
- Soft padded seat and arms provide postural support
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability



Code	BTS118030
Overall Height	670 - 790mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg



MAXI Over Toilet Aid

- Super heavy duty adjustable width over toilet aid
- Soft padded flip up seat for easy cleaning
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
- To create a MAXI bedside commode (BEC045120), simply add a padded backrest (as shown in MAXI Shower Chair above)



Code	BTT145020
Overall Height	670 - 790mm
Overall Width	645 - 845mm
Overall Depth	495mm
Seat Height	455 - 575mm
Sitting Width	575 - 775mm
Seat Depth	455mm
S.W.L	310kg





Super Heavy Duty Wheelchair

- Reinforced wheelchair for bariatric users
- Extra wide seat widths and reinforced construction
- Moulded composite mag wheels
- Swing-away detachable footrests
- Removable padded armrests
- Durable chrome-plated steel construction
- Vinyl wipe clean seating surfaces



Description	Super Heavy Duty
Seat Width	650 - 700mm
Seat Depth	500mm
Seat Height	500mm
S.W.L	318kg
HD S.W.L	318kg


Super Heavy Duty

MWS442901	650mm (26")
MWS442902	700mm (28")



Sorrento Chair

- Effective for seating clients with the following postural conditions: Cerebral palsy, spinal cord injury, head injury, stroke, multiple sclerosis, Parkinson's disease, motor neurone disease, dementia, Alzheimer's disease
- Tilt in space
- Adjustable back angle recline
- Seat depth adjustment
- Elevating legrest
- Optional lateral supports and supportive pillows can help to maintain or improve posture



Code	CHP197823
Description	Bariatric
Seat Depth	410 - 660mm
Seat Width	550 - 850mm
Total Height	1320mm
Backrest Height	750mm
Seat to Floor	600mm
Unit Weight	120kg
S.W.L	294kg

Bariatric Sorrento Chair



Bobath Therapy Couch

- Specifically designed for rehabilitation involving patients with neurological conditions
- Generous width therapy surface for mobilisation of patients up to 220kg SWL
- Height adjustable patient surface assists with therapy and examinations
- Gas assisted independent head section movement
- Durable medical grade vinyl upholstery with robust powder coated steel frame for easy cleaning and sanitising in support of infection control processes
- Medical grade foams for patient comfort and durability
- 125mm locking castors
- Linak electronics with intuitive hand control to provide simple and efficient patient handling



Code	WAT822640
Overall Length	2400mm
Overall Width	1200mm
Overall Height	560 - 820mm
Backrest Angle	0° to 69°
Castor Size	125mm
S.W.L	220kg





Viking XL Patient Lifter



- Ideal for bariatric lifts and transfers
- Universal TwinBar 670 Slingbar offers additional sling space for patients of size
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales

Code	LSS393170
Overall Height	2110mm
Overall Length	1420mm
Overall Width	815-1130mm
Lifting Height	1370-1910mm
Material	Aluminium
Leg Function	Electric
S.W.L	300kg



Soft Original Highback Sling with Head Support



- Soft pliable plate provides comfortable support when lifting from the recumbent position
- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester

Size	XL	XXL
S.W.L	500kg	500kg



Polyester (Loops)
 XL
 (LSS394590)
 XXL
 (LSS394600)



ComfortSling with Head Support



- This sling provides comfortable sitting posture with good back support
- Adapts to the body's contours without the need for individual adjustment
- Allows the patient to remain in the sling after application and transfer
- Enhanced width increases patient comfort
- Available in Mesh

Size	XL
S.W.L	300kg



Mesh (Loops)
 XL
 (LSS394730)





Original Highback Sling with Head Support


- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester or Mesh



Size XL
S.W.L 300kg



Polyester (Loops)
XL
(LSS394650)



Mesh (Loops)
XL
(LSS394700)



UltraSling Bariatric Sling

- Designed for safe and comfortable lifting of bariatric patients
- Supports the entire back up to the shoulders
- Leg supports are installed individually on the sling bar, which gives the stomach space and the legs optimal support
- Available in Short, Regular and Long lengths
- Available in Mesh & Polyester



Size SHORT REGULAR LONG
S.W.L 500kg 500kg 500kg



Polyester (Loops)
SHORT
(LSS394400)

REGULAR
(LSS394410)

LONG
(LSS394420)



Mesh (Loops)
REGULAR
(LSS394430)



Solo Highback Disposable Sling with Head Support

- Designed for single patient use
- Can be used in most lifting situations
- Designed to adapt to the patient without any individual adjustments
- Made of strong non-woven polypropylene fabric that is non-washable and engineered to be disposed of when it has become soiled or when the patient no longer requires it
- Supplied in boxes of 20 (Med-Lrg) & 8 (XL)



Size XL
S.W.L 500kg



Polyester (Loops)
XL
(LSS394770)



AmputeeSling with Head Support

- Designed for safe and secure lifting of patients with amputated legs
- Supports patient with single and double leg amputations as well as patients with contractures
- Can also be used for non-amputee patients who are slender and require a smaller sling
- Offers a sitting posture that is slightly reclined with leg supports applied parallel to each under both thighs
- Adjustable head support
- Available in Polyester



Size XL
S.W.L 300kg



Polyester (Loops)
XL (LSS394550)



HygieneSling

- Easy to apply sling that enables quick, smooth and secure transfers
- Equipped with a large opening that facilitates dressing and undressing
- Works well with patients regardless of what clothing they are wearing
- Recommended for lifting to sitting to/from close-fitting wheelchairs
- Suitable to be used during bathing
- Available in Polyester



Size XL
S.W.L 300kg



Polyester (Loops)
XL (LSS394510)



LiftPants Gait Training Sling

- Enables caregivers the ability to provide gait training to their patients while providing a sense of security for patients
- Offer freedom of movement while relieving the burden of body weight
- Provide necessary support to allow patient to move on his or her own with confidence
- Can be used on overhead or mobile patient lifters
- Available in Mesh



Size XL XXL XXXL
S.W.L 500kg 500kg 500kg



Mesh (Loops)
XL (LSS394370)
XXL (LSS394380)
XXXL (LSS394390)





A320 Aluminium Bariatric Lifter

Heavy duty yet manoeuvrable bariatric patient lifter ideal for the most demanding clinical environments.

- Ergonomic handles & support arms - provide multiple grip points. Folding support arms assist with stability or mobility
- Safety stop button - emergency stop button, emergency lowering and anti-entrapment technology
- Ergonomic handset - link technology fits neatly in the hand with an intuitive control panel
- Powered leg function - opens and closes the legs for lifting stability or accessibility
- Quick release 4 point yoke - one touch release and attach mechanism
- High quality castors - 2x 100mm front castors, 2x 125mm rear locking castors

Code	LSS390630
Overall Height	2165mm
Overall Length	1505mm
Overall Width	765mm - 1190mm
Lifter Weight	63kg
S.W.L	320kg



Deluxe General Purpose Sling

- Versatile sling offering body contoured fabric shaped for maximum patient comfort and support
- Available in Polyester or Mesh
- For FIM score 1 or 2
- Sizing colour reference on edging of slings for easy visual identification

Size	XL
S.W.L	300kg



Polyester (Loops)	Mesh (Loops)
XL (LSS390715)	XL (LSS390719)



Deluxe General Purpose Sling with Head Support

- This sling has a comfortable padded head support with rigidisers for posture management
- Body contoured fabric
- Available in Polyester or Mesh for longer use periods and Loop or Clip options.
- Sizing colour reference on edging of slings for easy visual identification

Size	XL - XXXL
S.W.L	300kg



Polyester (Loops)	Mesh (Loops)	Polyester (Clips)
XL (LSS390755)	XL (LSS390759)	XL (LSS390915)
XXL (LSS390751)	XXL (LSS390751)	
XXXL (LSS390752)	XXXL (LSS390752)	
		Mesh (Clips)
		XL (LSS390919)







Deluxe Hygiene Access Sling

Size XL - XXXL
S.W.L 300kg



- This sling offers generous underarm padding, quilted leg straps and a reinforced wide-cut back.
- Padded polyester fabric for increased patient comfort
- Security strap features velcro and buckle fastening for easy positioning and added safety
- Sizing colour reference on edging of slings for easy visual identification



Polyester (Loops)

XL (LSS390775)
XXL (LSS390776)
XXXL (LSS390777)




Deluxe Hygiene Access Sling with Head Support


Size XL
S.W.L 300kg



- Adjustable head support for additional patient safety and comfort.
- Padded poly fabric and generous under arm padding.
- Security belt with velcro and buckle
- For FIM score 2 or 3
- Available in padded polyester or polyester with clips
- Sizing colour reference on edging



Polyester (Loops)
XL (LSS390795)



Polyester (Clips)
XL (LSS390955)



Deluxe Hammock Sling with Head Support


Size XL
S.W.L 300kg



- This full body sling provides complete comfort and support for high needs patients.
- Available in Mesh for general use or Supersoft for long periods of use
- Parasilk leg strap lining for easy positioning
- For FIM score 1 or 2
- Sizing colour reference on edging



Mesh (Loops)
XL (LSS390875)



Supersoft (Loops)
XL (LSS390815)





Deluxe Low Friction Parasilk Sling

Size XL
S.W.L 300kg




S.W.L
300kg

- Parasilk material is low friction for easy positioning
- Supports the whole body and is suitable for general transfers, including: seated to seated and bed to chair transfers
- Clip option available
- For FIM score 1 or 2
- Sizing colour reference on edging


XL



Parasilk (Loops)
XL
(LSS390895)



Parasilk (Clips)
XL
(LSS390975)



Deluxe Standing Sling

Size XL - XXXL
S.W.L 300kg



S.W.L
300kg

- Generous underarm padding and a reinforced wide-cut back for comfort
- Quilted and padded under arm and on back to enhance patient comfort, stability and rigidity
- Security strap with velcro and buckle fastening for easy positioning and added safety.
- Available in padded polyester or polyester with clips
- For FIM score 3 or 4
- Sizing colour reference on edging


XL

XXL


XXXL



Polyester (Loops)
XL (LSS390735)
XXL (LSS390736)
XXXL (LSS390737)



Polyester (Clips)
XL (LSS390935)




Evacuation Sheet

Code TSE672526
Length 1900mm
Width 900mm
S.W.L 300kg



S.W.L
300kg

- Fits discreetly underneath the resident's mattress secured by cross and corner straps
- Allows bed making to occur as normal
- Pulling straps located at head and foot ends
- Remains readily accessible in case of emergency evacuation
- Is fire retardant and highly durable
- Low friction surface
- Suitable for use by one or two carers





C380 Bariatric Fixed Ceiling Hoist



- The C380 is specifically designed to lift clients that exceed a standard hoist's maximum weight capacity
- As the strongest Ceiling Track Hoist in its class, the C380 safely transfers bariatric clients with respect and dignity
- Low profile hoist unit is unobtrusive and provides greater lifting height
- Digital display screen provides essential information for the user and caregiver
- Return to Charge (RTC) option automatically returns the hoist to its charger when required
- Pneumatic handset is ideal for use in wet environments
- Smooth lifting/running tape with large lifting range
- Emergency Lowering allows for the safe lowering of the in emergency situations
- Client safety is assured with locking latches to secure slings in place for all lifts
- Available in manual or powered traverse options



Manual Traverse
(LIC350250)

Powered Traverse
(LIC350260)

Powered Traverse
& RTC (LIC350270)

Height Range	Up to 2080mm
Weight	10.4kg
S.W.L	380kg



LikoGuard XL Fixed Ceiling Hoist



- LikoGuard L & XL hoists combine powerful performance and enhanced safety design features
- Newly designed 450 SlingGuard slingbar encourages the correct placement of sling loops for greater confidence in safe patient mobilisation
- Innovative emergency stop and manual emergency lowering mechanism that is easy to reach, engage and operate with no tools, ladder or force needed
- Two choices of IPX7* certified hand controls, one LED and one LCD backlit display. Both easy-to-use hand controls features battery, service and overload-indication. Advanced and intuitive LCD hand control offers graphical user interface with tracking of user and service data
- The LikoGuard lift system offers lift intervals of 230cm and supports quick release hooks, fixed sling bars and stretchers
- Choice of on-wall or in-rail charger



Code	LIC359310
Size	XL
Height Range	Up to 2300mm
Weight	14kg
S.W.L	363kg



Multiline XXL Bariatric Treatment Chair




- 3 motor bariatric treatment chair (back recline, seat tilt and leg raise)
- True bariatric transport and treatment solution
- Full sit-to-stand functionality
- Easy access with fold down armrests and adjustable access height of 55 to 80cm
- Padded fold down armrests
- Adjustable push handles for ease of movement
- Soft upholstery promotes comfort whilst a large 70cm wide seat caters for bariatric users
- Battery backup for use in power outages or ward transport
- Full +/- Trendelenburg functionality
- Included IV Pole Holder can be pivoted left and right with just one hand (IV Pole is an optional accessory)

Optional Accessories



LED reading light



A range of custom colours are available upon request (Lead times may vary)



IV Pole



Sit-to-stand position



Seated position



Relax position



Trendelenburg & shock position



Bed/flat position



Code	CHT100300
Overall Height	1340 - 1600mm
Overall Width	890mm
Overall Depth	1300mm
Flat Position Length	2200mm
Seat Height	540 - 750mm
Seat Width	700mm
Seat Depth	700mm
Backrest Angle	-5° to 78°
Footrest Angle	0° to 90°
S.W.L	300kg