

SWEP : PAEDIATRIC WHEELCHAIR SCRIPT FORM

ACTION 3 JUNIOR

REFERRER DETAILS	CLIENT DETAILS
Referrer Name: _____	First Name: _____ Last Name: _____
Position: _____	Client DOB: ____ / ____ / ____
Phone: _____ Fax: _____	Aidacare Complex Rehab Representative: _____
Email: _____	

PLEASE FAX / EMAIL THIS FORM TO 9384 2088 OR SWEP@AIDACARE.COM.AU

ORDER EQUIPMENT DETAILS				OPTIONAL FEATURES - PLEASE COMPLETE FORM FOR ANY OPTIONAL FEATURES REQUIRED		
TYPE: Self Propelled	Width	Depth	Select	Swingaway Acrylic Tray (Includes T Armrest Support hardware)	RRP	Select
A3JRSP-EVO-1-AU	8-11"	8-10"	<input type="checkbox"/>	1513253	Tray for Seat Width 8" - 9"	\$384.50 <input type="checkbox"/>
A3JRSP-EVO-2-AU	8-11"	10-12"	<input type="checkbox"/>	1513254	Tray for Seat Width 10" - 11"	\$384.50 <input type="checkbox"/>
A3JRSP-EVO-3-AU	8-11"	12-14"	<input type="checkbox"/>	1513255	Tray for Seat Width 12" - 13"	\$384.50 <input type="checkbox"/>
A3JRSP-EVO-4-AU	12-15"	11-13"	<input type="checkbox"/>	1513256	Tray for Seat Width 14" - 15"	\$384.50 <input type="checkbox"/>
A3JRSP-EVO-5-AU	12-15"	13-15"	<input type="checkbox"/>	Elevating Legrest		RRP Select
TYPE: Transit	Width	Depth	Select	1506920	Left 330mm	\$275.00 <input type="checkbox"/>
A3JRTR-EVO-1-AU	8-11"	8-10"	<input type="checkbox"/>	1506483	Right 330mm	\$275.00 <input type="checkbox"/>
A3JRTR-EVO-2-AU	8-11"	10-12"	<input type="checkbox"/>	1513298	Left 355mm	\$275.00 <input type="checkbox"/>
A3JRTR-EVO-3-AU	8-11"	12-14"	<input type="checkbox"/>	1513297	Right 355mm	\$275.00 <input type="checkbox"/>
A3JRTR-EVO-4-AU	12-15"	11-13"	<input type="checkbox"/>	1512528-0601	Left 380mm	\$275.00 <input type="checkbox"/>
A3JRTR-EVO-5-AU	12-15"	13-15"	<input type="checkbox"/>	1512528-0601	Right 380mm	\$275.00 <input type="checkbox"/>
STANDARD FEATURES INCLUDED				80 Degree Hangers (Dark Red)		
Seating Type	Nylon w/velcro strips seat upholstery, width & depth adjustable	<input checked="" type="checkbox"/>	1428752-0607	Hangers 80 Degrees Short Right	\$57.50	<input type="checkbox"/>
Backrest Type	Fixed 7 degree backrest 3 heights - 15/16/17", delivered at 16"	<input checked="" type="checkbox"/>	1428753-0607	Hangers 80 Degrees Short Left	\$57.50	<input type="checkbox"/>
Backrest Upholstery	Tension adjustable nylon sling type	<input checked="" type="checkbox"/>	1514392-0607	Hangers 80 Degrees Long Right	\$57.50	<input type="checkbox"/>
Leg hangers	90 degree removable hanger	<input checked="" type="checkbox"/>	1514393-0607	Hangers 80 Degrees Long Left	\$57.50	<input type="checkbox"/>
Footplate Type	Flip up, angle & depth adjustable footplate	<input checked="" type="checkbox"/>	1511849	Foot Tube (x2 Required)	\$21.50	<input type="checkbox"/>
Armrest Type	Removable adjustable 15-25cm, delivered at 18cm	<input checked="" type="checkbox"/>	1434612	Tube Fixation Kit (x2 Required)	\$14.50	<input type="checkbox"/>
Frame Colour	Dark red	<input checked="" type="checkbox"/>	0572-06910	Locking block device to be fitted for each chair using 80 Degrees hangers (MUST order TWO of this code)		<input type="checkbox"/>
				Spoke Guards		
				1492138	Spoke Guard 20" Smiley	\$245.00 <input type="checkbox"/>
				1513257	Spoke Guard 20" Spider	\$245.00 <input type="checkbox"/>
				1513258	Spoke Guard 20" Flower	\$245.00 <input type="checkbox"/>
				1522862	Spoke Guard 22" Smiley	\$245.00 <input type="checkbox"/>
				1522863	Spoke Guard 22" Spider	\$245.00 <input type="checkbox"/>
				1522864	Spoke Guard 22" Flower	\$245.00 <input type="checkbox"/>



Swingaway Acrylic Tray



80 Degree Hangers



Spoke Guards

REFERRER SIGNATURE: _____

DATE: ____ / ____ / ____

Received by Aidacare: