



WORKSAFE EQUIPMENT CATALOGUE

OCTOBER 2019

VERSION: 1 | UPDATED: OCTOBER 2019



Quality Equipment
Enterprise Solutions
Home & Community





Aidacare's eTAC Ordering Portal

eTAC

Working together to achieve your goals!

- Intuitive dashboards - Full visibility of order status
- Integrated quote approval process - Includes automated notifications
- iPad App - Synchronised with offline processing
- Client Address Book - Manage your contacts
- Home Modification - Design tool

etac.aidacare.com.au

Our eTAC portal has been designed with the user in mind. It's simple to register and place orders. Everyday tasks such as placing orders, home mod / service requests and managing contacts are now at your fingertips with eTAC.

State Wide Distribution Network

AIDACARE

- | | |
|--|---|
| <p> MELBOURNE NORTH
82 Albion St,
Brunswick East VIC 3057
Tel: (03) 9981 2100
Fax: (03) 9384 2088</p> <p> MELBOURNE SOUTH
313 Warrigal Rd,
Cheltenham VIC 3192
Tel: (03) 9584 1644
Fax: (03) 9584 0562</p> <p> GEELONG
38 - 42 Mercer Street
Geelong VIC 3220
Tel: (03) 5223 2263
Fax: (03) 5223 2269</p> | <p> BENDIGO
51 Bridge St,
Bendigo VIC 3550
Tel: (03) 5441 8213
Fax: (03) 5441 7762</p> <p> MILDURA
18 Eighth Street
Mildura VIC 3500
Tel: (03) 5022 0755
Fax: (03) 5021 1845</p> <p> ALBURY
347 Urana Rd,
Lavington NSW 2641
Tel: (02) 6040 3444
Fax: (02) 6040 3222</p> |
|--|---|



SUB-CONTRACT

- | | |
|-------------------------------------|--------------------------------|
| Able Living Group - Mont Albert | Alter/VHHE |
| Able Living Group - Heidelberg West | Rehab Hire |
| Access Rehabilitation Equipment | South West Healthcare Supplies |
| Action Aids Australia | |

CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).



Product Index

Rehabilitation Equipment p 4

- How To Place An Order 4
- Equipment & Hire Orders 5
- Home Modification Orders 5
- Service & Repair Orders 6
- After Hours Repair Service 6

Bathroom & Toilet p 7

- Shower Commodes 7
- Shower Seating 11
- Toileting Aids 12
- Transfer Aids 13
- Other Bathing 14
- Bed Pans, Urinals & Bowls 16
- Grab Rails..... 17

Bedroom p 18

- Homecare Beds 18
- Foam Mattresses 20
- Alternating Air Mattresses 21
- Other Mattresses..... 23
- Foam Underlays & Overlays..... 24
- Pillows & Supports 25
- Bedside Commodes 26
- Overbed Tables 27
- Bed Sticks, Bed Poles & Other 28
- Falls Prevention & Alarms 30

Chairs p 31

- Adjustable Height Day Chairs 31
- Powerlift Recline Chairs 33
- Pressure Relief Chairs 35

Cushions & Seating p 36

- Cushions & Back Support 36
- Pressure Care Cushions 37

Daily Living Aids p 42

- Dressing Aids 42
- Kitchen Aids 44
- Gripping & Turning Aids..... 45
- Cutlery 47
- Crockery 49
- Reaching Aids 50
- General Household & Dementia Aids 51

Manual Handling p 53

- Lifters & Slings 53
- Transfer Boards & Sheets 56
- Transfer & Walk Belts 57
- Other Transfer Aids 58

Mobility Scooters p 59

- Mobility Scooters 59

Mobility Ramps p 61

- Modular & Rubber Ramps 61

Mobility Walking Frames & Aids p 62

- Seat Walkers/Rollators 62
- Walking Frames..... 64
- Walking Sticks 65
- Crutches 66
- General Accessories 67

Mobility Wheelchairs p 68

- Manual Wheelchairs 68
- Standard Power Wheelchair..... 76
- Backrests 77
- Accessories 78

Therapy & Exercise p 82

- Exercise & Resistance Training 82
- Physiotherapy 83
- Bracing Supports 84

Pressure Care & Circulation p 88

- Protection 88
- Compression 89



General Equipment Enquiries:

General Contact Details: 1300 133 120

WorkSafe Email Enquiries: worksafe@aidacare.com.au

WorkSafe Customer Service Team Members


Customer Service Manager: Michelle De Almeida

WorkSafe Team Leader: Helen Bakalis

Service Department Regular Business Hours:

Phone: 1300 133 120

Email: servicevic@aidacare.com.au





Ordering Equipment & Hire

- **Who can Order:** Clients/Claimants and Prescribing Therapists can place orders.

How to Order:

- Call (During business hours): 1300 133 120
 - Fax: (03) 9386 9170
 - Email: worksafe@aidacare.com.au.
 - Web Order: Use the WorkSafe Filter on our website to identify contracted products. Add products to cart and submit your quote.
- **Details Required when ordering:** The following Client/Claimant details are required when placing orders:
 - Type of Assistance Required: Equipment Purchase, Equipment Hire, Service, Collection, Transport, Home Modifications or Preventative Maintenance
 - Client Name
 - Claim Number
 - Contact Phone Number
 - Date of Birth
 - Date of Accident
 - Delivery Address
 - Details of Prescribing Clinician
 - The Equipment Item Number (refer to the Aidacare website or your catalogue)




Home Modification Orders



Modifications are based on an in-home assessment by a therapist and can range from something as simple as installing

grab-rails, hand-rails, lever taps and step-wedges to more complex home modifications such as installing fixed wheelchair ramps and widening doorways.


- **How to Order:** Requests can be placed by Clients/Claimants and Prescribing Therapists Fax or Email.
 - Email to: worksafe@aidacare.com.au
 - Fax to: (03) 9386 9170
- **What to send with your request:** Refer to Home Modification diagram
- **Once your order is placed:** When we receive the request via one of the methods above our Customer Service Team will coordinate the job with our Home Modification Provider who will contact the Client to arrange a time to do the Home Modification work.



HOME MODS

Our team of expert installers are on hand to meet your every needs

- Shower, Toilet & Bathroom Modifications
- Access Ramps
- Grab rails for bathrooms / toilets
- Handheld showers
- All modifications comply with the highest industry standards
- Emergency cases can be completed as soon as next day




Access Ramps




Hand Held Showers



Railings




Bath Rails




Shower Rails



Toilet Rails



Grab Rails



Stair Rails



Booking a Service/Repair

Our mobile-repair technicians can come to you to repair or maintain your equipment. Our qualified staff are industry trained, test and tag certified and Police checked.

How to request a repair:

- Call (during business hours): 1300 133 120.
- Email to: servicevic@aidacare.com.au
- Via Website: www.aidacare.com.au/WorkSafe



For Urgent After Hours Repairs: Please see the next section.

Maintaining Your Equipment: Certain equipment will automatically be scheduled for annual maintenance. You will be contacted by Aidacare’s maintenance department to book a time to conduct this maintenance.

The following items are covered by automatic maintenance:

- Electric hospital bed
- Electric standing hoist
- Electric ceiling hoist
- Electric scooter
- Motorised wheelchair
- Customised Manual wheelchair
- Tilt in space manual wheelchair
- Mobile Shower transporter/commode
- Alternating air pressure mattress
- Feed pump
- Respiratory - ventilator & suction unit (does not include CPAP or VPAP or nebulizers)

After Hours Repair Service



Aidacare offers an Urgent After-Hours Repair Service to WorkSafe clients in situations where “Client’s mobility or safety is at risk, requiring urgent equipment service or supply”:

- **Service Location:** Greater Metro Melbourne only
- How to request Urgent After Hours Repairs:
 - Call: 1300 133 128
 - Email to: serviceah@aidacare.com.au



Aidacare’s Urgent After-Hours Repair Service is available during the following times:

Aidacare Urgent After-Hours Service Availability			
Hours of Operation	Open	Closed	Phone
Weekdays	5.00pm	11.00pm	1300 133 128
Weekends	9.00am	5.00pm	
Public Holidays	9.00am	5.00pm	

Shower Commodes

Aspire Tilt in Space Shower Commode - Height Adjustable

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height
BTC066030	460mm - Attendant Propelled	1315 - 1485mm	560mm	460mm	900mm
BTC066040	530mm - Attendant Propelled	1315 - 1485mm	630mm	530mm	900mm



Shower Commodes

Description	Code	Image
-------------	------	-------

Attendant Propelled

Aspire - Aluminium Att Prop - 440mm (seat incl) - 136kg	BTC068010	
Aspire - SA Footrest - 460mm - S/Steel (seat sep) - 200kg	BTC065900	
Aspire - SA Footrest - 530mm - (seat sep) - 200kg	BTC066060	

Shower Commodes

Description	Code	Image
-------------	------	-------

Aspire - SA Footrest - 600mm - S/Steel (seat sep) - 200kg	BTC061300	
Aquatec Ocean - Att Prop - SA Footrest (130kg)	BTC083720	
Aspire - WB Platform - 460mm - S/Steel (seat sep) - Stainless Steel Castors & Safety Arms - 200kg	BTC066000	



Shower Commodes

Description	Code	Image
Aspire - WB Platform - 530mm - S/Steel (seat sep) - 200kg	BTC066050	
Aspire - WB Platform - 600mm - S/Steel (seat sep) - 400kg	BTC061400	

Self Propelled

Aspire - Aluminium Self Prop - 440mm (seat incl) - 136kg	BTC068020	
Aspire - Self Prop - SA Footrest - 460mm (seat sep) - 200kg	BTC066100	
Aspire - Self Prop - SA Footrest - 530mm - (seat sep) - 200kg	BTC066065	
Aspire - Self Prop - SA Footrest - 600mm - S/Steel (seat sep) - 200kg	BTC061305	
Aquatec Ocean - Self Prop - SA Footrest - 480mm - 150kg	BTC083730	
Aspire - Self Prop - WB Platform - 460mm (seat sep) - 200kg	BTC066200	
Aspire - Self Prop - WB Platform - 530mm (seat sep) - 200kg	BTC066055	

Tilt in Space

Tilt in Space w/ polyurethane back (seat sep)	BTC069500	
---	-----------	--

Shower Commodes

Description	Code	Image
Aquatec Ocean VIP - Tilt (150kg)	BTC083710	
Aquatec Ocean VIP -Dual - Adj	BTC083715	

Self Propelled

Aspire Open 460mm	BTC066600	
Aspire Open 530mm	BTC066750	
Aspire Open 600mm	BTC066800	
Polyurethane - Open front seat - fitted with pan carrier runners (for TIS Commode BTC069500)	BTC069510	
Polyurethane - Closed front seat - fitted with pan carrier runners (for TIS Commode BTC069500)	BTC069520	
Aspire - Pressure Reducing - 450x450 - Open Front	BTC089130	
Aspire - Pressure Reducing - 450x450 - Closed Front	BTC089140	
Aspire - Pressure Reducing - 450x450 - Open Left	BTC089150	
Aspire - Pressure Reducing - 450x450 - Open Right	BTC089160	



Shower Commodes

Description	Code	Image
Invacare - Universal Soft - for Aquatec Ocean Shower Commodes	BTC087640	
Invacare - for Aquatec Shower Commodes - Ergo Hygiene - Recess Seat	BTC087642	

Self Propelled

Aspire - Self Propelling Kit - 22in Wheel	BTC082900	
Aspire - Pan Carrier - (410mm wide) - For Aspire 460/530/600mm Commodes	BTC066400	
Aspire - Pan Carrier - (350mm wide) - For Aspire 410mm Commodes	BTC066300	
Aspire - Bed Pan - Standard	BTT151800	
Aspire Commode Bowl with Lid & Handle	BTC081610	
Aspire - Elevating Legrest - Left	BTC067110	
Aspire - Elevating Legrest - Right	BTC067115	
Aspire - Stump Support (Interchangeable) - Clamp On	BTC067105	
Aspire - Advanced Headrest (X-Shaped)	BTC067050	

Shower Commodes

Description	Code	Image
Aspire - Oxygen Bottle Carrier	BTC067200	
Aspire - Pelvic Strap - Adjustable	BTC067100	
Aspire - Adjustable IV Pole	BTC067103	
Aspire - Armrest - Padded Offset	BTC067000	
Invacare - Self Propelling Kit for Aquatec Commode	BTC083716	
Invacare - Safety Bar - for Aquatec Ocean & VIP Shower Commodes	BTC087606	
Invacare - Safety Bar - for Aquatec Ocean XI Shower Commodes	BTC087608	
Insert - Invacare - for Aquatec Shower Commodes - Recess Cover - Soft Seat	BTC087648	
Calf Support - Invacare - Below Knee - for Aquatec Ocean Shower Commodes - Left or Right	BTC087612	
Calf Support - Invacare - Above Knee - for Aquatec Ocean XL Shower Commodes - Left or Right	BTC087614	
Calf Support - Invacare - Above Knee - for Aquatec Ocean VIP Shower Commodes - Left	BTC087616	



Shower Commodes

Description	Code	Image
Calf Support - Invacare - Above Knee - for Aquatec Ocean VIP Shower Commodes - Right	BTC087618	
Calf Support - Invacare - Below Knee - for Aquatec Ocean VIP Shower Commodes - Left	BTC087620	
Calf Support - Invacare - Below Knee - for Aquatec Ocean VIP Shower Commodes - Right	BTC087622	
Armrest Lock - Invacare - for Aquatec Ocean VIP Shower Commodes	BTC087624	
Armrest Lock - Invacare - for Aquatec Ocean Shower Commodes	BTC087626	
Armrest - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087628	
Chest Belt - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087630	
Pelvic Belt - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087632	
Anti Tippers - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087636	
Plat Pin - Invacare - for Aquatec Shower Commodes - Height Adjustable Footrest	BTC087638	

Shower Commodes

Description	Code	Image
Splash Guard - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087650	
Safeguard - Invacare - for Aquatec Shower Commodes - S/Steel	BTC087656	
Frame Clip - Invacare - for Aquatec Ocean & VIP Shower Commodes	BTC087634	
Frame Clip - Invacare - for Aquatec Shower Commodes	BTC087644	
Footrest Clip - Invacare - for Aquatec Shower Commodes	BTC087646	



Product brochures also available!
Please contact us for a copy



Shower Seating


Shower Seating

Description	Code	Image
-------------	------	-------

Shower Chairs







Shower Chair w/ Arms - Aspire - 470mm - Treated Steel - 200kg	BTS118000	
Shower Chair w/ Arms - Aspire - 470mm - Aluminium - 200kg	BTS118005	
Shower Chair w/ Arms - Aspire - Wide - 525mm - Treated Steel - 200kg	BTS118010	
Shower Chair w/ Arms - Aspire - Wide - 525mm - Aluminium - 200kg	BTS118015	
Shower Chair w/ Padded Seat, Arms & Back - Aspire MAXI - Height Adjustable - Width Adjustable (550mm to 750mm) - Aluminium - 300kg	BTS118020	
Shower Chair w/ Arms - Invacare - Aquatec Pluto	BTS108611	
Shower Chair - Invacare - Aquatec Pico (Arms and Back independently removable) - 160kg	BTS108612	

Shower Stools

Shower Stool w/ Arms - Aspire - 495mm - Treated Steel - 175kg	BTS119000	
---	-----------	---

Shower Seating

Description	Code	Image
-------------	------	-------

Shower Stool w/ Arms - Aspire - 495mm - Aluminium - 175kg	BTS119005	
Shower Stool w/ Arms - Dayz - PVC Strapping Seat - Zinc**	BTS113600	
Shower Stool w/ Arms & Padded Seat & Back - Aspire - 495mm - Treated Steel - 175kg	BTS119025	
Shower Stool w/ Arms & Padded Seat - Aspire - 495mm - Treated Steel - 175kg	BTS119020	
Shower Stool w/ Padded Seat & Arms - Aspire MAXI (550mm to 750mm) - Aluminium - 300kg	BTS119030	
Shower Stool w/ Arms - Aspire - Wide - 525mm - Treated Steel - 175kg	BTS119010	
Shower Stool w/ Arms - Aspire - Wide - 525mm - Aluminium - 175kg	BTS119015	
Clip on Seat - Aspire Shower Stool	BTS119100	
Clip on Seat - Aspire Wide Shower Stool	BTS119105	
Clip on Padded Seat - Aspire Shower Stool	BTS119110	
Shower Stool w/ Arms - Clip on Padded Backrest only - Aspire - White	BTS119115	

Toileting Aids



Toileting Aids

Description	Code	Image
-------------	------	-------

Over Toilet Aids

Over Toilet Aid - Aspire - 460mm - Treated Steel - 175kg	BTT145000	
Over Toilet Aid - Aspire - 460mm - Aluminium - 175kg	BTT145005	
Over Toilet Aid - Aspire - Wide - 525mm - Treated Steel - 175kg	BTT145010	
Over Toilet Aid - Aspire - Wide - 525mm - Aluminium - 175kg	BTT145015	
Over Toilet Aid - Aspire MAXI -Extra Wide (550mm to 750mm) - Aluminium - 300kg	BTT145020	
Over Toilet Aid - Padded Toilet Seat Option - Aspire - Grey	BTT145110	
Over Toilet Aid - Plastic Toilet Seat Replacement - Aspire - White	BTT145100	

Toileting Aids


Description	Code	Image
-------------	------	-------

Over Toilet Aid - Plastic Splash Guard Replacement - Aspire - White	BTT145105	
Over Toilet Aid - U Frame - Adjustable	BTT156700	
Toilet Surround - Aspire - Height Adjustable 450mm - Aluminium - 175kg	BTT155100	

Aspire Toilet Seat Raisers - 200kg

50mm - With Lid	BTT149900	
100mm - With Lid	BTT150100	
150mm - With Lid	BTT150400	
50mm - No Lid	BTT149800	
100mm - No Lid	BTT150200	
150mm - No Lid	BTT150300	

Toilet Seat Raiser w/Arms

Raised Toilet Seat - w/ Swing Arms	BTT150450	
---------------------------------------	-----------	---



Transfer Aids

Transfer Aids

Description	Code	Image
-------------	------	-------

Bath / Shower Boards

Bath/Shower Board - Plastic w/ Handle & Adj Width **	BTS113700	
Bath/Shower Board - Regular	BTS111900	
Bath/Shower Board - Large	BTS111700	
Swivel Disk - Aquatec Disk - 130kg	BTS108603	
Swivel Transfer Board & Disk - Aquatec Trans - 130kg	BTS108604	

Bath Chairs

Swivel Bather - 710mm - 100kg	BTS106500	
Myco Swivel Bather - Corner	BTS106700	
Myco Ultra Swivel Bather - 535mm - 160kg	BTS106540	
Acc - Myco Spacer 3in (75mm)	BTS106800	
Acc - Myco Spacer 4in (100mm)	BTS106950	
Acc - Myco Spacer 6in (150mm)	BTS106900	

Transfer Aids

Description	Code	Image
-------------	------	-------

Shower Chair - Myco Swivel - 445mm Wide - 100kg Weight **	BTS107100	
Bath Lift - Aquatec Orca 170kg	BTS112251	

Bath Transfer Bench

Bath Transfer Bench - w/ Arm & Padded Seat & Backrest	BTS110900	
Bath Transfer Bench - w/ Arm & Plastic Seat & Backrest	BTS108200	
Bath Transfer Bench - Bariatric - w/ Arm & Plastic Seat & Backrest - 227kg	BTS105300	
Bath Transfer Bench - Leg Suction Cup - Screw on - for BTS110900, 108200	BTS108210	
Bath Transfer Bench - w/ Sliding Seat	BTS128600	

Other Bathing



Other Bathing

Description	Code	Image
Hand Held Shower		
Push On - Single Outlet - 1.25m - Plastic	BTS114400	
Push On - Single Outlet - 2.0m - Plastic	BTS114500	
Push On - Dual Outlet - 1.25m - Plastic	BTS114200	
Push On - Dual Outlet - 2.0m - Plastic	BTS114300	
1.5m Hose - Chrome Plated	BTS107300	
2.0m Hose - Chrome Plated	BTS107200	
Handheld Shower Set For low pressure system - 1.5m Chrome Plated Hose	BTS107495	
On Rail - 1.5m Hose - Chrome Plated	BTS107500	
On Rail - 2.0m Hose - Chrome Plated	BTS107490	
Component - Switchcock water diverter - Chrome plated	BTS107410	

Other Bathing

Description	Code	Image
Component - Non return backflow valve - Chrome plated	BTS107420	
Shower Curtain		
180cm L x 100cm W	BTS123100	
180cm L x 130cm W	BTS123300	
180cm L x 180cm W	BTS123400	
180cm L x 200cm W	BTS123500	
200cm L x 200cm W	BTS123600	
180cm L x 240cm W	BTS123700	
180cm L x 260cm W	BTS123750	
200cm L x 100cm W	BTS123760	
200cm L x 130cm W	BTS123770	
200cm L x 180cm W	BTS123780	
200cm L x 240cm W	BTS123800	
200cm L x 260cm W	BTS123820	
Shower Mat - 530x530mm	BTS109600	
Non-Slip Strips - 50mm wide (per meter)	YCH865100	



Other Bathing



Description	Code	Image
-------------	------	-------

Limb/Cast Protection






Adult Below Knee - X Sml	PTA534100	
Adult Below Knee - Lge	PTA534100AL	
Adult Below Knee - Sml	PTA534100AS	
Adult Below Knee - X Sml	PTA534100AXS	
Adult Below Knee - Lge	PTA534100SL	
Adult Below Knee - Sml	PTA534100SS	
Adult Below Knee - X Sml	PTA534100SXS	
Adult Below Knee - Lge	PTA534100TL	
Adult Below Knee - Sml	PTA534100TS	
Limb/Cast Protection - Limbo - Adult Below Knee - Waterproof - X Short - Sml	PTA534100XSS	
Limb/Cast Protection - Limbo - Adult Hand - Waterproof	PTA534101	
Adult Below Elbow - X Sml	PTA534102	
Adult Below Elbow - Lge	PTA534102L	
Adult Below Elbow - Sml	PTA534102S	
Adult Full Arm - Sml	PTA534103	
Adult Full Arm - Lge	PTA534103L	
Adult Elbow - Sml	PTA534104	
Adult Elbow - Lge	PTA534104L	

Other Bathing

Description	Code	Image
-------------	------	-------

Adult Full Leg - Short - Sml	PTA534105	
Adult Full Leg - Average - Lge	PTA534105AL	
Adult Full Leg - Average - Sml	PTA534105AS	
Adult Full Leg - Short - Lge	PTA534105SL	
Adult Full Leg - Tall - Lge	PTA534105TL	
Cast Protector / Wound Cover - Derma-Plas - Water Resistant	FBS281190	
Cast Protector - Long Leg	BTS109400	
Cast Protector - Short Leg	BTS109300	
Cast Protector - Foot/Ankle	BTS109500	

Washers & Wipes

Bath Sponge - Long Handled - Plastic Handle	BTS111100	
Bath Sponge - Long Handled - Wire Handle **	BTS112800	
Toe Washer **	BTS113000	
Bottom Wiper - Standard **	BTT150600	
Bottom Wiper - Folding	BTT150700	

Bed Pans, Urinals & Bowls



Bed Pans, Urinals & Bowls

Description	Code	Image
Urinal - Male with Screw Lid	BTT148400	
Bowl & Lid - Aspire - Classic/Deluxe Bedside & Shower Commodes	BTC081610	
Urinal - Spillproof - Female Adapter (adapter only)	BTT148300	
Urinal - Female	BTT148700	
Urinal - Female - Cygnet Design	BTT148800	

Bed Pans, Urinals & Bowls

Description	Code	Image
Urinal - Male with Lid	BTT148450	
Urinal - Male Clear **	BTT160900	
Urinal - Spillproof - Male**	BTT148200	
Urinal Bottle Holder - Nylon Coated **	BTT149000	



Grab Rails

Grab Rails

Description	Code	Image
-------------	------	-------

Concealed Flange 32mm

Brushed S/Steel		
Peened S/Steel		
White		
White ripple		
Almond ivory		
Almond ivory ripple		

Modular Tube - Grab Rail

Brushed S/Steel - Peened S/Steel - White - White ripple - Almond ivory - Almond ivory ripple	
--	--

Modular Tube

Brushed S/Steel	
40° Elbow RH - Brushed S/Steel	
40° Elbow LH - Brushed S/Steel	
90° Elbow RH - Brushed S/Steel	
90° Elbow LH - Brushed S/Steel	
Internal 90° - Brushed S/Steel	
90° vertical RH - Brushed S/Steel - White - Almond ivory	
90° vertical LH - Brushed S/Steel - White - Almond ivory	

Grab Rails

Description	Code	Image
-------------	------	-------

90° vertical centred LH - Brushed S/Steel- White - Almond ivory	
90° vertical centred RH - Brushed S/Steel - White - Almond ivory	
Elbow Ends (Pair) - Brushed S/Steel - White - Almond ivory	
Straight joiner - Brushed S/Steel - White - Almond ivory	
Straight end - Brushed S/Steel - White - Almond ivory	
Tee section - Brushed S/Steel - White - Almond ivory	
Connector Ring	
Hand rail support - Brushed S/Steel - White - Almond ivory	

Throne Toilet Grab Rails

3 in 1 - Powder Coated - S/Steel Wrap Around - S/Steel	
Spacer, Spacer Bracket - Powder Coated	

Bath Rails

Bath Rail - Clamp On	
Suction Grab Rail	

Homecare Beds



Aspire Premium Homecare Bed



Code	Description	External Length	External Width	Bed Height	Mattress Size
BEB046050K	Premium Homecare Bed - 4 Function	2140 - 2240 mm	1025 mm	220 - 710 mm	2000 x 900 mm



Aspire Community Care Bed




Code	Description	External Length	External Width	Bed Height	Mattress Size
BEB046050	Community Care Bed - 4 Function	2140 - 2240 mm	1025 mm	220 - 710 mm	2000 x 900 mm





Aidacare AC Bed Series

Code	Description	External Length	External Width	Bed Height	Mattress Size
BEB045900	AC3 Single Bed - 3 Function	2120 mm	920 mm	285 - 685 mm	2000 x 900 mm
BEB045905	AC3 King Single Bed - 3 Function	2140 mm	1090 mm	285 - 685 mm	2030 x 1050 mm



Contego Floorline Bed - Folding

Code	Description	Total Length	Total Width	Total Height
BEB049500	Contego Floorline Bed - Single	2250 mm	920 mm	110 - 710 mm



Homecare Beds

Description	Code	Image
-------------	------	-------

Aspire Community Care Bed

Side Rails - 3/4 Length Clip on	BEA016585	
Side Rails - Full Length Wooden	BEA029865	

Aidacare AC Bed

Bed Side Rails - AC3 Single & King-Single Beds	BEA029700	
Self Help Pole - AC3 Single & King-Single Beds	BEA029705	

Etude Plus - 180kg

Single 4 Section pull apart	BEB045360	
Folding Metal	BEA017600	

Homecare Beds


Description	Code	Image
-------------	------	-------

Bed Self Help Pole	BEA017900	
--------------------	-----------	--

Bed Rails

Bed Rail - Aspire - Fixed	BEA013600	
Bed Side Rails - Universal Folding - Aidacare AC/Ward Beds	BEA013700	
Bed Rail Protector - Aspire - 1260x670mm	BEA010775	
Bed Rail Protector - Aspire - Large - 1580x980mm	BEA010776	

Foam Mattresses



Foam Mattresses

Description	Code	Image
-------------	------	-------

Lifecomfort Mattresses

Aidacare Lifecomfort - AC Series - 2000 x 900 x130mm - 120kg	BEM046410	
Aidacare Lifecomfort - King Single 2030 x 1050 x130mm - 120kg	BEM046460	
Aidacare Lifecomfort - with Safety Sides - 1940x900x130mm - 120kg	BEM046475	
Aidacare Lifecomfort Deluxe - AC Series - 2000x900x150mm -150kg	BEM046510	
Aidacare Lifecomfort Deluxe - King Single 2030x1050x150mm -150kg	BEM046455	
Aidacare Lifecomfort Pressure Reducing - AC Series - 2000x900x150mm -150kg	BEM046610	

Foam Mattresses

Description	Code	Image
-------------	------	-------


Aidacare Lifecomfort Pressure Reducing - King Single - 2030x1050x150mm -150kg	BEM046660	
Aidacare Lifecomfort Bariatric - King Single - 2030x1050x200mm -300kg	BEM046700	
Aidacare Lifecomfort Pressure Reducing - w/ Roho Cutout - 1940x900x150mm -150kg	BEM046720	
Foam - Invacare - MSS Softform Maxi Glide Mattress (197x88x15.2cm)	PTM636380	
Aidacare Lifecomfort Acute Needs - Castellated Foam Top w/ Side-walls - 1980x880x150mm - 250kg	BEM046615	
Aidacare Lifecomfort Acute Needs - Castellated Foam Top w/ Side-walls - King Single - 2030x1050x150mm -250kg	BEM046616	



Alternating Air Mattresses



Aspire Lifecomfort Hybrid Mattress - Air & Foam Combination

Code	Description	Overall Length	Overall Width	Overall Height	Weight
PTM638010	Hybrid Single Mattress	1980 mm	895 mm	170 mm	20 kg
PTM638030	Hybrid Mattress Pump	320 mm	130 mm	230 mm	3.5 kg



Aspire Lifecomfort Hybrid Acute Mattress - Air & Foam Combination

Code	Description	Overall Length	Overall Width	Overall Height	Weight
PTM638020	Hybrid Acute Single Mattress	1980 mm	895 mm	170 mm	20 kg
PTM638030	Hybrid Mattress Pump	320 mm	130 mm	230 mm	3.5 kg





Alternating Air Mattresses

Description	Code	Image
Mattress Overlay - AirLo3 - Alt Air - 21cell - 2000x900x75mm - 120kg **	PTM636629	
Mattress Overlay -AirLo5 - Alt Air - 18cell - 2000x900x130mm - 150kg **	PTM636630	
Mattress Overlay - AirLo 5 King Single - Alt Air - 18cell - 2000x1070x130mm - 170kg **	PTM636631	

AirLo Mattress Replacements

Mattress - AirLo 8 - Alt Air - 20cell - 2000x900x205mm - 200kg **	PTM636635	
Mattress - AirLo 8 King Single - Alt Air - 20cell - 2000x1070x205mm - 200kg **	PTM636636	
AirLo - Thermal Hyperstretch Microfibre Coverlet	PTM636625	
AirLo - Thermal Hyperstretch Microfibre Coverlet - King Single	PTM636626	

Ultimate Overlays

Mattress Overlay - Ultimate 125 Digital - Alt Air - 18x18cm cell - 2000x880x130mm **	PTM636530	
Mattress Overlay - Ultimate 125 Digital - Alt Air - 18x18cm cell - 2000x880x130mm **	PTM636540	
Mattress - Ultimate 225 Digital - Alt Air - 21x27cm cell - 2000x880x275mm **	PTM636545	

Alternating Air Mattresses

Description	Code	Image
Mattress - Theraflow - Alt Air - Single - 2000x880x245mm - 200kg	PTM636640	
Mattress - Theraflow - Alt Air - King Single - 2000x1070x245mm - 200kg	PTM636645	
Mattress - Theraflow - Alt Air - Bariatric - 2000x1070x245mm - 350kg	PTM636646	



Product brochures also available!
Please contact us for a copy





Other Mattresses

Other Mattresses

Description	Code	Image
Mattress - Alt Air - Invacare - Softform Premier Active Mattress (197x88x15.2cm)	PTM636395	
Mattress Pump - Alt Air - Invacare Premier Active Pump Only	PTM636396	
Mattress Section - ROHO - Foam Insert Pad 21x34in	PTM636220	
Mattress Overlay - ROHO Levelling Pad - 34inx 80in	PTM636360	

Other Mattresses

Description	Code	Image
Mattress - ROHO (1 section) - 20inx34in **	PTM636320	
Mattress Overlay - ROHO (4 sections of 20x34in) complete with cover	PTM636335	
Mattress - ROHO - Mattress Attachment 4 Cell Wide Bariatric	PTM636700	


Foam Underlays & Overlays



Foam Underlays & Overlays

Description	Code	Image
-------------	------	-------

Underlays

Pressure Safety Underlay - Aidacare Lifecomfort - 2000x900x50mm - 135kg	BEM047050	
--	-----------	---



Foam Underlays & Overlays

Description	Code	Image
-------------	------	-------

Memory Foam - Naturelle Latex

Single - 1808x915x50mm	PTM636260	
Double - 1808x1370x50mm	PTM636270	
Queen - 2030x1520x50mm	PTM636280	

Overlays - Memory Foam

Topper Single	PTM636290	
Topper King Single	PTM636291	
Topper Double	PTM636292	
Topper Queen	PTM636293	
Topper King	PTM636294	
Naturelle Latex - Single - 1880x905x50mm - 100kg	BEM047181	
Naturelle Latex - Double - 1880x1355x50mm - 100kg	BEM047182	
Naturelle Latex - Queen - 2000x1505x50mm - 100kg	BEM047183	
Naturelle Latex - King - 2000x1805x50mm - 100kg	BEM047184	



Introducing our Aspire
Lifter & Sling Product Video



To view, visit the link below

<https://www.aidacare.com.au/knowledge/videos/>



Pillows & Supports




Pillows and Supports

Description	Code	Image
-------------	------	-------

Pillows and Supports

Description	Code	Image
-------------	------	-------

Pillows

Tranquillow	BEP053000	
Pillow - Tranquillow - Standard Deluxe	BEP053600	
Tranquillow - Small Soft	BEP053500	
Complete Sleeprrr	BEP055100	
Complete Sleeprrr Plus	BEP055400	
Family Pillow - Low Profile	BEP052300	
Family Pillow- Medium Profile	BEP052400	

Wedge Supports


Pyramid w/ Tailored Slip	BEP054200	
Leg Relaxer - w/ Quilted Slip	BEP056200	
Leg Relaxer - w/ Waterproof Slip	BEP056250	
Positioning Body Wedge	BEP055905	
BedWedge - 75x68x25cm - with Quilted Slip **	BEP054800	
BedWedge - 75x68x25cm - with Waterproof Slip	BEP055000	
Bed Backrest - Adjustable Deluxe	BEA009900	
Bed Backrest - Adjustable (Pillow Support)	BEA014700	

Bedside Commodes





Aspire Deluxe Bedside Commode

Code	Description	Overall Height	Overall Width	Seat Depth	Backrest Height
BEC045115	Champagne Vinyl	865 - 1015 mm	600 mm	450 mm	440 mm
BEC045116	Slate Vinyl	865 - 1015 mm	600 mm	450 mm	440 mm



Aspire Classic Bedside Commode

Code	Description	Overall Height	Overall Width	Seat Depth	Backrest Height
BEC045005	Champagne Vinyl	900 - 1040 mm	560 mm	440 mm	450 mm
BEC045006	Slate Vinyl	900 - 1040 mm	560 mm	440 mm	450 mm





Overbed Tables

Aspire Overbed Table

Code	Description	Height	Width	Length
BEA028000	Aspire Overbed Table Thermoform Recessed Top - Cream	730 - 1095 mm	430 mm	800 mm
BEA006700	Aspire Overbed Table Laminate Flat Top - White	730 - 1095 mm	430 mm	800 mm
BEA006701	Aspire Overbed Table Laminate Flat Top - Beech	730 - 1095 mm	430 mm	800 mm
BEA028100	Aspire Overbed Table Split Tilt Top	780 - 1140 mm	430 mm	800 mm



Thermoform Recessed Top



Laminate Flat Top - White



Laminate Flat Top - Beech



Split Tilt Top








Bed Sticks, Bed Poles & Other



Bed Sticks, Bed Poles & Other

Description	Code	Image
-------------	------	-------

Bed Poles

Bed Pole - Aspire - Single Clamp On - Folding - Right	BEA009300	
Bed Pole - Aspire - Single Clamp On - Folding - Left	BEA009310	
Bed Pole - Aspire - Adjustable Crook Handle	BEA014150	
Bed Pole - Lawson	BEA019785	
Bed Pole - Townsend - Tall - Clamp On (SA ADEonly)	BEA029536	

Bed Sticks, Bed Poles & Other

Description	Code	Image
-------------	------	-------

Bed Cradles

Bed Cradle - Aspire - Fixed - Treated Steel	BEA014605	
Bed Cradle - Aspire - Height Adjustable - Treated Steel	BEA014610	

Self Help Poles

Bed Self Help Pole - Aspire - Free Standing	BEA013100	
Bed Self Help Pole - Aspire - Handle only	BEA013400	

Security Poles

Security Pole and Curve Grab Bar - Black	BEA821860	
--	-----------	---






Bed Sticks, Bed Poles & Other

Description	Code	Image
Security Pole and Curve Grab Bar - White	BEA821870	
Security Pole - White	BEA821880	

Raisers

Blocks - 40mm Grey	BEA007400	
Blocks - 100mm Grey	BEA007401	

Bed Sticks, Bed Poles & Other

Description	Code	Image
Bed / Chair Raiser - Cones - 90mm (Pkt 4)	BEA009500	
Set - 140mm Raise, Plastic Elephant Feet, x4	BEA009600	
Bed / Chair Raiser - Blocks - Wooden 75mm (Pkt 4)	BEA009700	



Product brochures also available!
Please contact us for a copy



Falls Prevention & Alarms



Aidacare Lifecomfort Fall Safety Mat

Code	Description	Length	Width	Height
BEA008550	Fall Safety Mat	1750 mm	600 mm	50 mm





Adjustable Height Day Chairs

Aspire Adjustable High Back Day Chairs

Code	Description	Overall Height	Overall Width	Overall Depth	Backrest Height
CHP208025	Aspire Adjustable Day Chair - Latte Vinyl	1070 - 1220 mm	575 - 675 mm	760 - 860 mm	740 - 790 mm
CHP208030	Aspire Adjustable Day Chair - Ink Vinyl	1070 - 1220 mm	575 - 675 mm	760 - 860 mm	740 - 790 mm



Adjustable Height Day Chairs

Description	Code	Image
-------------	------	-------

Lowback Chairs

Chair Lowback - Aspire Classic Day Chair - Champ Vinyl - 460mm - 160kg	CHP208000	
Chair Lowback - Aspire Classic Day Chair - Slate Vinyl - 460mm - 160kg	CHP208005	
Chair Lowback - Harbour Bridge Chair - 450mm	CHP192700	

Adjustable Height Day Chairs

Description	Code	Image
-------------	------	-------

Highback Chairs

Chair Highback - Aspire Classic Day Chair - Champ Vinyl - 490mm - 160kg	CHP208050	
Chair Highback - Aspire Classic Day Chair - Slate Vinyl - 490mm - 160kg	CHP208055	



Adjustable Height Day Chairs

Description	Code	Image
-------------	------	-------




Adjustable Height Day Chairs

Description	Code	Image
-------------	------	-------

Waterfall Highback Chairs

Chair Highback - Aspire Waterfall Day Chair - Ink Vinyl - 535mm - 200kg	CHP208070	
Chair Highback - Aspire Waterfall Day Chair - Latte Vinyl - 535mm - 200kg	CHP208071	
Chair Highback - Aspire Waterfall Day Chair - Oatmeal Fab - 535mm - 200kg	CHP208075	
Chair Highback - Aspire Waterfall Day Chair - Mocha Fab - 535mm - 200kg	CHP208076	

Accessories

Chair Footrest - Adjustable - Padded - Aspire - Champagne Vinyl	CHP208080	
Chair Footrest - Adjustable - Padded - Aspire - Slate Vinyl	CHP208085	
Chair Rear Legs w/ Wheels - Adjustable - for Aspire Highback Day Chairs	CHP208090	



Aidacare have a range of Product Application Guides available for the following areas:






Powerlift Recline Chairs

Aspire Signature Lift Recline Chair - Space Saver

Code	Description	Overall Height	Overall Width	Overall Depth	Seat Width
CHP197700	Space Saver - Small - Mink	1010 mm	780 mm	690 mm	450 mm
CHP197705	Space Saver - Small - Lagoon	1010 mm	780 mm	690 mm	450 mm
CHP197710	Space Saver - Medium - Mink	1050 mm	780 mm	715 mm	480 mm
CHP197715	Space Saver - Medium - Lagoon	1050 mm	780 mm	715 mm	480 mm
CHP197716	Space Saver - Medium - Heritage Vinyl	1050 mm	780 mm	715 mm	480 mm
CHP197720	Space Saver - Large - Mink	1110 mm	780 mm	740 mm	500 mm
CHP197725	Space Saver - Large - Lagoon	1110 mm	780 mm	740 mm	500 mm
CHP197730	Space Saver - Large - Heritage Vinyl	1110 mm	780 mm	740 mm	500 mm



Aspire Signature Lift Recline Chair - Dual Action - Independent Leg Raise & Back

Code	Description	Overall Height	Overall Width	Overall Depth	Seat Width
CHP197735	Dual Action - Small - Mink	1010 mm	780 mm	710 mm	450 mm
CHP197740	Dual Action - Small - Lagoon	1010 mm	780 mm	710 mm	450 mm
CHP197745	Dual Action - Medium - Mink	1050 mm	780 mm	735 mm	480 mm
CHP197750	Dual Action - Medium - Lagoon	1050 mm	780 mm	735 mm	480 mm
CHP197751	Dual Action - Medium - Heritage Vinyl	1050 mm	780 mm	735 mm	480 mm
CHP197755	Dual Action - Large - Mink	1110 mm	780 mm	760 mm	500 mm
CHP197760	Dual Action - Large - Lagoon	1110 mm	780 mm	760 mm	500 mm
CHP197730	Space Saver - Large - Heritage Vinyl	1110 mm	780 mm	740 mm	500 mm





Powerlift Recline Chairs

Description	Code	Image
-------------	------	-------

Powerlift Recline Chairs

Description	Code	Image
-------------	------	-------

Chairs

Chair Powerlift Recline - Pride - C1 (Arctic Blue) - 460x510x660mm - 1 Motor **	CHP197209	
Chair Powerlift Recline - Pride - C1 (Chocolate) - 460x510x660mm - 1 Motor **	CHP197211	
Chair Powerlift Recline - Pride - 560 (Spruce) - 480x480x710mm - 1 Motor**	CHP193302	
Chair Powerlift Recline - Pride - 560 (Burgundy) - 480x480x710mm - 1 Motor**	CHP193304	
Chair Powerlift Recline - Pride - C5 (Chocolate) - 560x460x710mm - 1 Motor **	CHP197102	
Chair Powerlift Recline - Pride - C5 (Arctic Blue) - 560x460x710mm - 1 Motor **	CHP197104	
Chair Powerlift Recline - Pride - C101 (Sandal) - 510x510x700mm - 1 Motor**	CHP197213	
Chair Powerlift Recline - Pride - C101 (Cocoa) - 510x510x700mm - 1 Motor**	CHP197214	
Chair Powerlift Recline - Pride - D30 (Latte) - 560x560x740mm - 1 Motor **	CHP198002	
Chair Powerlift Recline - Pride - D30 (Spruce) - 560x560x740mm - 1 Motor **	CHP198004	
Chair Powerlift Recline - Pride - D30 (Royal Blue Vinyl) - 560x560x740mm - 1 Motor **	CHP198006	

Chair Powerlift Recline - Top Form - Ashley - Maxi 530x560x720mm - 140kg - 1 Motor	CHP195300	
Chair Powerlift Recline - Top Form - Ashley - LGE - 670x600x700mm - 170kg - 1 Motor	CHP190400	
Chair Powerlift Recline - Top Form - Ashley - XLGE - 860x650x720mm - 230kg - 1 Motor	CHP190500	
Chair Powerlift Recline - Pride - LC107 (Pecan) - 480x510x700mm - 2 Motor - 147kg	CHP197230	
Chair Powerlift Recline - Pride - LC107 (Burgundy) - 480x510x700mm - 2 Motor - 147kg	CHP197235	
Chair Powerlift Recline - Pride - LC107 (Ocean) - 480x510x700mm - 2 Motor - 147kg	CHP197240	
Chair Powerlift Recline - Pride - C6 (Chocolate) - 560x530x710mm - 2 Motor - 136kg SWL	CHP197202	
Chair Powerlift Recline - Pride - C6 (Arctic Blue) - 560x530x710mm - 2 Motor - 136kg SWL	CHP197204	
Chair Powerlift Recline - Pride - C6 (Royal Blue Vinyl) - 560x530x710mm - 2 Motor - 136kg SWL	CHP197206	

Head/Armrest Cover Set

Aspire Signature Recliners - Mink	CHP197770	
Aspire Signature Recliners - Lagoon	CHP197771	
Aspire Signature Recliners - Heritage	CHP197772	



Pressure Relief Chairs

Aspire Air Lift Chair

Code	Description	Seat Depth	Seat Width	Backrest Height	Seat Height
CHP198865	Space Saver - Petite	460 mm	450 mm	610 mm	480 mm
CHP198870	Space Saver - Large	510 mm	510 mm	610 mm	480 mm
CHP198875	Dual Action - Petite	460 mm	450 mm	610 mm	480 mm
CHP198880	Dual Action - Large	510 mm	510 mm	610 mm	480 mm



Cushions & Back Support



Cushions & Back Support

Description	Code	Image
-------------	------	-------

Lumbar Rolls

D Shape	GSS340240	
Demi D (small D shape)	GSS340248	
D Shape - Memory Foam	GSS340241	
Oval Shape	GSS340242	
Egg Shape	GSS340246	
Round	GSS340244	

Back Support

Back Hugger - Memory Foam	GSS340201	
Back Hugger - Durofab	GSS340200	
Contoured Large	GSS340217	
Contoured Medium	GSS340216	
Contoured Small	GSS340215	
Flobac Grey	GSS340230	
Flobac Grey with Beaded Seat	GSS340231	
Cushion - Back Support - Backform	GSS340205	
Total Spine Support Tall (Car/Office Chair)	GSS340220	
Total Spine Support Tall (Car/Office Chair)	GSS340221	

Cushions & Back Support

Description	Code	Image
-------------	------	-------

Orthopod - Lower Back	FBS285410	
Obus - Forme - Low	FBS285430	
Obus - Forme - Wide	FBS285440	
Therapod - Hight Adjustable	FBS285610	
Therapod - Cervical Pad	FBS285620	

Relief Ring Cushions

Cushion - Relief Ring **	GSS340150	
Cushion - Ring Silicone Fibre	PTC542250	

Tailbone Cushions

Cushion - Coccyx Tailbone Wedge	GSS340135	
Cushion - Coccyx Tailbone w/ W/Proof cover	GSS340136	

Keyhole Cushions

Cushion - Keyhole	GSS340130	
-------------------	-----------	--


Pressure Care Cushions



Pressure Care Cushions

Description	Code	Image
-------------	------	-------


JAY Basic Cushions

14x16in	PTC567600	
15x16in	PTC568400	
16x16in	PTC567400	
16x18in	PTC567500	
17x16in	PTC568450	
18x16in	PTC566800	
18x18in	PTC566900	
20x18in	PTC567700	
22x18in	PTC567750	
22x20in	PTC566950	

JAY Lite Cushions


14x14in	PTC571902	
14x16in	PTC571901	
15x15in	PTC571906	
15x16in	PTC571970	
16x16in	PTC571911	
16x17in	PTC585140	
16x18in	PTC571913	
17x17in	PTC585149	
18x17in	PTC585153	
18x18	PTC585155	
20x18in	PTC571925	
22x18in	PTC571950	

JAY Lite w/Positioning Inserts


14x14in	PTC571904	
14x16in	PTC571905	
15x15in	PTC571907	
15x16in	PTC571909	

Pressure Care Cushions


Description	Code	Image
-------------	------	-------

16x15in	PTC571908	
16x16in	PTC571910	
16x18in	PTC571912	
17x17in	PTC585150	
17x18in	PTC585151	
18x16in	PTC585152	
18x17in	PTC585154	
18x20in	PTC585160	
20x20in	PTC571930	


JAY Zip Cushions

10x10in	PTC585184	
12x12in	PTC585193	
12x16in	PTC585185	
14x14in	PTC585192	
14x16in	PTC585190	
15x15in	PTC585191	

Action Gel Pad Cushion

16x16in	PTC547200	
18x16in	PTC547100	
18x18in	PTC547000	
20x20in	PTC546990	
Cushion - Cover for Action Gel (*Specify Size*)	PTC547210	

JAY J2 Cushions

14x16in	PTC563400	
16x16in	PTC562500	
16inx18in	PTC562700	
16x20in	PTC563000	



Pressure Care Cushions

Description	Code	Image
-------------	------	-------

17x17in	PTC563800	
18x16in	PTC561100	
18x18in	PTC561700	
18x18in	PTC561300	
18x20in	PTC562100	
20x16in	PTC563700	
20x18in	PTC564100	
20x20in	PTC561900	
24x20in	PTC562400	

JAY J 2 Deep Contour Cushions

14x14in	PTC565300	
14x16in	PTC565400	
15x15in	PTC565700	
15x17in	PTC565800	
16x16in	PTC566000	
16x18in	PTC566100	
16x20in	PTC565900	
17x15in	PTC566300	
17x17in	PTC566400	
18x16in	PTC566600	
18x18in	PTC566700	
18x20in	PTC566500	
20x16in	PTC565150	
20x18in	PTC565000	
20x20in	PTC564900	
22x18in	PTC565200	

JAY J2 Deep Contour Cushions w/ Incontinence Cover

14x16in	PTC565600	
16x18in	PTC566200	
20x16in	PTC565156	
20x18in	PTC565100	

Ottobock Advantage Cushions

16x16in	PTC550050	
16x18in	PTC550075	
18x16in	PTC550100	
18x18in	PTC550150	

Diffuser Cushions

440x420mm	GSS340100	
490x460mm	GSS340104	
440x420mm	GSS340103	
490x460mm	GSS340105	

Multi-Purpose Cushions

Multipurpose Cushion	GSS340120	
Multipurpose Cushion w/ Waterproof Cover	GSS340121	

Pressure Care Cushions

Description	Code	Image
-------------	------	-------

JAY J2 Recline Cushions

15.5x18in	PTC562800	
15.5x20in	PTC563100	
15.5x20in	PTC563200	
16x16in	PTC562550	
17x17in	PTC563900	
18x18in	PTC561400	
18x18in	PTC561500	
18x18in	PTC561600	
18x20in	PTC562200	
20x18in	PTC564200	
20x20in	PTC562000	

JAY J2 Plus Bariatric Cushions

20x18in	PTC567000	
21x19in	PTC567050	
22x18in	PTC567100	
22x20in	PTC567200	
22x22in	PTC567250	
24x18in	PTC567300	

ROHO Harmony Cushions

14x16in	PTC554400	
16x16in	PTC554500	
16x18in	PTC554600	
17x16in	PTC554700	
17x18in	PTC554800	
18x16in	PTC554900	
18x18in	PTC555000	
20x18in	PTC555100	

ROHO Mini-Max Cushions

15x15in	PTC585800	
16x16in	PTC585900	
18x16in	PTC586000	

ROHO Mosaic Cushions

16x16in	PTC586100	
18x16in	PTC586200	
18x18in	PTC586300	

ROHO Quadtro Select Cushions

LP 11x11in (6x6cell)	PTC590300	
LP 11x13in (6x7cell)	PTC590500	
LP 11x15in(6x8cell)	PTC590600	
LP 13x13in (7x7cell)	PTC590800	
LP 13x15in (7x8cell)	PTC591100	
LP 15x15in (8x8cell)	PTC591900	
LP 16.75x15in (9x8cell)	PTC592900	
LP 16.75x20in (9x11cell)	PTC592500	
LP 18.5x15in (10x8cell)	PTC588400	



Pressure Care Cushions

Description	Code	Image
LP 15x13in (8x7cell)	PTC591700	
LP 15x16.75in (8x9cell)	PTC592100	
LP 15x18.5in (8x10cell)	PTC591400	
LP 16.75x16.75in (9x9cell)	PTC593200	
LP 16.75x18.5in (9x10cell)	PTC592300	
LP 18.5x16.75in (10x9cell)	PTC588700	
LP 18.5x18.5in (10x10cell)	PTC587900	
LP 18.5x20in (10x11cell)	PTC588100	
LP 20x16.75in (11x9cell)	PTC589500	
LP 20x18.5in (11x10cell)	PTC588900	
LP 20x20in (11x11cell)	PTC589100	
LP 22x18.5in (12x10cell)	PTC589700	
MP 16.75x16.75in (9x9cell)	PTC593205	
MP 16.75x18.5in (9x10cell)	PTC593220	
MP 18x18in (10x10cell)	PTC593250	
HP 11x13in (6x7cell)	PTC590400	
HP 12x12in (6x6cell)	PTC590200	
HP 13x13in (7x7cell)	PTC590700	
HP 13x15in (7x8cell)	PTC591000	
HP 15x13in (8x7cell)	PTC591600	
HP 15x15in (8x8cell)	PTC591800	
HP 15x16.75in (8x9cell)	PTC592000	
HP 15x17.75in (7x9cell)	PTC591200	
HP 15x18.5in (8x10cell)	PTC591300	
HP 15x20in (8x11cell)	PTC591500	
HP 16.75x15in (9x8cell)	PTC592800	
HP 16.75x16.75in (9x9cell)	PTC593000	
HP 16.75x18.5in (9x10cell)	PTC592200	
HP 16.75x20in (9x11cell)	PTC592400	
HP 17x15in (9x7cell)	PTC592700	
HP 18.25x22in (10x12cell)	PTC588200	
HP 18.5x15in (10x8cell)	PTC588300	
HP 18.5x16.75in (10x9cell)	PTC588500	
HP 18.5x18.5in (10x10cell)	PTC587700	
HP 18.5x20in (10x11cell)	PTC588000	
HP 20x15in (11x8cell)	PTC589300	
HP 20x16.75in (11x9cell)	PTC589400	
HP 20x18.5in (11x10cell)	PTC588800	
HP 20x20in (11x11cell)	PTC589000	
HP 20x22in (11x12cell)	PTC589200	
HP 22x18.5in (12x10cell)	PTC589600	
HP 22x20in (12x11cell)	PTC589800	
HP 22x22in (12x12cell)	PTC589900	
HP 24x18.5in (13x10cell)	PTC590000	
HP 24x20in (13x11cell)	PTC590100	

Pressure Care Cushions

Description	Code	Image
MP 16.75x16.75in (9x9cell)	PTC548600	
MP 16.75x18.5in (9x10cell)	PTC548601	
MP 16.75x20.25in (9x11cell)	PTC548602	
MP 18.5x16.75in (10x9cell)	PTC548603	
MP 18.5x18.5in (10x10cell)	PTC548604	
MP 18.5x20in (10x11cell)	PTC548605	
MP 20x18.5in (11x10cell)	PTC548607	
MP 20.25x20.25in (11x11cell)	PTC548608	
HP 11.5x10in (6x5cell)	PTC544900	
HP 11.5x8.5in (6x4cell)	PTC544800	
HP 13.25x13.25in (7x7cell)	PTC545100	
HP 13.25x16.75in (7x9cell)	PTC545200	
HP 15x11.5in (8x6cell)	PTC545375	
HP 15x13.25in (8x7cell)	PTC545400	
HP 15x15in (8x8cell)	PTC545500	
HP 15x16.75in (8x9cell)	PTC545700	
HP 15x18.5in (8x10cell)	PTC545301	
HP 16.75x13in (9x7cell)	PTC546450	
HP 16.75x15in (9x8cell)	PTC546500	
HP 16.75x16.75in (9x9cell)	PTC546600	
HP 16.75x18.5in (9x10cell)	PTC545900	
HP 16.75x20.25in (9x11cell)	PTC546100	
HP 18.5x16.75in (10x9cell)	PTC544000	
HP 18.5x18.5in (10x10cell)	PTC543700	
HP 18.5x20in (10x11cell)	PTC543900	
HP 20.25x16.75in (11x9cell)	PTC544400	
HP 20.25x20.25in (11x11cell)	PTC544300	
HP 20x18.5in (11x10cell)	PTC544200	
HP 22x18.75in (12x10cell)	PTC544500	
LP 11.5x13.25in (6x7cell)	PTC545000	
LP 15x11in (8x6cell)	PTC545377	
LP 15x15in (8x8cell)	PTC545600	
LP 15x16.75in (8x9cell)	PTC545800	
LP 15x18.5in (8x10cell)	PTC545300	
LP 16.75x16.75in (9x9cell)	PTC546700	
LP 16.75x18.5in (9x10cell)	PTC546000	
LP 16.75x20in (9x11cell)	PTC546200	
LP 18.5x16.75in (10x9cell)	PTC544100	
LP 18.5x18.5in (10x10cell)	PTC543800	
LP 18.5x20.25in (10x11cell)	PTC543910	

ROHO Single Compartment Cushions



Pressure Care Cushions

Description	Code	Image
LP 20.25x20.25in (11x11cell)	PTC544350	
LP 20x18.5in (11x10cell)	PTC544250	
LP 22x18.75in (12x10cell)	PTC544550	
LP 22x22in (12x12cell)	PTC544600	

ROHO Enhancer Cushions

15x15in (8x8cell)	PTC553300	
15x16.75in (8x9cell)	PTC553400	
16.75x16.75in (9x9cell)	PTC553600	
16.75x18.5in (9x10cell)	PTC553500	
18x16.75in (10x9cell)	PTC553200	
18.5x18.5in (10x10cell)	PTC553100	

ROHO Contour Select Cushions

15x15in (8x8cell)	PTC594100	
15x16.75in (8x9cell)	PTC594200	
16.75x15in (9x8cell)	PTC594600	
16.75x16.75in (9x9cell)	PTC594700	
16.75x18.5in (9x10cell)	PTC594300	
16.75x20in (9x11cell)	PTC594400	
18.5x16.75in (10x9cell)	PTC593800	
18.5x18.5in (10x10cell)	PTC593500	
18.5x20in (10x11cell)	PTC593700	
20x18.5in (11x10cell)	PTC593900	
20x20in (11x11cell)	PTC594000	

JAY Easy Cushions

14x14in Flat	PTC580002	
14x16in Flat Bottom	PTC580005	
15x15in Flat Bottom	PTC580100	
16x16in Flat Bottom	PTC580300	
17x17in Flat Bottom	PTC580350	
20x20in Flat Bottom	PTC580660	
14x16in Curved Bottom	PTC580004	
15x15in Curved Bottom	PTC580050	
16x16in Curved Bottom	PTC580200	
16x18in Curved Bottom	PTC580400	
18x16in Curved Bottom	PTC580500	
18x18in Curved Bottom	PTC580600	
19x16in Curved Bottom	PTC580640	

JAY Fusion Cushions

16x16in - Air Dual	PTC571201	
17x17in - Air Dual	PTC571203	
18x18in - Air Dual	PTC571202	
16x16in - w/fluid insert	PTC571204	
18x16in - w/fluid insert	PTC571208	
18x18in - w/fluid insert	PTC571210	

Pressure Care Cushions

Description	Code	Image
ROHO Nexus Spirit Cushions		
14x14in	PTC586500	
14x16in	PTC586600	
15x15in	PTC586700	
16x16in	PTC586800	
16x18in	PTC586900	
17x17in	PTC587000	
18x16in	PTC587100	
18x18in	PTC587200	
20x18in	PTC587300	

ROHO Airlite Cushions

16x18in	PTC551300	
18x18in	PTC551400	
20x18in	PTC551450	

ROHO Over Toilet Aid Cushions

Over Toilet Aid Toilet Seat - Pressure Care - ROHO Cushion	BTC067800	
--	-----------	--

ROHO Hybrid Elite Cushions

1R 15x15in	PTC594899	
1R 15x16in	PTC594900	
1R 16x16in	PTC594910	
1R 16x17in	PTC594911	
1R 16x18in	PTC594920	
1R 17x17in	PTC594930	
1R 17x18in	PTC594931	
1R 18x16in	PTC594935	
1R 18x18in	PTC594940	
1R 18x20in	PTC594945	
1R 20x18in	PTC594950	
1R 20x20in	PTC594951	
1R 22x20in	PTC594947	
2R 16x16in	PTC594952	
2R 16x18in	PTC594953	
2R 17x18in	PTC594971	
2R 18x18in	PTC594954	
2R 20x18in	PTC594955	

ROHO Dual Compartment Cushions

HP 18.5x18.5in (10x10cell)	PTC547300	
HP 18.5x16.75in	PTC547500	
HP 15x13.25in (8x7cell)	PTC547550	
HP 15x15in	PTC547600	
HP 15x16.75in	PTC547800	
HP 16.75x18.5in	PTC548000	
HP 16.75x15in	PTC548200	



Pressure Care Cushions

Description	Code	Image
HP 16.75x16.75in (9x9cell)	PTC548400	
LP 18.5x18.5in (10x10cell)	PTC547400	
LP 15x15in	PTC547700	
LP 15x16.75in	PTC547900	
LP 16.75x18.5in	PTC548100	
LP 16.75x15in (9x8cell)	PTC548300	
LP 16.75x16.75in (9x9cell)	PTC548500	

JAY Cushion Covers

18x18in - Stretch	PTC575700	
18x18in - Air Exch	PTC575800	
18x18in - Air Exch	PTC573600	
18x18in - Incontinence	PTC575900	
Deep Contour 18x18in	PTC579000	
Deep Contour 18x18in - Air Exch	PTC579100	
Deep Contour 18x18in - Incontinence	PTC579200	

ROHO Cushion Covers

HP - 18x18in	PTC555300	
HP - 18.5x18.5in	PTC555500	
LP - 18.5x18.5in	PTC556000	
HP - 18.5x20in	PTC555650	
HP - 18.5x16.75in	PTC555700	
LP - 20x18.5in	PTC555740	
LP - 18.5x16.75in	PTC555750	
HP - 18.5x15in	PTC555760	
HP - 20x18.5in	PTC555800	
HP - 20x20in	PTC555900	
HP - 20x18in	PTC556000	
HP - 13.5x13.5in	PTC556100	
HP - 15x15in	PTC556200	
LP - 15x15in	PTC556300	
HP - 15x16.75in	PTC556400	
LP - 15x16.75in	PTC556500	
HP - 16.75x18.5in	PTC556600	
LP - 16.75x18.5in	PTC556700	
HP - 16.75x20in	PTC556705	
LP - 16.75x20in	PTC556710	
LP - 16.75x16.75in	PTC556750	
HP - 16.75x16.75in	PTC556800	

ROHO Nexus Spirit Cushion Cover

ROHO Nexus Spirit 18x18in	PTC587400	
---------------------------	-----------	--

Pressure Care Cushions

Description	Code	Image
ROHO Universal Cushion Cover		
LP 18.5x18.5in	PTC596300	
HP 18.5x18.5in	PTC596200	
HP 18.5x22in	PTC596500	
HP 22x18.5in	PTC596750	
HP 22x20in	PTC596800	

Equagel General Cushions

41x41cm (16x16in)	PTC551043	
41x46cm (16x18in)	PTC551045	
41x51cm (16x20in)	PTC551047	
46x41cm (18x16in)	PTC551049	
46x46cm (18x18in)	PTC551051	
46x51cm (18x20in)	PTC551053	
51x41cm (20x16in)	PTC551055	
51x46cm (20x18in)	PTC551057	
51x51cm (20x20in)	PTC551059	
56x41cm (22x16in)	PTC551061	
56x46cm (22x18in)	PTC551063	
56x51cm (22x20in)	PTC551065	

Equagel Adjustable Cushions

41x41cm (16x16in)	PTC551011	
41x46cm (16x18in)	PTC551013	
46x41cm (18x16in)	PTC551015	
46x46cm (18x18in)	PTC551017	
51x41cm (20x16in)	PTC551019	
51x46cm (20x18in)	PTC551021	
56x41cm (22x16in)	PTC551023	
56x46cm (22x18in)	PTC551025	

Equagel Protector Cushions

41x41cm (16x16in)	PTC551093	
41x46cm (16x18in)	PTC551095	
46x46cm (18x18in)	PTC551097	
51x41cm (20x16in)	PTC551099	
51x46cm (20x18in)	PTC551101	
56x41cm (22x16in)	PTC551103	
56x46cm (22x18in)	PTC551105	
61x46cm (24x18in)	PTC551109	

Cushion Accessories

Cushion Acc - ROHO Adaptor Pad	PTA507400	
Cushion Acc - ROHO - Pump	PTA509400	

Dressing Aids



Dressing Aids

Description	Code	Image
-------------	------	-------

Notch Cord Sock Aid

Notch - Cord Style	DLD269100	
--------------------	-----------	---

Stocking Aid - Cord

Cord Style	DLD269275	
------------	-----------	---

Stocking Aid - Dorking




Dorking (Homecraft)	DLD269260	
Dorking - Side Slot- (091095462)	DLD269270	

EzyAs Compression Stocking Applicator


Compression Stocking/ Garment Aid - Ezy As - Large Blue	DLD269003	
---	-----------	---

Dressing Aids

Description	Code	Image
-------------	------	-------

Compression Stocking/ Garment Aid - Ezy As - Medium Yellow	DLD269002	
Compression Stocking/ Garment Aid - Ezy As - Small Red	DLD269001	
Compression Stocking/ Garment Aid - Handle Attachment	DLD269004	

Compression Stocking/Garment Aid

Compression Stocking Aid - Eureka	DLD269010	
--------------------------------------	-----------	---

Sock Aid

Terry Cloth Sock Aid	DLD269090	
----------------------	-----------	---



Dressing Aids

Description	Code	Image
-------------	------	-------

Easy-Slide Applicator - Donning Aid

O/Toe Stockings - Med	GMS340952	
O/Toe Stockings - Lge	GMS340953	
O/Toe Stockings - XL	GMS340954	
Magnide for C/Toe Stockings - Lge	GMS340956	

JOBST Stocking Donner

Stocking Aid - JOBST Stocking Donner	PTS641090	
--------------------------------------	-----------	--

Button Hook

Button Hook - Kings	DLD269300	
---------------------	-----------	--

Shoe Horn

Rigid Plastic	DLD269160	
---------------	-----------	--

Dressing Aids

Description	Code	Image
-------------	------	-------

Metal	DLD269170	
-------	-----------	--

Dressing Stick

Dressing Stick	DLD269310	
----------------	-----------	--

Shoelace - Spring Coilers

Black 600mm	DLD269140	
White 600mm	DLD269150	
Royal Blue 600mm	DLD269153	

Leg Lifter

Leg Lifter	DLD269390	
------------	-----------	--



Product brochures also available! Please contact us for a copy






Kitchen Aids






Kitchen Aids

Description	Code	Image
-------------	------	-------

Meal Trolley/Walker

Meal Trolley - Aspire - Height Adj - Plastic Trays	DLG272060	
Meal Tray Walker - Aspire - Hand Brakes - Height Adj - Plastic Trays - 125kg	DLG272061	
Aspire Walker Spares - Plastic Trays - Meal Tray Walker/Trolley	WAA686926	

Counter Top




Counter Top - Spread Board Plastic - Reflex	DLK275730	
Counter Top - Bread Board with Spikes - Reflex	DLK275710	
Counter Top - Food Preparation System - Reflex	DLK275720	

Kitchen Aids



Description	Code	Image
-------------	------	-------

Counter Top - Kitchen Workstation - Reflex	DLK275700	
--	-----------	---

Peeler

Peeler - Good Grips - Y Peeler	DLK275100	
Peeler - Good Grips - Swivel Peeler	DLK275110	
Peeler - Gordon	DLK275630	

Kettle Tipper

Kettle Tipper - Cordless	DLK276610	
Kettle Tipper	DLK276600	



Gripping & Turning Aids

Gripping & Turning Aids

Description	Code	Image
-------------	------	-------

Jar Opener

Jar Opener - Twister	DLK275520	
Jar Opener - Strongboy	DLK275550	
Jar & Bottle Opener - Undo-It	DLK275530	
Jar Opener - Good Grips **	DLK275120	
Jar Opener - One-Touch Automatic (Batteries incl)	DLG272340	

Bottle Holder

Bottle Holder - Belliclamp	DLK275580	
----------------------------	-----------	--

Gripping & Turning Aids

Description	Code	Image
-------------	------	-------

Strap Handle

Strap-handle - Multi-purpose	DLK276430	
------------------------------	-----------	--

Can Opener

Can Opener - Good Grips	DLK275090	
Bottle Opener - One-Touch Automatic (Batteries incl)	DLG272310	
Can Opener - One-Touch Automatic II (Batteries incl)	DLG272330	
Opener - Can Pull	DLK275610	



Gripping & Turning Aids

Description	Code	Image
Rectangular 25x18cm - Blue	DLG271390	
Rectangular 35x25cm - Blue	DLG271400	
Round 14cm - Blue - Small	DLG271370	
Round 19cm - Blue - Large	DLG271380	

Key Turner

One - White	DLG271310	
Two - Blue	DLG271320	

Tap Turner

Tap Turner - Red & Blue	DLG271280	
-------------------------	-----------	--

Gripping & Turning Aids

Description	Code	Image
Tap Turner - Contour	DLG271300	

Uno Grips - Red & Blue

Clear	DLG271100	
Semi Clear	DLG271110	
Red/Blue	DLG271120	

Utensil Handles - Foam Tubing

White 6mmx1m (19mm)	DLG271210	
White 6mmx1m (31mm)	DLG271220	
White 9mmx1m (31mm)	DLG271240	
White 12mmx1m (31mm)	DLG271260	
Grey 6mm x 1m (31mm)	DLG271230	
Grey 9mmx1m (31mm)	DLG271250	
Grey 12mmx1m (31mm)	DLG271270	



Aidacare have a range of Product Application Guides available for the following areas:





Cutlery

Cutlery

Description	Code	Image
-------------	------	-------

Supergrip Cutlery

Tablespoon	DLK275320	
Teaspoon	DLK275310	

Supergrip Bendable Cutlery

Teaspoon	DLK275360	
Tablespoon	DLK275370	
Fork	DLK275380	

Reflex Cutlery

Chefs Knife 20.3cm	DLK275670	
Carving Knife 20.3cm	DLK275660	
Preparation Knife 12.7cm	DLK275650	
Slicing Knife 25cm	DLK275680	

Nelson Knife

Nelson Knife - Queens	DLK275810	
-----------------------	-----------	--

Queens Cutlery

Fork	DLK275770	
Knife	DLK275760	
Cutlery Set	DLK275870	

Cutlery

Description	Code	Image
-------------	------	-------

Queens Angled Cutlery

Angled Built Up Spoon Left	DLK275850	
Angled Built Up Spoon Right	DLK275860	
Angled Built Up Fork Left	DLK275830	
Angled Built Up Fork Right	DLK275840	
Angled Built Up Knife	DLK275820	

Kings Softcoated Head Spoon

Softcoated Spoon	DLK275880	
Softcoated Teaspoon	DLK275890	

Kings Junior Spoon

Junior Spoon	DLK275930	
--------------	-----------	--

Kings Modular Cutlery

Modular Fork	DLK275910	
Modular Spoon	DLK275920	
Soup Spoon	DLK275950	
Teaspoon	DLK275940	
Nelson Knife	DLK275970	
Rocker Knife Large	DLK276040	
Rocker Knife Small	DLK275980	



Cutlery

Description	Code	Image
-------------	------	-------

Kings Angled Cutlery

Spoon Left Hand	DLK276020	
Spoon Right Hand	DLK276030	
Fork Left Hand	DLK276000	
Fork Right Hand	DLK276010	

Kings Handles

Standard Handle	DLK276070	
Heavy Handle	DLK276100	
Slim Handle	DLK276120	

Lightweight Cutlery

Cutlery Set with Foam Handle (Set)	DLK276220	
Teaspoon w/ Foam Handle	DLK276210	
Spoon w/ Foam Handle	DLK276200	
Knife w/ Foam Handle	DLK276180	
Fork w/ Foam Handle	DLK276190	

Kings Oval Built Up Sleeves

Handles - Oval	DLK276140	
Handles - Large	DLK276160	
Handles - Small	DLK276150	

Cutlery

Description	Code	Image
-------------	------	-------

Lightweight Angled Cutlery

Spoon w/ Foam Handle Angled Left	DLK276260	
Spoon w/ Foam Handle Angled Right	DLK276270	
Fork w/ Foam Handle Angled Left	DLK276240	
Fork w/ Foam Handle Angled Right	DLK276250	
Knife w/ Foam Handle Angled	DLK276230	

Kings Lightweight Handles

Lightweight Foam Handle - Large	DLK276080	
Lightweight Foam Handle - Slim	DLK276090	

Kings Mushroom Handle

Cutlery - Kings - Handles - Mushroom Handle	DLK276130	
---	-----------	--

Caring Cutlery

Spoon	DLK276320	
Fork	DLK276310	
Knife	DLK276300	
Spoon - Angled Right	DLK276330	



Introducing our Aspire Lifter & Sling Product Video



To view, visit the link below

<https://www.aidacare.com.au/knowledge/videos/>



Crockery

Crockery

Description	Code	Image
Manoy Plates		
Large	DLK276460	
Small	DLK276450	
Inner Lip Plate		
Plates & Dishes - Inner Lip Plate	DLK275290	
Economy Scoop Dish		
Plates & Dishes - Scoop Dish Economy	DLK275300	
Plate Guard		
Plates & Dishes - Plate Guard Large Clear (A7062)	DLK275450	
Incurve Plate Surround		
Plates & Dishes - Incurve - Plate Surround **	DLK276490	
Caring Mug		
Cups & Mugs - Caring - Mug Clear	DLK276510	

Crockery

Description	Code	Image
Two Handled Mug		
Cups & Mugs - Two Handled Mug	DLK276520	
Anti Spill Mug		
Cups & Mugs - Anti Spill Mug	DLK276500	
Nosey Cup		
115ml Blue Tint	DLK275480	
225ml Blue Tint	DLK275490	
340ml Blue Tint	DLK275500	
225ml White	DLK275510	
Feeding Cup		
8mm	DLK276560	
4mm	DLK276550	
Feeding Cup - Adjustable Spout		
4mm	DLK276570	
8mm	DLK276580	
Dysphagia Cup		
Almond	DLK275460	
Green	DLK275470	

Reaching Aids



Reaching Aids

Description	Code	Image
-------------	------	-------

Pickup Reacher

Short 340mm	DLR279220	
Standard 600mm	DLR279150	
Long 750mm	DLR279170	
Long w/ Accessories 750mm	DLR279180	
Extra Long 900mm	DLR279230	

Handireacher

Standard 610mm	DLR279110	
Long 760mm	DLR279120	
Extra Long 890mm	DLR279130	


Reaching Aids

Description	Code	Image
-------------	------	-------

Economy Reacher

660mm	DLR279010	
838mm	DLR279020	

Easireach 2 Reacher

Folding - Short	DLR279080	
380mm	DLR279040	
530mm	DLR279050	
660mm	DLR279060	
Folding - Long	DLR279090	
810mm	DLR279070	



General Household & Dementia Aids

General Household & Dementia Aids

Description	Code	Image
-------------	------	-------



Vehicle Handy Bar

Vehicle - Handy Bar **	DLG271670	
------------------------	-----------	---

Easigrip Scissors

Chiropracist	DLG271360	
Pointed 45mm	DLG271340	
Round 45mm	DLG271330	

Long Handled Comb/Brush

Long Handled Comb	DLD269340	
Hairbrush - Long Handled Brush	DLD269350	
Cleaning - Long Handle Toilet Brush	LAC500565	

Long Handled Dustpan & Brush

Dust Pan & Broom - Long Handled	DLG271500	
---------------------------------	-----------	---

General Household & Dementia Aids

Description	Code	Image
-------------	------	-------

Wool Duster

Cleaning - Telescopic Wool Duster	LAC500505	
-----------------------------------	-----------	---


Powerfibre Duster

Cleaning - Powerfibre Duster	LAC500520	
------------------------------	-----------	---


Enviro Dust Wand

Cleaning - Oates Enviro Dust Wand	LAC500510	
-----------------------------------	-----------	---


Window Washer

Cleaning - Swivel Head	LAC500570	
------------------------	-----------	---

Carpet Sweeper

Cleaning - Carpet Sweeper	LAC500585	
---------------------------	-----------	---

Enviro Dusting Floor Mop

Cleaning - Enviro Dusting Floor Mop, 350mm, Oates	LAC500555	
---	-----------	---



General Household & Dementia Aids

Description	Code	Image
-------------	------	-------

Magic Mop

Cleaning - Vileda - Magic Mop	LAC500531	
-------------------------------	-----------	---

Microfibre Flat Mop

Cleaning - Flat Mop: Microfibre Flat Mop (600mm)	LAC500545	
Cleaning - Flat Mop Refill: Microfibre Flat Mop Pad (600mm)	LAC500550	

Super Swish Spray Mop

Cleaning - Sabco Super Swish Spray Mop	LAC500539	
Cleaning - Sabco Super Swish Mop - Refills	LAC500541	

General Household & Dementia Aids

Description	Code	Image
-------------	------	-------

Eager Beaver Floor Mop

Floor Tool Kit (inc3 Pads) (Plus Handle)	LAC500577	
Floor Pad - Red - 10 Pack	LAC500578	
Floor Pad - Blue - 10 Pack	LAC500579	

Clothes Line - Handyline

Standard Large - 16.1m line	DLG272255L	
Standard Small - 11.2m line	DLG272255S	



Aidacare have a range of Product Application Guides available for the following areas:





Lifters & Slings

Aidacare A150F Folding Aluminium Patient Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390600	Lifter A150F - Folding Lifter - Aluminium	1873 mm	1280 mm	660 mm	32.9 kg



Aidacare A205 Aluminium Compact Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390610	Lifter A205 - Patient Lifter - Aluminium	2045 mm	1135 mm	620 mm	38.1 kg



Aidacare 200 Steel Patient Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390620	Lifter 200 - Patient Lifter - Steel	1940 mm	1280 mm	665 mm	43.9 kg





Aidacare A320 Aluminium Patient Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390630	Lifter A320 - Bariatric Lifter - Aluminium	2165 mm	1505 mm	765 mm	63 kg



Aidacare A200S Aluminium Standing Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390640	Lifter A200S - Patient Stander - Aluminium	1455 mm	985 mm	590 mm	30.4 kg



Aidacare 200S Steel Standing Lifter

Code	Description	Overall Height	Overall Length	Overall Width	Lifter Weight
LSS390650	Lifter 200S - Patient Stander - Steel	1790 mm	1150 mm	655 mm	44.9 kg



Lifters & Slings

Description	Code	Image
-------------	------	-------

Aidacare Manual Pivot Frame

Manual Pivot Frame - Aidacare Aspire Lifters - 227kg	LSS390670	
--	-----------	---

Invacare Birdie Lifter

Lifter - Patient Lifter - Invacare Birdie Compact Lifter - 150kg SWL	LSS364100	
Lifter - Patient Lifter - Invacare Birdie Lifter - 170kg SWL - Manual leg	LSS364111	

Lifters & Slings

Description	Code	Image
-------------	------	-------

Linak Battery

2.9 Amp/Hour without Port	MBT401190	
2.9 Amp/Hour with Port	MBT401200	

Linak Battery Charger

Battery Charger Plug Pack - Linak - with Port	MBT401220	
---	-----------	---



Lifters & Slings			Lifters & Slings		
Description	Code	Image	Description	Code	Image
Aspire Deluxe General Purpose Sling - Polyester			Aspire Deluxe Padded Standing Sling - Polyester		
Small - 200kg	LSS390700		Small - 200kg	LSS390720	
Medium - 200kg	LSS390705		Medium - 200kg	LSS390725	
Large - 300kg	LSS390710		Large - 300kg	LSS390730	
XLarge - 300kg	LSS390715		XLarge - 300kg	LSS390735	
Aspire Deluxe General Purpose Sling - Mesh			Kerry General Purpose Sling		
Small - 200kg	LSS390716		Lifter Sling - Gen Purp Polyester - Kerry - Medium	LSS381000	
Medium - 200kg	LSS390717				
Large - 300kg	LSS390718				
XLarge - 300kg	LSS390719				
Aspire Deluxe General Purpose Sling w/ Head Support - Polyester					
Small - 200kg	LSS390740				
Medium - 200kg	LSS390745				
Large - 300kg	LSS390750				
XLarge - 300kg	LSS390755				
Aspire Deluxe General Purpose Sling w/ Head Support - Mesh					
Small - 200kg	LSS390756				
Medium - 200kg	LSS390757				
Large - 300kg	LSS390758				
XLarge - 300kg	LSS390759				
Aspire Deluxe General Purpose Sling w/ Head Support & Clips - Polyester					
Small - 200kg	LSS390900				
Medium - 200kg	LSS390905				
Large - 300kg	LSS390910				
XLarge - 300kg	LSS390915				
Aspire Deluxe Hygiene/Access Sling - Polyester					
Small - 200kg	LSS390760				
Medium - 200kg	LSS390765				
Large - 300kg	LSS390770				
XLarge - 300kg	LSS390775				
Aspire Deluxe Hygiene/Access Sling w/Head Support - Polyester					
Small - 200kg	LSS390780				
Medium - 200kg	LSS390785				
Large - 300kg	LSS390790				
XLarge - 300kg	LSS390795				
Aspire Deluxe Hygiene/Access Sling w/Head Support & Clips - Polyester					
Small - 200kg	LSS390940				
Medium - 200kg	LSS390945				
Large - 300kg	LSS390950				
XLarge - 300kg	LSS390955				

Product brochures also available!
Please contact us for a copy

aidacare
aspire

PATIENT LIFTERS & SLINGS

*Mobile Patient Lifters & Slings
Designed for a Range of Environments,
from Community to Facility*

LIFTERS

- Aidacare Aspire A100F
- Aidacare Aspire A205
- Aidacare Aspire 200
- Aidacare Aspire A320

STANDSSES

- Aidacare Aspire A200S
- Aidacare Aspire 200S
- Aspire Swift Patient Mover

OPTIONS

- Aidacare Aspire Inline Weigh Scale
- Aidacare Aspire Manual & Powered Post Frames

SLINGS

- Aspire Deluxe General Purpose
- Aspire Deluxe General Purpose with Head Support
- Aspire Deluxe Hygiene Access
- Aspire Deluxe Hygiene Access with Head Support
- Aspire Deluxe Hammock with Head Support
- Aspire Deluxe Low Friction
- Aspire Deluxe Standing

www.aidacarelifters.com.au
1300 133 120

Transfer Boards & Sheets



Transfer Boards & Sheets

Description	Code	Image
-------------	------	-------

Transfer Board - The Slider

610x230mm No Cut Out	TSE672810	
540x360mm - w/cut out	TSE672800	
645x297mm	TSE672790	
645x360mm	TSE672780	

Transfer Board - Curved

Transfer - Board - Curved	TSE672920	
Transfer - Board - Curved Q2	TSE672192	

Transfer Board - Sallyboard

Large - 180x58cm	TSE672196	
Medium - 150x58cm	TSE672195	

Transfer Boards & Sheets

Description	Code	Image
-------------	------	-------

Curved	TSE672191	
--------	-----------	---

Slide Board

Slide Board - Wheelchair Cut Out	TSE672460	
----------------------------------	-----------	---

Aspire Transfer Slide Sheet

Small - 1000x1450mm	TSE673070	
Large - 2000x1450mm	TSE673080	

Patterson Transfer Board

Patterson Transfer Board - Long	TSE672910	
---------------------------------	-----------	--



Transfer & Walk Belts


Transfer & Walk Belts

Description	Code	Image
-------------	------	-------


Walk Belt - Tedi Grip

Small	TSE673200	
Medium	TSE673190	
Large	TSE673180	

Transfer Belt - SallyStrap

130cm - Small	TSE672193	
170cm - Large	TSE672194	

Transfer Walk Belt w/ Buckle

630-1140mm - Small	TSE672040	
1010-1570mm - Large	TSE672030	



Product brochures also available!
Please contact us for a copy




Other Transfer Aids



Other Transfer Aids

Description	Code	Image
-------------	------	-------


Aspire Swivel Cushion

Transfer - Swivel Cushion - Aspire	TSE673056	
------------------------------------	-----------	---

Aspire Slide Sheet

Transfer - One Way Slide - Non Slip - Aspire	TSE673050	
Transfer - One Way Slide - Foot Operated - Non Slip - Aspire	TSE673052	
Transfer - Positioning Slide - Aspire	TSE673054	

Roller Sheet

Transfer - Sheet - Roller Sheet	TSE672510	
---------------------------------	-----------	---



Introducing our Aspire Lifter & Sling Product Video




To view, visit the link below

<https://www.aidacare.com.au/knowledge/videos/>


Mobility Scooters




Mobility Scooters

Description	Code	Image
Scooter Cane Holder - (Requires L Piece Bracket SCA668310)		
Single Rear	SCA668230	
Single Side	SCA668250	
Double Rear	SCA668220	
Double Side	SCA668240	

Scooter L Piece

Scooter Accessory - L Piece for Accessories	SCA668310	
---	-----------	---



Scooter Walking Frame Holder

Scooter Walking Frame Holder (Requires L-Piece bracket SCA668310 & Anti-Rattle Clamp SCA668312)	SCA668390	
---	-----------	---


Scooter Safety Flag

Scooter Safety Flag	SCA668370	
---------------------	-----------	---


Mobility Scooters

Description	Code	Image
Scooter Seat Sheepskin Cover		
Scooter Seat Sheepskin Cover	SCA668410	
Scooter Wind Shield (Requires SCA668380 Sun Canopy)		
Side Wind Shield	SCA668395	
Rear Wind Shield	SCA668396	
Shield Front Section Including Frame	SCA668397	

Scooter Sun Canopy

Scooter Sun Canopy Complete (Includes all Mounting Hardware)	SCA668380	
--	-----------	---

Scooter Rain Cover


Scooter Rain Cover (Requires SCA668380 Sun Canopy)	SCA668392	
--	-----------	---




Mobility Scooters

Description	Code	Image
-------------	------	-------


Scooter Mirror

Scooter Mirror - Celebrity & Pathrider Scooters	SCA668300	
---	-----------	---


Scooter Rear Bag

Scooter Rear Bag (Includes all Mounting Hardware)	SCA668340	
---	-----------	---


Scooter Rear Basket

Scooter Rear Basket	SCA668180	
---------------------	-----------	--


Scooter Lockable Rear Pod

Scooter Lockable Rear Pod	SCA668399	
---------------------------	-----------	---


Scooter Oxygen Bottle Holder

Scooter Oxygen Bottle Holder	SCA668330	
------------------------------	-----------	---

Scooter Cup Holder


Scooter Cup Holder	SCA668170	
--------------------	-----------	---

Scooter Anti Rattle Clamp

Scooter Anti Rattle Clamp for inLin Piece	SCA668311	
---	-----------	---

Mobility Scooters

Description	Code	Image
-------------	------	-------

Scooter Acc - Anti Rattle Clamp for Walking Frame Holder	SCA668312	
--	-----------	---

Dust Cover

Dust Cover	SCA668210	
------------	-----------	---



Product brochures also available!
Please contact us for a copy



WHEELCHAIR RANGE

aspire for...

Comfort
Posture
Mobility











www.aspirecare.com.au
1300 133 120





Modular & Rubber Ramps

Modular & Rubber Ramps

Description	Code	Image
-------------	------	-------


Rubber Ramps - Threshold

Small - 1010x200x25mm	RAM656211	
Large - 1140x400x50mm	RAM656212	

Folding Ramps - Silver Star

610mm	RAM656370	
910mm	RAM656380	
1220mm	RAM656390	
1520mm	RAM656400	
1830mm	RAM656410	


Decpac Senior Ramp

1350mm	RAM656320	
1650mm	RAM656330	
2000mm	RAM656340	
2500mm	RAM656345	
3500mm	RAM656350	
900mm	RAM656360	
Multi Purpose **	RAM656180	
Multi Purpose Ramp Carry Bag only	RAM656190	

Modular & Rubber Ramps

Description	Code	Image
-------------	------	-------

Decpac EBL Ramp

900mm	RAM656170	
1200mm	RAM656090	
1350mm	RAM656100	
1650mm	RAM656110	
2000mm	RAM656120	
2500mm	RAM656130	
3000mm	RAM656140	
3500mm	RAM656150	




Introducing our Aspire Lifter & Sling Product Video



<https://www.aidacare.com.au/knowledge/videos/>

Seat Walkers/ Rollators




Aspire Vogue Seat Walker

Code	Description	Overall Height	Overall Length	Overall Width	Total Weight
WAF705300RE	Vogue Walker - Red	695 - 920 mm	680 mm	600 mm	6.7 kg
WAF705300WH	Vogue Walker - White	695 - 920 mm	680 mm	600 mm	6.7 kg



Red White




Aspire Vogue Carbon Fibre Seat Walker

Code	Description	Overall Height	Overall Length	Overall Width	Total Weight
WAF705350	Vogue Walker - Carbon Fibre	695 - 920 mm	680 mm	600 mm	5 kg



Carbon Fibre




Aspire Vogue Indoor Walker

Code	Description	Overall Height	Overall Length	Overall Width	Total Weight
WAF705100RE	Vogue Indoor Walker - Red / Black	860 - 970 mm	650 mm	540 mm	5.2 kg
WAF705100CH	Vogue Indoor Walker - Champagne / Black	860 - 970 mm	650 mm	540 mm	5.2 kg



Red / Black Red / Champagne





Seat Walkers/Rollators

Description	Code	Image
Red - HBrakes	WAF705700RE	
Blue - HBrakes	WAF705700BL	

Aspire Classic 6 Walker - 6" Wheels - 130kg

Blue - HBrakes	WAF709600BL	
Red - HBrakes	WAF709600RE	

Aspire Mini Walker - 6" Wheels - 130kg

Blue - HBrakes	WAF705800BL	
Red - HBrakes	WAF705800RE	

Aspire Deluxe Walker - 8" Wheels - 130kg

Blue - Inline HBrakes	WAF705750BL	
Red - Inline HBrakes	WAF705750RE	

Aspire Push Down Walker - 7" Wheels - 130kg

Blue	WAF705600BL	
Red	WAF705600RE	

Aspire XL Walker

Seat Walker - Aspire XL - Silver - HBrakes - 8in Wheels - 180kg	WAF750020SL	
---	-------------	--

Knee Walker

Knee Walker - for Lower Leg Injuries	WAF712000	
--------------------------------------	-----------	--

Lets Go Indoor Walker

Indoor Walker - Lets Go Indoor - Wood Grain w/ Tray & Bag	WAF706150	
---	-----------	--

Lets Go Outdoor Walker

Seat Walker - Let's Go - Outdoor Walker - Beige (with backrest & bag)	WAF706160	
---	-----------	--

Seat Walkers/Rollators

Description	Code	Image
Oxygen Bottle Holder Bag - Aspire VOGUE Walkers	WAA695600	
Crutch & Cane Holder - Aspire VOGUE Walkers	WAA695610	
Universal LED Night Lamp - Aspire VOGUE Walkers	WAA695630	
Single Hand Brake System - Aspire VOGUE Walkers	WAA695640	
IV Pole - Height Adjustable - Aspire VOGUE Walkers	WAA695650	

Aspire Walker Accessories

Oxygen Bottle Holder - Aspire Walkers - Classic/Mini/Deluxe/Pdown/XL	WAA686905	
Clip on Carry Bag - Classic/Mini/Deluxe/PDown	WAA686915	
Clip on Carry Bag - XL Walker	WAA686920	
Walking Stick Holder - Aspire Walkers	WAA686925	
Underseat Basket - Classic/Mini/Deluxe/Pdown	WAA686910	
Underseat Basket - XL	WAA686911	
Vogue Forearm Attachment Kit - Aspire Vogue Walkers	WAA695660	




Product brochures also available!





Walking Frames

Description	Code	Image
Alum Folding - Small	WAF706700K	
Alum Folding - Medium	WAF706800K	
Alum Folding - Large	WAF706900K	
Alum Non Folding - Small	WAF704900K	
Alum Non Folding - Medium	WAF704800K	
Alum Non Folding - Large	WAF704700K	


Walking Frame - 4 Wheels - Side Folding

Adult - 100kg SWL	WAF706920	
Youth - 100kg SWL	WAF706930	

Walking Frame - 2 Front Wheels

Walking Frame w/ 2 Front Wheels - Folding Dual Release - Heavy Duty - 227kg	WAF704750	
---	-----------	---

Walking Frames

Description	Code	Image
Small	WAF706700	
Medium	WAF706800	
Large	WAF706900	

Walking Frame - Folding

Silver	WAF705000	
--------	-----------	---

Walking Frame - Aluminium Non Folding

Small	WAF704900	
Medium	WAF704800	
Large	WAF704700	

Walking Frame - Folding Dual Release

Walking Frame - Folding Dual Release - Heavy Duty - 227kg	WAF704745	
---	-----------	---

Aspire Vogue Forearm Seat Walker

Code	Description	Overall Height	Overall Length	Overall Width	Total Weight
WAF705400	Vogue Forearm Walker - White	1020 - 1270 mm	660 mm	620 mm	9 kg



Forearm Walkers

Description	Code	Image
-------------	------	-------

Forearm Walker - Walking Tutor - Folding - 100kg

Adult w/ Padded Rest - 850-1350mm	WAF709100	
Adult w/ Rear Glides - 850-1350mm	WAF709000	

Forearm Walker - Gutter Frame - Folding

Forearm Walker - Gutter Frame - 2 Wheeled - Folding	WAF705950	
---	-----------	---


Forearm Walkers

Description	Code	Image
-------------	------	-------

Forearm Walker - Days - Gutter Frame

Forearm Walker - DAYS - Gutter Walker w/ Handbrakes - 150kg	WAF705960	
---	-----------	---

Forearm Walker - Gutter Frame HD

Forearm Walker - Gutter Frame Heavy Duty w/ Hand Brakes - 180kg	WAF699200	
---	-----------	---




Walking Sticks

Walking Sticks

Description	Code	Image
-------------	------	-------


Walking Stick - T Handle

Walking Stick - Ht Adj - T Handle - Alum **	WAS743600	
---	-----------	---

Walking Stick - Height Adjustable - Swan Neck

w/ Strap - Alum	WAS743700	
T Handle - Alum	WAS743800	
Walking Stick - Ht Adj - T Handle Swan Neck - Heavy Duty	WAS744750	

Walking Stick - Crook Handle

Walking Stick - Ht Adj - Crook Handle (foam) - Alum	WAS744200	
---	-----------	---

Walking Stick - Height Adjustable

Palm Grip Right - Alum	WAS744300	
Palm Grip Left - Alum	WAS744400	

Walking Sticks

Description	Code	Image
-------------	------	-------


Walking Stick - T Handle - Folding

Walking Stick - Ht Adj - T Handle - Folding - Alum	WAS744700	
--	-----------	---

Walking Stick - Quad - Height Adjustable

Large Base	WAS744500	
Small Base	WAS744600	

Walking Stick Holder

Walk Stick Acc - Walking Stick Holder - Table/Chair attachable	WAA686880	
--	-----------	---



Product brochures also available!








Crutches







Crutches			Crutches		
Description	Code	Image	Description	Code	Image
Underarm Crutches - Adjustable Handle - Aluminium			Forearm Crutches - Forearm - Aluminium - Child		
Child/Paediatric - 800-950mm	WAC692300		4-7 Years Small - 406-559mm	WAC691700	
Small - 940-1145mm	WAC692500		6-10 Years - 457-656mm	WAC691600	
Medium - 1120-1320mm	WAC692200				
Tall Large - 1320-1520mm	WAC692400				
Underarm Crutches - Adjustable Handle - Bariatric			Forearm Crutches - Forearm - Aluminium - Youth		
249kg - 1550-1750mm	WAC688300		Small - 530-735mm	WAC692000	
249kg - 1750-1980mm	WAC688200				
Aspire Underarm Crutches - Wing Nut - Aluminium			Underarm Crutches - Forearm - Aluminium - Adult		
136kg - Youth Small	WAC692503		Medium - 635-910mm	WAC691400	
136kg - Medium	WAC692502		Tall / Large - 820-1075mm	WAC691800	
136kg - Adult Tall	WAC692501				
Underarm Crutches - Gutter			Underarm Crutches - Forearm - Ergo Grip - Aluminium		
Crutches Gutter - Alum - 980-1235mm	WAC694200		Medium - 580-835mm	WAC691500	
			Large - 680-935mm	WAC691900	



Crutches			Crutches		
Description	Code	Image	Description	Code	Image
Crutch Cushion			Crutch Protector		
Padding Grey Axilla	WAA686240		Gel & Foam	PTA509020	
Padding Grey Axilla Bariatric	WAA686250		Gel	PTA509021	
Padding Axilla - Wing Nut Crutches	WAA686245				
Underarm Crutch Pad - Sheepskin			Crutch Hand Grip w/ Padding		
Axilla	PTA497500		Hand Grip	WAA686220	
Hand Grip	PTA497510		Hand Grip - Grey - Wing Nut Crutches	WAA686225	
			Hand Grip - Grey - Bariatric Crutches	WAA686230	



General Accessories

General Accessories			Crutches		
Description	Code	Image	Description	Code	Image
Rubber Ferrule			Glides		
Quad Cane - 13mm Black	WAA686280		White 19mm Black Dot	WAA686690	
Quad Cane - 16mm Black	WAA686290		White 22.2mm Green Dot	WAA686700	
Type Z 16mm Grey	WAA686300		White 25mm Red Dot	WAA686710	
Type Z 18mm Grey	WAA686310		White 25.4mm Yellow Dot	WAA686720	
Type Z 22mm Grey	WAA686320		White 28.6mm Blue Dot	WAA686730	
Type Z 25mm Grey	WAA686330				
Type D 18mm Black	WAA686332				
Aspire Walker Meal Tray			Wheels		
Aspire - Meal Tray Walker/Trolley	WAA686926		Walker Accessory - Walker Front Castors (PAIR)	WAA686450	

Manual Wheelchairs



Aspire Evoke 2 & 2 HD Wheelchair

Features & Benefits:

- Seat & Backrest Height Adjustable
- Customisable Backrest
- Seat Depth Adjustable
- Height and angle adjustable front castors (8")
- Folding Lightweight Aluminium Design
- Height Adjustable back cranes 17" - 19"
- Optional Contoured Padded Overlay
- Swing-back Armrests
- Swing Away and removable legrests
- Height and angle adjustable footrests
- Quick release 24" solid rear wheels



Code	MWS449500	MWS449510	MWS449520	MWS449530	MWS449540 (HD)	MWS449550 (HD)
Seat Width	350 mm	400 mm	450 mm	500 mm	570 mm	610 mm
Seat Depth	375 - 525 mm	375 - 525 mm	375 - 525 mm	375 - 525 mm	375 - 525 mm	375 - 525 mm
Seat Height	440 - 480 mm	440 - 480 mm	440 - 480 mm	440 - 480 mm	440 - 480 mm	440 - 480 mm
Overall Width	560 mm	610 mm	660 mm	710 mm	780 mm	820 mm
Backrest Height	430 - 490 mm	430 - 490 mm	430 - 490 mm	430 - 490 mm	430 - 490 mm	430 - 490 mm
Armrest Height	205 - 225 mm	205 - 225 mm	205 - 225 mm	205 - 225 mm	205 - 225 mm	205 - 225 mm
SWL	140 kg	140 kg	140 kg	140 kg	180 kg	180 kg

Contoured Padded Overlays

Code	WCA848330	WCA848340	WCA848350	WCA848360	WCA848370	WCA848380
Overall Width	350 mm	400 mm	450 mm	500 mm	570 mm	610 mm



Aspire Transit 2 Manual Folding Wheelchair

Code	MWS449410	MWS449420	MWS449430
Seat Width	400 mm	450 mm	500 mm
Seat Depth	400 mm	400 mm	400 mm
Seat Height	460 - 500 mm	460 - 500 mm	460 - 500 mm
Overall Width	540 mm	590 mm	640 mm
Backrest Height	440 mm	440 mm	440 mm
Armrest Height	210 mm	210 mm	210 mm
SWL	140 kg	140 kg	140 kg



Features & Benefits:

- Folding Steel Frame
- Padded Nylon Upholstery
- Optional Contoured Padded Overlay
- Quick Release Solid 12" Rear Wheels & Height Adjustable, 8" Castors
- Dacron Fire Retardant Fabric
- Swing-Away, Removable Armrests
- Attendant Push-to-lock Brakes
- Removable Swing Away Legrests & Height Adjustable Footrest



Aspire Assist 2 Manual Folding Wheelchair

Code	MWS449300	MWS449310	MWS449320	MWS449330
Seat Width	350 mm	400 mm	450 mm	500 mm
Seat Depth	400 mm	400 mm	400 mm	400 mm
Seat Height	460 - 500 mm	460 - 500 mm	460 - 500 mm	460 - 500 mm
Overall Width	570 mm	620 mm	670 mm	720 mm
Backrest Height	440 mm	440 mm	440 mm	440 mm
Armrest Height	210 mm	210 mm	210 mm	210 mm
SWL	140 kg	140 kg	140 kg	140 kg



Features & Benefits:

- Foot Steel Frame
- Optional Contoured Padded Overlay
- Quick Release Solid 24" Rear Wheels
- Height Adjustable, 8" Castors
- Dacron Fire Retardant Fabric
- Swing-back Removable Armrests
- Push-to-lock Brakes
- Removable Swing Away Legrests (Elevating / Articulating optional accessory)
- Height Adjustable Footrest



Aspire Rehab RS Tilt in Space Wheelchair

Code	MWS449600	MWS449610	MWS449620
Seat Width	410 - 460 mm	410 - 460 mm	410 - 460 mm
Seat Depth	440 - 530 mm	440 - 530 mm	440 - 560 mm
Seat Height (Front)	410 - 440 mm	410 - 440 mm	410 - 440 mm
Seat Height (Rear)	400 - 475 mm	400 - 475 mm	400 - 475 mm
Overall Width	690 mm	740 mm	790 mm
Backrest Height	460 - 710 mm	460 - 710 mm	460 - 710 mm
Armrest Height	310 - 460 mm	310 - 460 mm	310 - 460 mm
SWL	180 kg	180 kg	180 kg



Features & Benefits:

- Attendant Operated Drum Brakes & Push-to-lock Handbrakes
- Removable, Height & Depth Adjustable Armrest
- Tilt-in-Space Functionality
- Quick release, Height & Depth Adjustable Solid 24" Rear Wheels & Height Adjustable 6" CastorsAdjustable, 8" Castors
- Height, Depth, Angle & Width Adjustable Removable Headrest
- Depth Adjustable Seat, Height & Depth Adjustable Backrest
- Swing-away, Removable, Padded & Elevating Legrest
- Height & Angle Adjustable Footrest



Aspire Lite Transit Manual Wheelchair

Code	MWS442301
Seat Width	460 mm
Seat Depth	400 mm
Seat Height	460 mm
Overall Width	600 mm
Backrest Height	465 mm
Armrest Height	235 mm
SWL	120 kg



Features & Benefits:

- Ultra-lightweight wheelchair comes with attendant brakes fitted as standard
- Folding frame for easy storage and transport
- Weighs just 11.5kg
- Ergonomic handles – provides comfort for the caregiver
- Push button, height adjustable arm rests
- Swing away, flip up foot plates – this facilitates user transfers
- 12" rear wheels – with easy to access brakes, for improved riding over outdoor terrain



Aspire Rehab RX JNR Wheelchair

Code	MWS449630
Seat Width	320 - 370 mm
Seat Depth	440 - 490 mm
Seat Height (Front)	410 - 440 mm
Seat Height (Rear)	400 - 475 mm
Overall Width	600 mm
Backrest Height	460 - 710 mm
Armrest Height	310 - 460 mm
SWL	180 kg



Features & Benefits:

- Attendant Operated Drum Brakes & Push-to-lock Handbrakes
- Removable, Height & Depth Adjustable Armrest
- Quick release, Height & Depth Adjustable Solid 20" Rear Wheels & Height Adjustable 6" Castors
- Tilt-in-Space Functionality
- Height, Depth, Angle & Width Adjustable Removable Headrest
- Depth Adjustable Seat, Height & Depth Adjustable Backrest
- Swing-away, Removable, Paddeed & Elevating Legrest
- Height & Angle Adjustable Footrest



Aspire Transit 2 Manual Folding Wheelchair

Code	MWS442300
Seat Width	460 mm
Seat Depth	400 mm
Seat Height	460 mm
Overall Width	550 mm
Backrest Height	465 mm
Armrest Height	235 mm
SWL	120 kg



Features & Benefits:

- An ultra-light, folding frame – for easy storage and transport
- Weighs just 9.5kg
- Ergonomic handles – provides comfort for the caregiver
- Push button, height adjustable arm rests
- Swing away, flip up foot plates – this facilitates user transfers.
- Rear wheel locks – for stability and security.



ottobock.

Ottobock Start Manual Wheelchair - Start M2S



Code	MWC449214 - Standard	MWC449215 - Standard	MWC449216 - Standard
Description	16in (405mm) w/ seat depth	18in (455mm) w/ seat depth	20in (480mm) w/ seat depth
Seat Width	405mm	455mm	505mm
SWL	130kg	130kg	130kg

Code	MWC449204 - Adjustable	MWC449220 - Adjustable	MWC449210 - Adjustable
Description	15in (380mm) w/ seat depth	16in (405mm) w/ seat depth	18in (455mm) w/ seat depth
Seat Width	405mm	455mm	505mm
SWL	130kg	130kg	130kg

Features & Benefits:

- Patented, tool-free, quick seat depth adjustment)
- Lightweight 14.7kg total, 9.6 transport weight.
- Antibacterial push handles
- Excellent component compatibility
- Adjustable economy back upholstery
- Aluminium frame
- Height of side panel can be adjusted (without tools) and adjustable depth set (with tools)



Breezy Basix 2 Fixed Manual Wheelchair



Code	MWS449091	MWS449097	MWS449108	MWS449114
Description	14in - Solid Wheels - RED	20in - Solid Wheels - RED	16in - Solid Wheels - RED	18in - Solid Wheels - RED
Seat Depth	410 - 460 mm	410 - 460 mm	410 - 460 mm	410 - 460 mm
Seat Width	520 mm	520 mm	400 mm	460 mm
Seat Height	425 - 510 mm	425 - 510 mm	425 - 510 mm	425 - 510 mm
Backrest Height	410 - 460 mm	410 - 460 mm	410 - 460 mm	410 - 460 mm
SWL	125 kg	125 kg	125 kg	125 kg

Code	MWS449106	MWS449113	MWS449120
Description	16in - Pneumatic Wheels - SILVER	18in - Pneumatic Wheels - SILVER	20in - Pneumatic Wheels - SILVER
Seat Depth	410 - 460 mm	410 - 460 mm	410 - 460 mm
Seat Width	400 mm	460 mm	520 mm
Seat Height	425 - 510 mm	425 - 510 mm	425 - 510 mm
Backrest Height	410 - 460 mm	410 - 460 mm	410 - 460 mm
SWL	125 kg	125 kg	125 kg

Features & Benefits:

- Lightweight chair
- Available in self-propelled or transit
- Adjustable push handle height
- Tension adjustable backrest upholstery
- Swing away legrests as standard
- Drum break option available (16, 18 & 20")
- 41 - 46 cm adjustable seat depth
- Castor angle adjustment 0, 3 & 6 degrees



Breezy Basix 2 Folding Manual Wheelchair



Code	MWS449104	MWS449112	MWS449119
Description	Paediatric - 14in	18in - Att Prop.	20in - Att Prop.
Seat Depth	410 - 460 mm	410 - 460 mm	410 - 460 mm
Seat Width	360 mm	460 mm	520 mm
Seat Height	425 - 510 mm	425 - 510 mm	425 - 510 mm
Backrest Height	410 - 460 mm	410 - 460 mm	410 - 460 mm
SWL	125 kg	125 kg	125 kg

Code	MWS449111	MWS449118	MWS449123
Description	16in - Transit	18in - Transit	20in - Transit
Seat Depth	410 - 460 mm	410 - 460 mm	410 - 460 mm
Seat Width	400 mm	460 mm	520 mm
Seat Height	425 - 510 mm	425 - 510 mm	425 - 510 mm
Backrest Height	410 - 460 mm	410 - 460 mm	410 - 460 mm
SWL	125 kg	125 kg	125 kg



Features & Benefits:

- Folding lightweight chair
- Available in self-propelled or transit
- Adjustable push handle height (+/- 3")
- Available with folding backrest
- Tension adjustable upholstery on backrest
- Swing away legrests as standard
- Drum break option available (16, 18 & 20")
- 41 - 46 cm adjustable seat depth
- Castor angle adjustment 0, 3 & 6 degrees



Breezy Rubix 2 XL Heavy Duty Wheelchair



Code	MWS449130	MWS449135	MWS449140
Seat Depth	410 - 510 mm	410 - 510 mm	410 - 510 mm
Seat Width	520 mm	520 mm	520 mm
Seat Height	375 - 510 mm	375 - 510 mm	375 - 510 mm
Backrest Height	410 - 460 mm	410 - 460 mm	410 - 460 mm
SWL	170kg	170kg	170kg



Features & Benefits:

- This lightweight wheelchair designed for individual user requirements
- Adjustable back upholstery, adjustable seat depth and angle, angle adjustable footrests and comfort seating
- Castor fork with increased adjustment (0, 3 & 6 degrees) and added stability
- Seat width up to 600mm, with a 170kg user weight
- Axle plate for seat height and centre of gravity adjustment
- Elevating footrest made of aluminium cast housing for greater stability
- Wider, rounded front armrest for greater comfort and safety



Heavy Duty & Super Heavy Duty Folding Wheelchair

Code	MWS442910 - HD	MWS442911 - HD	MWS442912 - HD	MWS442901 - Super HD	MWS442902 - Super HD
Seat Depth	400 mm	450 mm	450 mm	500 mm	500 mm
Seat Width	500 mm	550 mm	600 mm	650 mm	700 mm
Seat Height	500 mm	500 mm	500 mm	500 mm	500 mm
Backrest Height	400 mm	400 mm	400 mm	500 mm	500 mm
SWL	136kg	204kg	204kg	318kg	318kg



Features & Benefits:

- Reinforced wheelchair for bariatric users
- Extra wide seat widths and reinforced construction accommodates larger users
- Moulded composite mag wheels
- Swing-away detachable footrests and removable padded armrests



Action 3 Lever Drive Folding Manual Wheelchair

Code	MWS442600	MWS442700	MWS442800	MWS442900
Description	Left Hand Drive	Right Hand Drive	Left Hand Drive	Right Hand Drive
Seat Depth	410mm	410mm	410mm	410mm
Seat Width	440mm	440mm	500mm	500mm
Seat Height	400 - 500mm	400 - 500mm	400 - 500mm	400 - 500mm
Backrest Height	410mm	410mm	410mm	410mm
SWL	125kg	125kg	125kg	125kg



Features & Benefits:

- Single arm propelled wheelchair for homes or hospitals
- Push pull lever designed for users with adequate upper limb strength in one arm
- Folding frame for easy storage and transportation
- Swing-away detachable footrests and armrests
- 24" quick-release pneumatic rear wheels plus 8" front castors
- Durable steel construction, with a 14.2kg weight



Action 3NG Folding Light Manual Wheelchair

Code	MWS449205	MWS449210	MWS449215
Description	16in - Self Prop.	18in - Self Prop.	20in - Self Prop.
Seat Depth	405, 455, 510 mm	405, 455, 510 mm	405, 455, 510 mm
Seat Width	405, 455, 510 mm	405, 455, 510 mm	405, 455, 510 mm
Total Width	405 mm	455 mm	505 mm
Backrest Height	400, 430, 510 mm	400, 430, 510 mm	400, 430, 510 mm
SWL	125 kg	125 kg	125 kg



Features & Benefits:

- Adaptable and modular to accommodate progressive changes in user needs
- Low transfer weight, easy to fold and to dismantle, making it easily transportable in a small car
- 20 rear wheel positions.
- Brake and Brake extension, fold the brakehandle for easier side transfers
- Magnesium castor attachment
- Simple to refurbish



Action Folding Manual Wheelchair - Paediatric

Code	MWS449180	MWS449181
Seat Width	12 - 15 in	8 - 11 in
Seat Depth	11 - 13 in	10 - 12 in
SWL	80 kg	80 kg

Code	MWS449182	MWS449183	MWS449192
Seat Width	8 - 11 in	8 - 11 in	12 - 15 in
Seat Depth	12 - 14 in	8 - 10 in	13 - 15 in
SWL	80 kg	80 kg	80 kg

Features & Benefits:

- 4 seat width adjustment positions - 205 - 380mm
- 20 different rear wheel mounting positions for improved performance height options
- 2 tray options - Slide-on tray table or swing-away tray
- A front stabiliser kit is available to increase forward stability for active children
- 4 different types of spoke guards are offered to match the different wheelchair colours



ottobock.

Discovery T-Max - Tilt in Space Wheelchair

Code	MWC427105
Seat Depth	300 - 400 mm
Seat Width	300 - 400 mm
Seat Tilt	-5° - +50°
Backrest Height	530 / 610 mm
Backrest Tilt	85° - 120°
SWL	90kg

Features & Benefits:

- Suitable for use indoors and outdoors
- Approved for transport in a car (wheelchair accessible vehicle)
- Air-filled or puncture-proof front and rear wheels available
- Width, depth and height of seat adjustable in 2,5 cm / 1 inch increments
- Cable-free, maintenance-friendly seat angle release with foot pedal and release lock
- Smooth seat angle adjustment without pneumatic spring thanks to adjustable centre of gravity



Standard Power Wheelchairs



Pride Jazzy Select Elite 2 Power Wheelchair



Code	PWS654255
Total Length	873 mm
Total Width	584 mm
Max Speed	6.5 km/h
SWL	90kg

Features & Benefits:

- The Jazzy Select Elite features in-line, front-wheel drive technology, which gives excellent stability and maneuverability.
- The Jazzy Select Elite comes standard with Comfort high-back seat with headrest
- Reliable blend of power, performance and style
- Excellent for indoor and outdoor performance
- Available in Jazzy Red or Jazzy Blue




Backrests


Backrests

Description	Code	Image
-------------	------	-------

JAY Wheelchair Basic Backrests

20x18in	PTB532505	
18x18in	PTB532510	
14x18in	PTB532503	
16x20in	PTB532504	

JAY J2 Wheelchair Backrests

12in	PTB532260	
13in	PTB532265	
14in	PTB532270	
14in Tall	PTB532280	
15in Tall	PTB532300	
16in	PTB532310	
16in Tall	PTB532320	
17in	PTB532330	
17in Tall	PTB532340	
18in	PTB532350	
18in Tall	PTB532360	
19in	PTB532370	
19in Tall	PTB532380	
20in	PTB532390	
20in Tall	PTB532400	

Backrests

Description	Code	Image
-------------	------	-------


JAY J2 Deep Contour Backrests

14in	PTB532180	
15in	PTB532190	
16in	PTB532200	
17in	PTB532210	
18in	PTB532220	
19in	PTB532230	
20in	PTB532240	


JAY - Lumbar Support Combi

Wheelchair Back - JAY - Lumbar Support Combi - Adult	PTB532250	
--	-----------	---


Synergy Shape Backrest

16in	PTB533180	
18in	PTB533190	
20in	PTB533200	

Tru Comfort Backrest

22x22in	PTB533202	
Plus - 18x22in	PTB533205	
Plus - 18x25in	PTB533207	

Curved Seat Backrest

17in w x 17in h	PTB533155	
18in w x 19in h	PTB533160	
21in w x 21in h	PTB533170	

Accessories



Accessories

Description	Code	Image
-------------	------	-------

Wheelchair Wheels - 24inch for Aspire ASSIST/ EVOKE/RELAX

Pneumatic	MWT469865	
Solid	MWT469840	

Wheelchair Wheels - 12inch for Aspire TANSIT

Pneumatic	MWT469851	
Solid	MWT469850	

Rex Tyres

24x1 3/8 - Primo - Orion Grey	MWT450230	
24x1 3/8 - Solid	MWT450221	
24x1 - Solid	MWT450222	
200x50 - Grey Rib Tread	MWT450030	
12.5x2 1/4 - Grey Stud Tread	MWT450170	
300x8 - Grey Stud Tread - Solid	MWT450195	
12 1/2x2 1/4 - Straight Valve	MWT450710	
24x1 3/8 - Straight Valve	MWT450750	
200x50 - 90x90 Degree	MWT450650	
4.10/3.5x4	MWT450670	
300x8	MWT450700	

QHDC Tube

Tube - QHDC - 24x1 - Straight Valve	MWT450760	
-------------------------------------	-----------	--

Accessories

Description	Code	Image
-------------	------	-------

Wheelchair Lap Belt

Wheelchair Lap Belt - Metal Seatbelt Buckle - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX	WCA848300	
W'Chair Belt - Generic	WCA858560	
W'Chair Belt w/ Buckle - Attachment	WCA858590	
W'Chair Belt w/ Anchor Strap	WCA858600	
WChair Belt - 50in Lap Belt	WCA858200	
WChair Belt - 60in Lap Belt	WCA858205	

Wheelchair Chest Strap - Elastic

Small - 75mmx1520mm	WCA860539	
Lge - 115mmx1520mm	WCA860540	
XLge - 115mmx1520mm	WCA860541	

Wheelchair Shoulder Harness - Rear Pull

Flat Mount - Sml	WCA860468	
Flat Mount - Med	WCA860470	
Flat Mount - Lge	WCA860472	

Wheelchair Trunk Support - Monoflex

XSml - 76mmx318mm	WCA860546	
Small - 89mmx397mm	WCA860547	
Med - 102mmx476mm	WCA860548	
Lge - 114mmx556mm	WCA860549	
Extra Lge - 127mmx635mm	WCA860550	



Accessories

Description	Code	Image
Wheelchair Leg Harness - Padded		
Sml - 230mm	WCA860419	
Med - 330mm	WCA860420	
Lge - 380mm	WCA860421	

Wheelchair Calf Strap - Padded

Small - 32in Long	WCA860464	
Med - 36in Long	WCA860465	
Lge - 40in Long	WCA860466	

Wheelchair Ankle Huggers - Velcro Attach

X-Sml	WCA860422	
Sml	WCA860423	
Med	WCA860424	
Lge	WCA860425	
XL	WCA860426	

Wheelchair Foot Positioners

Small	WCA859280	
Medium	WCA859290	
Large	WCA859300	
Extra Large	WCA859310	

Wheelchair Tray Table

Aspire Rehab RS/Rehab RX	WCA848560	
16in for Breezy Basix	WCA858020	
18in for Breezy Basix	WCA858021	
20in for Breezy Basix	WCA858022	

Wheelchair Channel Forearm Pad - Ottobock

Large	WCA858980	
Medium	WCA858990	
Small	WCA859000	
One Piece	WCA859010	

Wheelchair Height Adjustable Armrest - Aspire Evoke 2

Left	WCA848390	
Right	WCA848400	

Wheelchair Gutter Armrest Pad - Aspire Evoke 2

Left	WCA848410	
Right	WCA848420	

Wheelchair Elevating Legrests Other

Alu f/plates - Heavy Duty WChair (20-24in)	WCA858005	
Alu f/plates - Super Heavy Duty WChair (26-30in)	WCA858006	
Breezy Basix2	WCA858016	

Accessories

Description	Code	Image
Wheelchair Elevating Legrests - Standard		
Left - 320/350/400/410mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848100	
Left - 450/460/500/510mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848120	
Left - 570/610mm - Aspire Assist2/Transit2/Evoke2	WCA848140	
Right - 320/350/400/410mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848110	
Right - 450/460/500/510mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848130	
Right - 570/610mm - Aspire Assist2/Transit2/Evoke2	WCA848150	

Wheelchair Elevating Legrests - Advanced

Left - 320/350/400/410mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848160	
Left - 450/460/500/510mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848180	
Left - 570/610mm - Aspire Assist2/Transit2/Evoke2	WCA848200	
Right - 320/350/400/410mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848170	
Right - 450/460/500/510mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848190	
Right - 570/610mm - Aspire Assist2/Transit2/Evoke2	WCA848210	

Wheelchair Stump Support

Left - Silver - 320/350/400/450mm - Aspire Assist2/Transit2/ Evoke2	WCA848220	
Left - Silver - 500/550/570/610mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848240	
Right - Silver - 320/350/400/450mm - Aspire Assist2/Transit2/ Evoke2	WCA848230	
Right - Silver - 500/550/570/610mm - Aspire Assist2/Transit2/ Evoke2/Rehab RS/Rehab RX	WCA848250	



Accessories

Description	Code	Image
-------------	------	-------


Wheelchair Stump Support Other

Left - Heavy Duty (20-24in)	WCA858007	
Right - Heavy Duty (20-24in)	WCA858008	
Left - Super Heavy Duty (26-30in)	WCA858009	
Right - Super Heavy Duty (26-30in)	WCA858010	
Left for Breezy Basix Wheelchair	WCA858014	
Right for Breezy Basix Wheelchair	WCA858015	
Extension Bracket - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX	WCA848260	

Wheelchair Pommel - Flip Down

2in Pad	WCA858340	
4in Pad	WCA858330	

Wheelchair Foot Plate Protectors

WChair Foot Plate Protectors	WCA859275	
------------------------------	-----------	---


Wheelchair Gel Foot Plate Protector Covers

Footplate - Gel Footplate Protector Covers 20.32 x 15.24cm	PTA508600	
--	-----------	---

Wheelchair Gel Armrest Pads

3.5x12in	PTA508400	
3.5x16in	PTA508450	

Wheelchair Armrest Cover

WChair Armrest Cover - Sheepskin - 335mmx18-22 cm circumf.	WCA858195	
--	-----------	---

Wheelchair Gel Calf Support Panel

Calf - Gel Calf Support Panel	PTA508500	
-------------------------------	-----------	---

Wheelchair Knee Protector Buttons

Knee Protector - Gel Knee Protector Buttons 3.18cm - Swing Away Legrests	PTA508800	
--	-----------	---

Accessories

Description	Code	Image
-------------	------	-------

Wheelchair Legrest - Hanger Protector Gel

Legrest - Hanger Protectors Gel	PTA509000	
---------------------------------	-----------	---

Wheelchair Abductor Button Gel Padded

Abductor Button Gel Padded	PTA514100	
----------------------------	-----------	---


Wheelchair Crutch Holder

Crutch Holder for Breezy Basix Wheelchair	WCA858012	
---	-----------	---


Wheelchair Brake Extension Lever

Left - Aspire Assist2/Transit2/Evoke2	WCA848310	
Right - Aspire Assist2/Transit2/Evoke2	WCA848320	


Wheelchair Deluxe Headrest

Wheelchair Headrest Assembly - Deluxe Style - Aspire Evoke2	WCA848480	
---	-----------	---


Wheelchair IV Pole

Wheelchair IV Pole - Aspire Evoke2/Transit2/Assist2/Rehab RS/Rehab RX	WCA848280	
---	-----------	---


Wheelchair Anti Tip Bar

WChair Anti Tip Bar - Aspire Evoke2/Assist2/Transit2/Rehab RS/Rehab RX	WCA848270	
--	-----------	---

Wheelchair Oxygen Bottle Holder

Wheelchair Oxygen Bottle Holder - Aspire Assist2/Transit2/Evoke2	WCA848290	
--	-----------	---

Wheelchair Oxygen Bottle Holder Bag

Oxygen Bottle Holder Bag - Aspire Vogue Walkers & Rehab Wheelchairs	WAA695600	
---	-----------	---



Accessories

Description	Code	Image
-------------	------	-------

Wheelchair Camber Kit

Wheelchair Camber Kit - Aspire Evoke2	WCA848490	
---------------------------------------	-----------	--

Wheelchair Tie Downs

Wheelchair Tie Downs (set of 4) - Aspire Evoke2/Rehab RS/Rehab RX	GSP493290	
---	-----------	--

Wheelchair Drum Brake Kit

Aspire Evoke2	WCA848500	
16in for Aspire Rehab RX	WCA848510	
24in for Aspire Rehab RX	WCA848520	

Wheelchair Therapist Kit for Hand Adjustment

Wheelchair Therapist Kit for Hand Adjustment - Aspire Rehab RS/Rehab RX	WCA848550	
---	-----------	--

Wheelchair Extended Back Cane

Wheelchair Extended Back Cane - Pair - 19/21in - Aspire Evoke2	WCA848470	
--	-----------	--

Wheelchair Lateral Chest Support Assembly

Left - Aspire Evoke2	WCA848430	
Right - Aspire Evoke2	WCA848440	

Wheelchair Back Cushion - 4 Way Stretch PU

410mm - Aspire Rehab RS	GSP492620	
460mm - Aspire Rehab RS	GSP492630	
510mm - Aspire Rehab RS	GSP492640	

Wheelchair Back Cushion - Woven Breathable Mesh

410mm - Aspire Rehab RX	GSP492660	
460mm - Aspire Rehab RX	GSP492670	
510mm - Aspire Rehab RX	GSP492680	

Wheelchair Seat Cushion - 4 Way Stretch PU

410mm - Aspire Rehab RS	GSP492480	
460mm - Aspire Rehab RS	GSP492490	
510mm - Aspire Rehab RS	GSP492500	

Accessories

Description	Code	Image
-------------	------	-------

Wheelchair Seat Cushion - Woven Breathable Mesh

410mm - Aspire Rehab RX	GSP492520	
460mm - Aspire Rehab RX	GSP492530	
510mm - Aspire Rehab RX	GSP492540	

Wheelchair Depth Extension Kit

50mm - Aspire Rehab RS/Rehab RX	WCA848530	
100mm - Aspire Rehab RS/Rehab RX	WCA848540	

Smart Drive

Smart Drive - inc batt & mounting hardware (As per Script)	WCA858508	
--	-----------	--

Power Assist

Power Drive Unit - Pride	WCA858501	
Power Stroll Std P/Pack with Reverse	WCA858505	
Stroll HD - (As per Script)	WCA858511	

iXpress Power Add On

iXpress - Sunrise - Power Add On - For Manual Wheelchairs (As per script)	WCA858506	
---	-----------	--

Regal Battery - Gel

12V 33Ah	MBT403014	
----------	-----------	--

Pride Battery Charger

3 Amp/Hour	MBT401575	
------------	-----------	--

Exercise & Resistance Training



Exercise & Resistance Training

Description	Code	Image
Theraband Tubing with Soft Grip Handles		
Soft Grip Handles 121cm Yellow Thin	FPP309910	
Soft Grip Handles 121cm Red Medium	FPP309911	
Soft Grip Handles 121cm Green Heavy	FPP309912	
Soft Grip Handles 121cm Blue Extra Heavy	FPP309913	
Soft Grip Handles 121cm Black Special Heavy	FPP309914	
Soft Grip Handles 121cm Silver Super Heavy	FPP309915	
Padded Cuffs Red Beginner/Intermediate	FPP309920	
Padded Cuffs Green Intermediate	FPP309921	
-Padded Cuffs Blue Intermediate/Advanced	FPP309922	
Padded Cuffs Black Advanced	FPP309923	

Exercise & Resistance Training

Description	Code	Image
Theraband		
152cm Yellow Thin	FPP308450	
152cm Red Medium	FPP308451	
152cm Green Heavy	FPP308452	
152cm Blue Extra Heavy	FPP308453	
152cm Black Special Heavy	FPP308454	
Exercise - Theraband - Door Anchors	FPP310100	
Weights - Hand Grips		
Medium	FPF290440	
Strong	FPF290450	
Extra Strong	FPF290460	
Exercise System - Digi Flex		
Yellow 0.7kg	FPP298700	
Red 1.4kg	FPP298800	
Green 2.3kg	FPP298900	
Blue 3.2kg	FPP299000	
Black 4.1kg	FPP299100	



Exercise & Resistance Training

Description	Code	Image
Putty - 85g		
Extra Soft - Tan	FPP298100	
Soft - Yellow	FPP298200	
Medium - Red	FPP298300	
Extra Firm - Blue	FPP298500	
Firm - Green	FPP298400	

Exercise & Resistance Training

Description	Code	Image
Weights - Soft Ankle / Wrist		
Wrist Weights 2x0.5kg	FPF290480	
Wrist Weights 2x1.0kg	FPF290470	
Wrist Weights 2x1.5kg	FPF290490	
Massager		
Mini Massager	FPP301600	
Backnobber	FPP303900	



Physiotherapy

Description	Code	Image
Weights - Soft Ankle / Wrist		
(Analogue) - PROTENS - Supplied w/ Electrodes + Leads	FPP315200	
(Digital) - Supplied w/ Electrodes + Leads + Batteries	FPP315210	
TENS Electrodes - For PROTENS Tens unit - 5cmx5cm (Pkt 4)	GSP342245	

Physiotherapy

Description	Code	Image
TENS Batteries - Protens	GSP342230	
TENS	GSP342240	

Bracing & Supports



Bracing & Supports

Description	Code	Image
-------------	------	-------

Walker Boots - Fixed Ankle

Standard Walker - Small	FBS301584	
Standard Walker - Medium	FBS3015842	
Standard Walker - X-Large	FBS3015844	

Walker Boots - Airselect Elite

X-Small	FPP303703	
Small	FPP303701	
Medium	FPP303700	
Large	FPP303699	
X-Large	FPP303702	
Short - X-Small	FPP303708	
Short - Small	FPP303706	
Short - Medium	FPP303705	
Short - Large	FPP303704	
Short - X-Large	FPP303707	

Walker Boots - XcelTrax Air


Tall - X-Small	FPP303709	
Tall - Small	FPP303710	
Tall - Medium	FPP303711	
Tall - Large	FPP303712	
Tall - X-Large	FPP303713	
Short - X-Small	FPP303714	

Bracing & Supports

Description	Code	Image
-------------	------	-------

Short - Small	FPP303715	
Short - Medium	FPP303716	
Short - Large	FPP303717	
Short - X-Large	FPP303718	
Standard Tall - X-Small	FBS303432	
Standard Tall - Small	FBS303433	
Standard Tall - Medium	FBS303434	
Standard Tall - Large	FBS303435	
Standard Tall - X-Large	FBS303436	
Standard Short - X-Small	FBS303437	
Standard Short - Small	FBS303438	
Standard Short - Medium	FBS303439	
Standard Short - Large	FBS303440	
Standard Short - X-Large	FBS303441	

Cast Shoe - Darco - Body Armor

X-Small	FBS301582	
Small	FBS3015821	
Medium	FBS3015822	
Large	FBS3015823	
X-Large	FBS3015824	

Cast Shoe - Darco - Slim Line

X-Small	FBS301583	
Small	FBS3015831B	
Medium	FBS3015832B	
Large	FBS3015833B	
X-Large	FBS3015834B	




Bracing & Supports

Description	Code	Image
-------------	------	-------

Cast Shoe - Rocker Sole - Standard

Large	FBS301583L	
Medium	FBS301583M	
Small	FBS301583S	

Knee Braces - Donjoy OA Reaction Web Brace

Right Medial/Left Lateral - XS	FBS303444	
Right Medial/Left Lateral - S	FBS303445	
Right Medial/Left Lateral - M	FBS303446	
Right Medial/Left Lateral - L	FBS303447	
Right Medial/Left Lateral - XL	FBS303448	
Right Medial/Left Lateral - XXL	FBS303449	
Right Medial/Left Lateral - XXXL	FBS303450	
Left Medial/Right Lateral - XS	FBS303451	
Left Medial/Right Lateral - S	FBS303452	
Left Medial/Right Lateral - Medium	FBS303453	
Left Medial/Right Lateral - Large	FBS303454	
Left Medial/Right Lateral - X-Large	FBS303455	
Left Medial/Right Lateral - XX-Large	FBS303456	
Left Medial/Right Lateral - XXX-Large	FBS303457	

Knee Brace Hinged - Oppo - Motion Control ROM

Open Knee - Length 46 cm	FBS301536	
Open Knee - Length 56 cm	FBS30153622	
Open Knee - Length 66 cm	FBS30153626	

Knee Brace - BodyAssist - Elastic Cross-Cut w/ Rods

Small	FBS284180	
Medium	FBS284190	
Large	FBS284200	
X-Large	FBS284210	
XX-Large	FBS284215	

Knee Brace - BodyAssist - Elastic Slip-on Support

Small	FBS284260	
Medium	FBS284270	
Large	FBS284280	
X-Large	FBS284290	

Bracing & Supports

Description	Code	Image
-------------	------	-------

Knee Brace - BodyAssist - Elastic Wrap Cartilage Support

Small	FBS284255	
Medium	FBS284256	
Large	FBS284257	
X-Large	FBS280765	

Knee Immobiliser - Oppo

Tri Panel (UNIVERSAL)	FBS301539	
-----------------------	-----------	---


Knee Brace - Oppo - Elastic Slip-on

O/Patella - Small	FBS301512	
O/Patella - Large	FBS301512L	
O/Patella - Medium	FBS301512M	
O/Patella - X-Large	FBS301512XL	
Extra Length - Medium - Beige	FBS280790	


Knee Brace - McDavid - Neoprene Ligament Support

Black - X-Large	FBS284881	
Black - XX-Large	FBS284882	
Black - Large	FBS284880	
Black - Small	FBS284860	
Black - Medium	FBS284870	

Calf Support - Oppo - Elastic Slip-on Shin

Brace - Small	FBS301548	
Support - Small	FBS301547	
Support - Medium	FBS301547M	
Support - Large	FBS301547L	
Support - X-Large	FBS301547XL	

Thigh Support - Oppo - Slip-on

Elastic - Small	FBS301545	
Neoprene - Small	FBS301544	
Neoprene - Medium	FBS301544M	
Neoprene - Large	FBS301544L	
Neoprene - X-Large	FBS301544XL	

Wrist Wrap

Wrist Wrap - Oppo - Elastic - (UNIVERSAL)	FBS301369	
Wrist Wrap - Oppo - Elastic w/ Silicone - (UNIVERSAL)	FBS301371	
Wrist Wrap - Oppo - Neoprene - (UNIVERSAL)	FBS301367	



Bracing & Supports

Description	Code	Image
Wrist Wrap - Oppo - Thumb Loop - Neoprene - (UNIVERSAL)	FBS301368	
Wrist Wrap - Oppo - Thumb Loop - Cotton Blend - (UNIVERSAL)	FBS301370	

Roylan Palm Protector w/ Finger Separator

Left (UNIVERSAL Size)	FBS282830	
Right (UNIVERSAL Size)	FBS282840	

Roylan Palm Shield - Foam

Left (UNIVERSAL Size)	FBS283010	
Right (UNIVERSAL Size)	FBS283020	

Roylan Palm Protector - Foam

Left (UNIVERSAL Size)	FBS282860	
Right (UNIVERSAL Size)	FBS282880	

Finger Splints - Body Assist

Multi Purpose - Rigid Splint - Small	FBS284070	
Curved or Straight - Small	FBS284060	

Baseball Finger Splints - Oppo

Small	FBS301466	
Medium	FBS3014662	
Large	FBS3014663	

Curved Finger Splints - Oppo

Small	FBS301467	
Medium	FBS3014672	
Large	FBS3014673	

Finger Splints - Aluminium - Oppo

Small	FBS301465	
Medium	FBS3014652	
Large	FBS3014653	

Bracing & Supports

Description	Code	Image
Open Thimble Finger Splints - Oppo		
Small	FBS301469	
Medium	FBS3014692	
Large	FBS3014693	

Handeze Support Gloves - Pair

X-Small - Beige	FBS301170	
Small - Beige	FBS301170S1	
Medium - Beige	FBS301170S2	
Large - Beige	FBS301170S3	

Oppo - Positioning splint

Right - Small	FBS301470	
Left - Large	FBS301470LL	
Left - Medium	FBS301470LM	
Left - Small	FBS301470LS	
Right - Large	FBS301470RL	
Right - Medium	FBS301470RM	

Roylan - Carve-it Hand Positioning splint

Large/XL	FBS283070	
Small/Medium	FBS283075	

Elbow Sleeves - Body Assist - Slip-on Elastic

Beige - Small	FBS284640	
Beige - Medium	FBS284650	
Beige - Large	FBS284660	
Beige - XL	FBS284670	

Arm Slings

Pouch - Body Assist - Small	FBS284580	
Pouch - Body Assist - Medium	FBS284590	
Pouch - Body Assist - Large	FBS284600	
Sling - Arm - Pouch - Body Assist - Deluxe & Swathe - Small-Medium	FBS284610	



Bracing & Supports

Description	Code	Image
Strap - Body Assist - Small	FBS284520	
Strap - Body Assist - Medium	FBS284530	
Strap - Body Assist - Large	FBS284540	
Tube - Body Assist - Small	FBS284550	
Tube - Body Assist - Medium	FBS284560	
Tube - Body Assist - Large	FBS284570	

Abdominal and Lumbar Support

Belt - BodyAssist - Small	FBS283680	
Belt - BodyAssist - Medium	FBS283690	
Belt - BodyAssist - Large	FBS283700	
Belt - BodyAssist - Extra Large	FBS283710	
Abdominal Support - Belt - BodyAssist - XXL	FBS280550	
Surgical Support Belt - Back Brace - Post Surgical Belt - Medium	FBS283740	
Back Brace - OAPL - Thermovent Neoprene - Large	FBS281890	
BodyAssist - Lumbo Sacral Support - Medium	FBS280525	
BodyAssist - Lumbo Sacral Support - Large	FBS280526	
BodyAssist - Lumbo Sacral Support - XL	FBS280527	

Back Brace - Black

Betabak - Lycra - Medium	FBS280470	
Betabak - Mesh - Medium	FBS280460	
Betabak - Lycra - Large	FBS280473	
Betabak - Lycra - X-Large	FBS280475	
Betabak - Lycra - X-Large	FBS280476	

Neck Support

Soft Foam Collar - Cervical Collar - Small	FBS301201	
Soft Foam Collar - Cervical Collar - Small / Long	FBS301202	
Soft Foam Collar - Cervical Collar - Medium	FBS301203	
Soft Foam Collar - Cervical Collar - Medium / Long	FBS301204	
Soft Foam Collar - Cervical Collar - Large -	FBS301205	

Bracing & Supports

Description	Code	Image
Oppo - Cervical Collar - Rigid Adjustable Collar - w/ chin - Small	FBS301210	
Oppo - Cervical Collar - Rigid Adjustable Collar - w/ chin - Large	FBS301210L	
Aspen Collar - Cervical Collar - Adult Short	FBS301211	
Aspen Collar - Cervical Collar - Adult Reg	FBS301212	
Aspen Collar - Cervical Collar - Adult Tall	FBS301213	
Aspen Collar - Cervical Collar - Adult X-Tall	FBS301214	
Aspen Collar - Cervical Collar - Set (one complete Collar + spare pads) - Adult Short	FBS301215	
Aspen Collar - Cervical Collar - Set (one complete Collar + spare pads) - Adult Reg	FBS301216	
Aspen Collar - Cervical Collar - Set (one complete Collar + spare pads) - Adult Tall	FBS301217	
Aspen Collar - Cervical Collar - Set (one complete Collar + spare pads) - Adult X-Tall	FBS301218	
Aspen Collar - Cervical Collar - Replacement pad Set - Adult - Small	FBS301219	
Aspen Collar - Cervical Collar - Replacement pad Set - Adult Regular	FBS301219AR	
Aspen Collar - Cervical Collar - Replacement pad Set - Adult tall	FBS301219AT	
Aspen Collar - Cervical Collar - Replacement pad Set - Adult X tall	FBS301219AXT	



Protection







Protection

Description	Code	Image
Heel Protection		
Heel - Bed Heel Elevator	PTA509050	
Heel - ROHO Heel Pad Single	PTC541600	
Heel / Elbow - Pad Convuluted Foam	PTA501300	
Heel Protectors - Silicone Fibre (Pair) **	PTA503350	
Heel Protectors Pelican	PTA504200	
Heel Protector - Heelift Suspension Boot - Smooth	PTA509200	
Heel Protector - Heelift Suspension Boot - Convuluted	PTA509100	

DermaSaver - Heel Elevator

Small	PTA533316	
Large	PTA533316L	

Protection

Description	Code	Image
Sheepskin Overlays		
Sheepskin Hi Temp - Classic Overlay 75cmx125cm - Blue	PTA521400	
Sheepskin Hi Temp - Medical Overlay XD1900 65cmx100cm - Blue	PTA521500	
Sheepskin Hi Temp - Paediatric Overlay 75cmx45cm	PTA520301	
Sheepskin Hi Temp - Cushion Seat or Back 45x50cm	PTA521300	
Sheepskin - Natural Shaped Gold 75cm x 100cm **	PTA502500	



Compression

Description	Code	Image
-------------	------	-------

Stocking Applicators

Stocking Accessories – "Compression Assist" limb lubricant - 60ml bottle	GMS336050	
Stocking Accessories - JOBST - Body Adhesive - 59ml Plastic Bottle	PTS641100	
Jobst Donning Gloves - Small	FBS303490	
Jobst Donning Gloves - Stocking Applicator - Medium	FBS303491	
Doff N Donner - Donning Cone (Black) for Doff N Donner	PTS642661	
Doff N Donner	PTS642660	
Doff N Donner & Cone	PTS642662	

Compression

Description	Code	Image
-------------	------	-------

TED Thigh Stockings w/ Belt - Anti-embolic - Toe Hole

Reg / X Small - White	PTS641425	
Reg / Small - White	PTS641430	
Reg / Medium- White	PTS641460	
Reg / Large - White	PTS641470	
Reg / XL - White	PTS641500	
Long / X Small - White	PTS641501	
Long / Small - White	PTS641502	
Long / Medium- White	PTS641503	
Long / Large - White	PTS641504	
Long / XL - White	PTS641505	

TED Thigh Stockings w/ Grip Top - Anti-embolic - Toe Hole

Short / Small - White	PTS641471	
Short / Large - White	PTS641472	
Short / XL - White	PTS641473	
Short / XXL - White	PTS641474	
Short / Medium- White	PTS641475	
Reg / Small - White	PTS641476	
Reg / Medium- White	PTS641477	
Reg / Large - White	PTS641478	
Reg / XL - White	PTS641479	
Reg / XXL - White	PTS641481	
Long / Small - White	PTS641482	
Long / Medium- White	PTS641483	
Long / Large - White	PTS641484	
Long / XL - White	PTS641485	
Long / XXL - White	PTS641486	



Compression

Description	Code	Image
Reg / Small - White	PTS641570	
Reg / Medium- White	PTS641580	
Reg / Large - White	PTS641590	
Reg / XL - White	PTS641591	
Reg / XXL - White	PTS641592	
Reg / XXXL - White	PTS641593	
Long / Small - White	PTS641600	
Long / Medium- White	PTS641610	
Long / Large - White	PTS641620	
Long / XL - White	PTS641621	
Long / XXL - White	PTS641622	
Long / XXXL - White	PTS641623	

TED Continuing Care Knee Stockings - Anti-embolic - C/Toe

Reg / Small - White	PTS641630	
Reg / Medium- White	PTS641631	
Reg / Large - White	PTS641632	
Reg / XL - White	PTS641633	
Long / Small - White	PTS641634	
Long / Medium- White	PTS641635	
Long / Large - White	PTS641636	
Long / XL - White	PTS641637	

Sigvaris Stockings - Cotton

AD Calf Open Or Closed Toe	PTS642310	
Extra Ad Calf Open or Closed Toe	PTS642315	
Calf Open/Closed Toe w/ Grip Top (ID CLASS, SIZE & COLOUR)	PTS642320	
Thigh Open/Closed Toe w/ Grip Top (ID CLASS, SIZE & COLOUR)	PTS642325	
Panty (ID CLASS, SIZE & COLOUR)	PTS642330	
Medium Chaps Open Toe (ID CLASS, SIZE & COLOUR)	PTS642335	

Sigvaris Stockings - Traditional

Traditional Calf Open/ Closed Toe (ID CLASS, SIZE & COLOUR)	PTS642340	
Traditional Half Thigh Open/ Closed Toe (ID CLASS, SIZE & COLOUR)	PTS642345	
Traditional Half Thigh Open/ Closed Toe w/ Grip Top (ID CLASS, SIZE & COLOUR)	PTS642350	
Traditional Thigh Open/ Closed Toe (ID CLASS, SIZE & COLOUR)	PTS642355	

Compression

Description	Code	Image
Traditional Thigh Open/ Closed Toe w/ Grip Top (ID CLASS, SIZE & COLOUR)	PTS642360	
Traditional Chaps Open Toe (ID CLASS, SIZE & COLOUR)	PTS642365	
Traditional Panty (ID CLASS, SIZE & COLOUR)	PTS642370	
Traditional Panty For Men (ID CLASS, SIZE & COLOUR)	PTS642375	

Sigvaris Stockings - Magic

Magic Calf Open/Closed Toe (ID CLASS, SIZE & COLOUR)	PTS642380	
Magic Thigh Open/Closed Toe w/ Grip Top (ID CLASS, SIZE & COLOUR)	PTS642385	
Magic Panty (ID CLASS, SIZE & COLOUR)	PTS642390	

Sigvaris Stockings - Ulcrex

Kit Calf Open Toe (ID CLASS, SIZE & COLOUR)	PTS642395	
Underliners Calf Open Toe (ID CLASS, SIZE & COLOUR)	PTS642400	

Sigvaris Stockings - Calf - Open Toe

Post Surgery Kit (ID CLASS, SIZE & COLOUR)	PTS642405	
Underliners (ID CLASS, SIZE & COLOUR)	PTS642410	

JOBST Relief Stockings - Knee High - Open Toe - 15-20mmHg

Small/Beige	PTS642516	
Med/Beige	PTS642517	
Large/Beige	PTS642518	
XLarge/Beige	PTS642519	
Large-Full Calf/Beige	PTS642520	
XLarge-Full Calf/Beige	PTS642521	

JOBST Relief Stockings - Knee High - Closed Toe - 15-20mmHg

Small/Beige	PTS642522	
Med/Beige	PTS642523	
Large/Beige	PTS642524	
XLarge/Beige	PTS642525	
Large-Full Calf/Beige	PTS642526	
X-Large-Full Calf/Beige	PTS642527	



Compression

Description	Code	Image
-------------	------	-------

JOBST Knee Stockings - For Men - C/Toe - 15-20mmHg

Small - Black	PTS640040	
Medium- Black	PTS640050	
Large - Black	PTS640060	
XL - Black	PTS640070	

JOBST Knee Stockings - UltraSheer - C/Toe - 15-20mmHg

Small - Black	PTS640140	
Medium- Black	PTS640150	
Large - Black	PTS640160	
XL - Black	PTS640170	
Small - Natural	PTS640220	
Medium- Natural	PTS640230	
Large - Natural	PTS640240	
XL - Natural	PTS640250	

JOBST Stockings - UltraSheer - Thigh w/ Grip Top - 15-20mmHg

Medium- Black	PTS640300	
Large - Black	PTS640310	
XL - Black	PTS640315	
Small - Natural	PTS640330	
Medium- Natural	PTS640335	
Large - Natural	PTS640350	
XL - Natural	PTS640360	

JOBST Stockings - UltraSheer - Waist - C/Toe - P/hose - 15-20mmHg

Small - Black	PTS640380	
Medium- Black	PTS640390	
Large - Black	PTS640400	
XL - Black	PTS640410	
Small - Natural	PTS640450	
Medium- Natural	PTS640460	
Large - Natural	PTS640470	
XL - Natural	PTS640480	

JOBST Knee Stockings - For Men - C/Toe - 20-30mmHg

Small - Black	PTS640530	
Medium- Black	PTS640540	
Large - Black	PTS640550	
XL - Black	PTS640560	
Small - Black	PTS640601	
Medium- Black	PTS640602	
Large - Black	PTS640603	
XL - Black	PTS640604	
Small - Black	PTS640605	
Medium- Black	PTS640606	
Large - Black	PTS640607	
XL - Black	PTS640608	

Compression

Description	Code	Image
-------------	------	-------

JOBST Knee Stockings - UltraSheer - C/Toe - 20-30mmHg

Small - Black	PTS640610	
Medium- Black	PTS640620	
Large - Black	PTS640630	
XL - Black	PTS640640	
Small - Natural	PTS640650	
Medium- Natural	PTS640660	
Large - Natural	PTS640670	
XL - Natural	PTS640680	

JOBST Stockings - UltraSheer - Thigh w/ Grip Top - C Toe - 20-30mmHg

Small - Black	PTS640710	
Medium- Black	PTS640720	
Large - Black	PTS640730	
XL - Black	PTS640735	
Small - Natural	PTS640740	
Medium- Natural	PTS640750	
Large - Natural	PTS640760	
XL - Natural	PTS640761	

JOBST Stockings - UltraSheer - Waist - C/Toe - 20-30mmHg

Small - Black	PTS640779	
Medium- Black	PTS640780	
Large - Black	PTS640781	
XL - Black	PTS640782	
Small - Natural	PTS640789	
Medium- Natural	PTS640790	
Large - Natural	PTS640800	
XL - Natural	PTS640810	

JOBST Knee Stockings - For Men - C/Toe - 30-40mmHg

Small - Black	PTS640830	
Medium- Black	PTS640840	
Large - Black	PTS640850	
XL - Black	PTS640860	

JOBST Knee Stockings - UltraSheer - C/Toe - 30-40mmHg

Small - Black	PTS640870	
Medium- Black	PTS640880	
Large - Black	PTS640890	
XL - Black	PTS640900	
Small - Natural	PTS640910	
Medium- Natural	PTS640920	
Large - Natural	PTS640930	
XL - Natural	PTS640940	



Compression

Description	Code	Image
JOBST Stockings - UltraSheer - Thigh w/ Grip Top - C/Toe - 30-40mmHg		
Small - Black	PTS640945	
Medium- Black	PTS640946	
Large - Black	PTS640947	
XL - Black	PTS640948	
Small - Natural	PTS640949	
Medium- Natural	PTS640950	
Large - Natural	PTS640955	
XL - Natural	PTS640960	

JOBST Stockings - UltraSheer - Waist - P/hose - C/Toe - 30-40mmHg

Small - Natural	PTS640969	
Medium- Natural	PTS640970	
Large - Natural	PTS640980	
XL - Natural	PTS640981	
Small - Black	PTS640983	
Medium- Black	PTS640984	
Large - Black	PTS640985	
XL - Black	PTS640986	

Stockings - JOBST UlcerCare - Starter Kit (Stocking + 2 Liners) - 40mmHg

Right Zip - Small	PTS640990	
Right Zip - Med	PTS641000	
Right Zip - Large	PTS641010	
Right Zip - XL	PTS641020	
Right Zip - XXL	PTS641021	
Replacement Liners (pkt 3) - 17mmHg - Small	PTS641025	
Replacement Liners (pkt 3) - 17mmHg - Med	PTS641030	
Replacement Liners (pkt 3) - 17mmHg - XL	PTS641031	
Replacement Liners (pkt 3) - 17mmHg - XXL	PTS641032	
Replacement Liners (pkt 3) - 17mmHg - Large	PTS641040	

Stockings - JOBST Relief - Knee - 20-30mmHG - Unisex

C/Toe - Small	PTS641110	
C/Toe - Med	PTS641120	
C/Toe - Large	PTS641130	
C/Toe - XL	PTS641140	
C/Toe - Large Full Calf	PTS641141	
C/Toe - XL/Full Calf	PTS641142	
O/Toe - Small	PTS641150	
O/Toe - Med	PTS641160	

Compression

Description	Code	Image
O/Toe - Large	PTS641170	
O/Toe - XL	PTS641180	
O/Toe - Large/ Full Calf	PTS641181	
O/Toe - XL/ Full Calf	PTS641182	

Stockings - JOBST Relief - Knee - 30-40mmHG - Unisex

O/Toe - Small	PTS641190	
O/Toe - Med	PTS641200	
O/Toe - Large	PTS641210	
O/Toe - XL	PTS641220	
O/Toe - Large Full Calf	PTS641221	
O/Toe - XL /Full Calf	PTS641222	
C/Toe - Small	PTS641230	
C/Toe - Med	PTS641240	
C/Toe - Large	PTS641250	
C/Toe - XL	PTS641260	
C/Toe - Large Full Calf	PTS641261	
C/Toe - XL Full Calf	PTS641262	

Stockings - JOBST Relief - Thigh w/Grip Top - 20-30mmHG - Unisex

O/Toe - Small	PTS641265	
O/Toe - Med	PTS641270	
O/Toe - Large	PTS641280	
O/Toe - XL	PTS641290	
C/Toe - Small	PTS641300	
C/Toe - Med	PTS641310	
C/Toe - Large	PTS641319	
C/Toe - XL	PTS641320	

Stockings - JOBST Relief - Waist - P/hose - 20-30mmHG - Unisex

C/Toe - Med	PTS641330	
C/Toe - XL	PTS641335	
O/Toe - Small	PTS641339	
O/Toe - Med	PTS641340	
O/Toe - Large	PTS641341	
O/Toe - XL	PTS641342	
C/Toe - Small	PTS641938	
C/Toe - Large	PTS641939	

Stockings - JOBST Relief - Thigh w/ Belt (Chaps) - 20-30mmHG - Unisex


Left - O/Toe - Small	PTS641349	
Left - O/Toe - XL	PTS641350	
Left - O/Toe - Med	PTS641351	
Left - O/Toe - Large	PTS641352	
Right - O/Toe - Small	PTS641359	
Right - O/Toe - Med	PTS641360	
Right - O/Toe - Large	PTS641365	
Right - O/Toe - XL	PTS641370	



Compression

Description	Code	Image
-------------	------	-------


Stockings - JOBST Relief - Thigh w/ Belt (Chaps) - 30-40mmHG - Unisex

Left - O/Toe - Small	PTS641371	
Left - O/Toe - Large	PTS641372	
Left - O/Toe - XL	PTS641373	
Right - O/Toe - Small	PTS641374	
Right - O/Toe - Med	PTS641375	
Right - O/Toe - Large	PTS641376	
Right - O/Toe - XL	PTS641377	
Left - O/Toe - Med	PTS641380	

Stockings - JOBST Relief - Thigh w/ Grip - 30-40mmHG - Unisex

O/Toe - Small	PTS641940	
O/Toe - Med	PTS641941	
O/Toe - Large	PTS641942	
O/Toe - XL	PTS641943	
C/Toe - Small	PTS641944	
C/Toe - Med	PTS641945	
C/Toe - Large	PTS641946	
C/Toe - XL	PTS642264	


Stockings - JOBST Relief - Waist - P/hose - 30-40mmHG - Unisex

O/Toe - Small	PTS641947	
O/Toe - Med	PTS641948	
O/Toe - Large	PTS641949	
O/Toe - XL	PTS641950	

Stockings - JOBST UlcerCare - Starter Kit - 40mmHG

(Stocking + 2 Liners) - Left Zip - Small	PTS641951	
(Stocking + 2 Liners) - Left Zip - Med	PTS641952	
(Stocking + 2 Liners) - Left Zip - Large	PTS641953	
(Stocking + 2 Liners) - Left Zip - XL	PTS641954	
(Stocking + 2 Liners) - Left Zip - XXL	PTS641955	


JOBST Bella Strong - Gauntlet 20 - 30mmHg

Size 1	PTS642534	
Size 2	PTS642535	
Size 3	PTS642536	
Size 4	PTS642537	
Size 5	PTS642538	
Size 6	PTS642539	
Size 7	PTS642540	
Size 8	PTS642541	


Compression

Description	Code	Image
-------------	------	-------

Sigvaris - Traditional Sleeves

Traditional - w/out Mitten	PTS642415	
Traditional - w/ Mitten	PTS642420	
Mitten - Traditional Mitten	PTS642445	

Sigvaris - Advance Sleeves

Advance - w/o Mitten	PTS642425	
Advance - w/ Mitten	PTS642430	
Advance - w/o Mitten Grip Top	PTS642435	
Advance - w/ Mitten Grip Top	PTS642440	


JOBST Elvarex - Armsleeve - 20-30mmHg

Right - Small	PTS642149	
Left - Small	PTS642150	
Right - Med	PTS642151	
Left - Med	PTS642152	
Right - Large	PTS642153	
Left - Large	PTS642154	


JOBST Elvarex - Armsleeve with Shoulder Strap - 20-30mmHg

Right - Small	PTS642155	
Left - Small	PTS642156	
Right - Med	PTS642157	
Left - Med	PTS642158	
Right - Large	PTS642159	
Left - Large	PTS642160	

JOBST Elvarex - Armsleeve with Bra Attachment - 20-30mmHg

Right - Small	PTS642161	
Left - Small	PTS642162	
Right - Med	PTS642163	
Left - Med	PTS642164	
Right - Large	PTS642165	
Left - Large	PTS642166	

JOBST Elvarex - Handpiece with Thumb - 20-30mmHg

Right - Small	PTS642167	
Left - Small	PTS642168	
Right - Med	PTS642169	
Left - Med	PTS642170	
Right - Large	PTS642171	
Left - Large	PTS642172	

JOBST Elvarex - Handpiece with Full Fingers - 2 - 30mmHg

Right - Small	PTS642173	
Left - Small	PTS642174	
Right - Med	PTS642175	
Left - Med	PTS642176	
Right - Large	PTS642177	
Left - Large	PTS642178	



1300 133 120
www.aidacare.com.au

