

MULTILINE XXL

THE BARIATRIC TREATMENT CHAIR FOR WARD, THEATRE AND CONSULTING SUITES



300kg

Safe working load



Exclusively Distributed by:



GREINER



MULTILINE XXL

MULTILINE XXL.

This heavy duty couch meets the bariatric needs of oncology day hospitals, infusion medicine and transfusion medicine. It combines the comfort and security of oncology couches with the ability to facilitate bariatric treatment. Patients and those caring for them enjoy the full range of amenities and features, with the couch quickly conveying patients into the most suitable position for therapy. With a wide seat width of 70cm and ability to transfer patients from a sitting to standing position, the XXL is a true bariatric solution. The comfortable, extra soft upholstery, patient friendly armrests that support therapy and wide range of fittings make the MULTILINE XXL a next-generation couch.



A TRUE BARIATRIC SOLUTION.

The development of the XXL chair ensures that bariatric treatment in all areas of care and application is simplified. As a result, the XXL chair can be universally used for transport, examinations and treatments, as a nursing chair and for rest or recovery. The armrests can be easily removed for patient transfer, side entry, treatment and relocation.

GREINER for oncology:

- ◆ Totally new concept in terms of design, technology and features; smooth-running mobility with central locking if required
- ◆ Version with three motors for the stepless adjustment of backrest, seat and leg section into required positions, including the emergency Trendelenburg position
- ◆ Comfortable cushioning in a wide range of colours, high levels of seated and reclining comfort, maximum hygiene with a minimum of seams
- ◆ Unique design with 90° angle of leg section – perfect for chair-like access and exit, minimal space requirement
- ◆ Flexible design variants and configurations adapt to every practical purpose
- ◆ Independently certified product safety and superb reliability



MULTILINE XXL. BARIATRIC TREATMENT SOLUTION.

The new therapy couch developed by Greiner packs the performance of technically sophisticated couches into a bariatric solution. Designed for patients who may be very weak or afflicted by extreme symptoms, the couch helps them through the difficult times. It features ease of access, smooth rollers, clear functions, stepless motorised movement to the required position, comfortable armrests and cosy upholstery in soothing colours. The emphasis is on patient comfort, ergonomics and efficiency to ensure the best treatment conditions for patients over hours of chemotherapy.



GREINER hand control unit with intelligent operation: simple, internationally understandable symbols; keys for basic settings for seated, relax, flat bed and Trendelenburg positions are easily accessed in an emergency. Safety is assured thanks to a magnetic key lock.

Comfort, functionality and ergonomics.

With stepless motorised adjustment and intelligent and intuitive operation via the selection and control keys of the GREINER control element, the MULTILINE XXL for infusion and transfusion therapy promises high levels of seated and reclining comfort in all situations. The three-motor couch quickly and gently moves to any required or preferred position – from seated, relax and bed positions to the emergency Trendelenburg position – at the touch of a button. Access and exit is considerably eased by the vertical leg section and adjustable access height of 55-80cm.

The GREINER design concept makes it easy for carers and patients alike to use the powerful motors for the backrest, seat and leg sections efficiently. The couch is mobile, stable and secure.



Sit-to-stand position



Seated position



Relax position



Bed/flat position



Trendelenburg and shock position



MULTILINE NEXT XXL. NEW LEVEL OF COMFORT FOR INFUSION AND TRANSFUSION.

GREINER treatment couches are supplied with many inclusions as standard, however can be tailored to your precise needs. Premium covering materials and upholstery offer extra soft, wear resistant surfaces; a dynamically flexible suspension system within the seating surface supports freedom of movement for the spine, pelvis and hips, thus relieving pressure to the benefit of patients.



Electric sit-to-stand function



Electric patient positioning



Head rest as standard




Fold-down armrests



Available in Black Vinyl.
A range of custom colours are available upon request. (Lead times may vary)



Central locking castors



Adjustable push handles

Hygiene

Design compliant with stringent regulations, ease of access, clean surface finish, tested surfaces, minimum of seams

Stability of movement

Fixing of backrest in any position makes couch easier to push

Comfort

Comfortable upholstery with wide 70cm hand made seat

Flexibility

Infusion pole holder can be pivoted left and right with just one hand

Sit-to-Stand

Ability for patients to be easily moved from a seating to standing position

Access

Ease of access from front or side with fold down armrest and adjustable access height of 55-80cm

Seated position

Genuine chair position with leg section at 90° angle



Mobility and braking

10 cm twin-wheel castors with central locking if required

MULTILINE XXL: Equipment and accessories

Product Codes (treatment chair & inclusions - listed below)

- ◆ CHT100300 - XXL Treatment Chair - Black Vinyl

Inclusions

- ◆ Battery backup
- ◆ Adjustable push handles
- ◆ Padded armrests
- ◆ IV pole holder
- ◆ Head cushion

Additions

- ◆ Premium vinyl
- ◆ Patient table
- ◆ Seat with pressure relieving springs
- ◆ Protective legrest cover

Dimensions and technical data

◆ Chair upright (floor space)	110 x 78 cm
◆ Recline position (floor space)	196 x 78 cm
◆ Total width	89 cm
◆ Surface in reclined position (L x W)	185 x 70 cm
◆ Thickness of padding	7 cm
◆ Armrest (L x W)	70 x 7 cm
◆ Access height	55 - 80 cm
◆ Castor diameter	15 cm
◆ Backrest	- 5° to 78°
◆ Seat section	0° to 22°
◆ Leg section	0° to 90°
◆ Maximum patient weight	300 kg
◆ Voltage	230 V/50-60 Hz, 110 V/50-60 Hz
◆ Weight of couch	148 kg



GERMAN DESIGN AND QUALITY

Made by GREINER means made in Germany. GREINER has been supplying products for the medical field since the 1970s. A consistently high standard is applied in the development and production of GREINER couches. The same applies to our range of services. We achieve this by translating our flair for innovation into fully automated processes and handmade quality. High standards in design, materials, engineering and ergonomics combined with client-oriented flexibility make GREINER a partner to quality-conscious customers at the highest level.



we design quality.

Our quality commitment is verified by an independently audited quality assurance system: GREINER is certified under DIN EN ISO 9001 and 13485.



Test mark for
quality management system



Test mark for
product certification

Exclusively Distributed by:



Aidacare P/L

4 Noonan Rd,
Ingleburn NSW 2565
T – 1300 133 120

W – www.aidacare.com.au

